



Jun 27, 2022

Isolation of ECs from Brain tissue for scRNAseq on the 10x Chromium

Guido Krähenbühl¹¹MDC Berlin*In Development*

This protocol is published without a DOI.

AG Gerhardt



Guido Krähenbühl

ABSTRACT

Protocols for endothelial cell isolation from mouse tissues: brain, adapted from

Open Access • DOI: 10.1016/j.xpro.2021.100508

PROTOCOL CITATION

Guido Krähenbühl 2022. Isolation of ECs from Brain tissue for scRNAseq on the 10x Chromium. **protocols.io**
<https://protocols.io/view/isolation-of-ecs-from-brain-tissue-for-scrnaseq-on-cbh2sj8e>



LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jun 17, 2022

LAST MODIFIED

Jun 27, 2022

PROTOCOL INTEGER ID

64794

Preparation

- 1 Day before experiment check Antibodies, Enzyme Mix 1 and 2 Reagents to be used and prepare the following Buffer:

A	B	C
BSA	0.5%	2.5 g
EDTA, stock: 0.5M	2 mM	2 mL
PBS	-	498 mL
Total		500 mL

MACS buffer, store at 4C for up to 3 months (B is Final Concentration)

A	B	C
BSA	2%	10 mL
EDTA, stock: 0.5M	2 mM	2 mL
PBS	-	488 mL
Total		500 mL

FACS buffer, store at 4C for up to 3 months

- Before starting prepare:
Dissection Tools
20G needle and thin syringe
If Brain tissue is also taken, prepare perfusion kit
10mL DMEM Genta

A	B	C	D	E
	Enzyme Mix 1		Enzyme Mix 2	
	Buffer X	Enzyme P	Buffer Y	Enzyme A
1 Brain	1900 ul	50 ul	20 ul	10 ul
2 Brains	3800 ul	100 ul	40 ul	20 ul
5 Brains	9500 ul	250 ul	100 ul	50 ul

Keep all reagents on ice

Tissue Collection


- Sacrifice mice using Isoflurine
(do not use cervical dislocation as this might affect perfusion of the brain)
- Make an incision in the skin with fine scissors about a couple of millimeters above the orifice.
From this opening the incision must continue on the mid ventral side up to the chin.


- 5 Perform transcardial perfusion via the left ventricle with **8 mL Ice-cold PBS** / mouse at a perfusion rate of 2 mL/min for 5 min.
- 6 If cervical lymph nodes are needed dissect them now, otherwise remove the Head and harvest Brain regions of interest and collect all Brains in **10 mL Ice-cold DMEM, genta**
- 7 If lymph nodes are also analysed **go to step #3** and repeat for up to 5 mice.
If only brain ECs are analysed a total of 2 mice should suffice


Tissue processing 15m





- 8 Mince Brain samples with a scalpel in a 10cm dish until they may pass through a cut 1000P pipette
- 9 Transfer minced sample to a 15 mL conical tube and add **1950 µL Enzyme mix 1 per Brain** at **10 rpm, 37°C, 00:17:00** on the tube rotator
- 10 Add **30 µL enzyme mix 2 per Brain**
- 11 Pipette 10 times up and down with a pasteur pipette (until the suspension passes without clogging)
- 12 Incubate at **10 rpm, 37°C, 00:12:00** on the tube rotator
- 13 Pass the cell suspension 10 times through a 20G needle on a 1mL syringe (dont create bubbles or foam)
- 14 Incubate at **10 rpm, 37°C, 00:10:00** on the tube rotator



15 Transfer cell suspension to a 50 mL conical tube, add 40 mL HBSS and filter through a 70 µm cell strainer.


16 Centrifuge at  **300 x g, 4°C, 00:05:00** 5m

16.1 Transfer supernatant to new 50 mL conical tube, keep pellet on ice and centrifuge supernatant at  **300 x g, 4°C, 00:05:00** 5m

16.2 Transfer supernatant to new 50 mL conical tube, keep pellet on ice and centrifuge supernatant at  **300 x g, 4°C, 00:05:00** 5m


17 Discard Supernatant and pool 3 pellets in  **3600 µL MACS buffer** per Brain, add  **400 µL Myelin Removal MicroBeads** per Brain and incubate at  **4 °C** for  **00:15:00** 15m

18 Add a total volume of  **40 mL MACS buffer** and centrifuge at  **300 x g, 4°C, 00:05:00** 5m
Prepare 6 times the number of LS columns per Brain by rinsing each with 3mL MACS buffer





19 Resuspend in  **6 mL MACS buffer** per Brain sample and run 1mL per LS column (collect elluent)

19.1 Wash each column twice with 1mL MACS buffer (collect elluent)

20 Centrifuge total eluent at  **300 x g, 4°C, 00:05:00** , remove supernatant 5m

21 Pool Pellets in  **2 mL MACS buffer** per Brain and centrifuge at  **300 x g, 4°C, 00:05:00** 5m

FACS sorting 15m

- 22 Resuspend the pellet in  **500 µL FACS buffer** per brain
- 23 Prepare Controls and Samples with each antibody at 1/200:
Cells only
CD31
CD45
Viability dye
PDPN
- 24 Stain Cells for  **00:30:00** at  **4 °C** 30m
- 25 Add 3mL of FACS Buffer, centrifuge at  **300 x g, 4°C, 00:05:00** and resuspend in 2 mL of ^{5m}FACS Buffer