



Upload image

Jul 20, 2020

Nickel-NTA Protein Purification

Nay Chi Khin¹, Diep R Ganguly¹, Suyan Yee¹, Kai Xun Chan¹¹The Australian National University

1 Works for me dx.doi.org/10.17504/protocols.io.bihhkb36

Pogson Genomics Group



Diep Ganguly

The Australian National University



ABSTRACT

Protocol to purify 6xHis-tagged recombinant proteins expressed in *E. coli* using Ni-NTA pull-down.

DOI

dx.doi.org/10.17504/protocols.io.bihhkb36

PROTOCOL CITATION

Nay Chi Khin, Diep R Ganguly, Suyan Yee, Kai Xun Chan 2020. Nickel-NTA Protein Purification. **protocols.io**
dx.doi.org/10.17504/protocols.io.bihhkb36

KEYWORDS

Ni-NTA, Protein purification

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jul 13, 2020

LAST MODIFIED

Jul 20, 2020

PROTOCOL INTEGER ID

39177

GUIDELINES

Try to perform all steps at 4 °C as best as possible and keep protein samples on ice. For all spin-concentrate and buffer exchange steps, use slow speed (e.g. 3500-5000 rcf) to prevent precipitation.

MATERIALS

NAME	CATALOG #	VENDOR
Sodium dihydrogen phosphate		P212121
Disodium hydrogen phosphate		
Triton-X100		
1.5 mL Eppendorf tubes		
Glycerol	GB0232.SIZE.500ml	Bio Basic Inc.
Poly-Prep Chromatography Columns	731-1550	BioRad Sciences
DTT	43816-10ML	Sigma Aldrich
Sodium chloride	S3014	Sigma Aldrich
Imidazole	I5513	Sigma Aldrich

NAME	CATALOG #	VENDOR
Bovine Thrombin Protein	1996000JL	Thermo Fisher
HisPur [®] ; Ni-NTA Resin	88221	Thermo Fisher
Pierce [®] ; HRV 3C Protease Solution Kit (2 units/µL)	88946	Thermo Fisher
TEV protease	P8112S	New England Biolabs
Tris Hydrochloride (Tris-HCl)	RES3098T-B7	Sigma Aldrich
Amicon [®] Ultra-4 Centrifugal Filter Unit	UFC8003	Merck Millipore
Amicon [®] Ultra-15 Centrifugal Filter Unit	UFC9003	Merck Millipore

BEFORE STARTING

<https://www.protocols.io/private/3797D176057F0C37477116D259C75682>

Prepare buffers

1h

1h

1 Buffer recipes

Resuspension buffer

- 50 mM Tris-HCl, pH 8.0
- 2 mM EDTA

1X Binding buffer

- 50 mM sodium phosphate buffer, pH 8.0 (ensure >1 pH units away from pI of expressed protein)
- 500 mM sodium chloride
- 0.5% Triton X-100
- 10% glycerol
- 10 mM imidazole

N.b. This should reflect the cell lysis buffer (e.g.

<https://www.protocols.io/private/3797D176057F0C37477116D259C75682>).

1X Wash buffer

- 50 mM sodium phosphate buffer, pH 8.0
- 500 mM sodium chloride
- 0.5% Triton X-100
- 10% glycerol
- 20 mM imidazole

1X Elute buffer

- 50 mM sodium phosphate buffer, pH 8.0
- 500 mM sodium chloride
- 0.5% Triton X-100
- 10% glycerol
- 250 mM imidazole

N.b. DTT cleaves/interferes with Ni²⁺-binding, use low DTT concentrations (<2 mM), in the above buffers, if needed.

Digestion buffer

- 50 mM Tris-HCl, pH 8.0
- 150 mM NaCl

For *storage*, supplement *digestion buffer* with:

- 1 mM DTT
- 20% glycerol

Ni-NTA Affinity Purification

35m

15m

- Keep solubilised protein samples (from completion of <https://www.protocols.io/private/3797D176057F0C37477116D259C75682>) on ice to thaw completely.

- 3 Centrifuge samples (14,000 rcf for 10-20 minutes) and syringe filter (0.45 µm) supernatant into clean microfuge tube^{30m}. Keep on ice.
- 4 Add an appropriate amount of Ni-NTA resin slurry (50% slurry in 20% ethanol) into gravity-flow column. For 15 mL bacterial cultures resuspended in 2 mL lysis buffer, equating to ~1-4 mg protein (1-2 mg/mL), 500 µL resin slurry is added.

Top-up column with 5x volume of *1X binding buffer* (e.g. 500 µL slurry = 250 µL resin = ~2 mL 1X binding buffer).

Allow binding buffer to drop into a waste bottle to equilibrate the resin.
- 5 Transfer syringe filtered soluble protein (i.e. supernatant from step 3) into equilibrated resin, cap column tightly at both ends (including parafilm to prevent leakage) and incubate with shaking for 1-2 hours @ 4°C (max 10 mL per column)^{2h}.
- 6 Keep column up-right and allow resin to settle. Once settled, open cap and capture flow-through (FT).
- 7 Wash resin three-times with *1X Wash buffer* (2x resin volume, i.e. 1 mL) and collect separately (i.e. W1, W2, W3).
- 8 Elute four-times with *1X Elute buffer* (1X resin volume, i.e. 500 µL), collect separately (i.e. E1, E2, E3, E4).
- 9 Run SDS-PAGE of pre-FT (20 µL), FT (20 µL), W1-W3 (20 µL), and E1-E4 (8 µL) to quality-check purification.
- 10 Store all samples @ 4 °C for the short-term. Keep columns containing resin for second Ni-NTA step.

Protease digestion 16h

- 11 Pool elutes containing protein fragment of expected size.
- 12 Perform buffer exchange via centrifuge filtration (or dialysis) into *Digestion buffer* to reduce [imidazole] < 1mM.

For example, if *1X Elution buffer* contains 250 mM imidazole, you will need 4x centrifugation at a dilution of 1:4 per spin (4 °C @ 3,500 rcf for 20 minutes).
- 13 Add appropriate protease (e.g. TEV, Thrombin, HRV 3C) and incubate overnight (or > 16hrs) at 4 °C. Ensure to add any required cofactors (e.g. citrate for TEV)^{16h}.
- 14 Equilibrate Ni-NTA resin (from previously used and stored columns) with 10-20 mL *Digestion buffer* and allow resin to settle. Let buffer flow-through.

- 15 Add digested protein to resin and collect flow-through (dig-FT). Can perform multiple elutes with *Digestion buffer* to obtain as much protein as possible.
- 16 Run elutes on SDS-PAGE gel to test for cleavage.
- 17 Pool and concentrate the pure fraction(s) using buffer exchange.
- 18 Prepare aliquots of purified protein supplemented with 10-20% glycerol and 1mM DTT. Store at -20 °C.