



Jun 04, 2025

# 🌐 Transfecting COS-1 cells using FuGENE® 4K Transfection Reagent and cell lysis (M-PER) in a 6-well plate

DOI

[dx.doi.org/10.17504/protocols.io.261ged1rov47/v1](https://dx.doi.org/10.17504/protocols.io.261ged1rov47/v1)

Jailyn Izu<sup>1</sup>

<sup>1</sup>NIMH

TIU



Jailyn Izu

NIMH

OPEN  ACCESS



**DOI:** [dx.doi.org/10.17504/protocols.io.261ged1rov47/v1](https://dx.doi.org/10.17504/protocols.io.261ged1rov47/v1)

**Protocol Citation:** Jailyn Izu 2025. Transfecting COS-1 cells using FuGENE® 4K Transfection Reagent and cell lysis (M-PER) in a 6-well plate. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.261ged1rov47/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**Created:** January 03, 2024

**Last Modified:** June 04, 2025

**Protocol Integer ID:** 92930

## Abstract


General protocol for transfecting COS-1 cells with FuGene.




## Materials

- 6-well plate (Corning, #3506)
- FuGene 4K Transfection Reagent (Promega, #E5911)
- Opti-MEM Reduced Serum Medium (Thermo Fisher Scientific, #31985062)
- Plasmid constructs
- Dulbecco's Modified Eagle Medium (DMEM) with GlutaMax (Thermo Fisher Scientific, #10569010)
- FBS (Gibco, #26140-079)
- Penicillin-streptomycin (Thermo Fisher Scientific, #15140122)
- PBS, pH 7.4 (Thermo Fisher Scientific, #10010023)
- M-PER mammalian protein extraction reagent (Thermo Fisher Scientific, #78501)
- Halt Protease Inhibitor Cocktail (Thermo Fisher Scientific, #78429)


## Protocol materials

 FuGENE® 4K Transfection Reagent **Promega Catalog #E5911**

 FuGENE® 4K Transfection Reagent **Promega Catalog #E5911**











## Seed cells

- 1 Cos1 cell line were purchased from (ATCC # CRL-1650) and cultured in DMEM with GlutaMax **DMEM with GlutaMax (Thermo Fisher Scientific, # 10569010)** supplemented with 10% FBS **Fetal Bovine Serum, certified # 16000044** and 1% **penicillin-streptomycin (Thermo Fisher Scientific, # 15140122)**.
- 2 Cells were maintained at 37 °C in 5% CO<sub>2</sub> and split 1:5 every 3-5 days based on the confluence (80-100 %) reached through visual monitoring. Cells were not passaged more than 30 times.
- 3 Seed COS-1 cells in a 6-well plate (*see specific protocols for seeding density and total volume amounts*). Incubate the plated cells in  37 °C with 5% CO<sub>2</sub> (see specific protocols for duration).

## Transfection

2d

- 4 Once cells have grown to 50-80% confluence, transfect. In a biosafety cabinet, combine FuGene and Opti-Mem (referred to in the Promega protocol as just "medium," volumes according to the manufacturer's protocol  
 FuGENE® 4K Transfection Reagent **Promega Catalog #E5911** ) in a  1.5 mL eppendorf tube. Vortex gently for 1-2 seconds and incubate at room temperature for  00:05:00 .
- 5 Add plasmid constructs to each tube containing FuGene:Opti-MEM (ratio according to the manufacturer's protocol  
 FuGENE® 4K Transfection Reagent **Promega Catalog #E5911** ). Vortex briefly and incubate at room temperature for  00:15:00 .
- 6 While incubating, the DNA with FuGene:Opti-MEM, aspirate the old media from the cells. Add  1 mL of fresh DMEM + 10% FBS + 1% pen-strep to each well.
- 7 After incubation, pipette  100 µL of the reaction mixture into the corresponding wells. Tilt the plate(s) back and forth to mix.
- 8 Incubate the plate(s) in a  37 °C incubator with 5% CO<sub>2</sub> according to the manufacturer's protocol.

5m








15m

10m



## Lysis

10m

- 9 Aspirate the media from all wells.
- 10 Wash each well 3 times with  2 mL ice cold 1X PBS (filtered).
- 11 Add M-PER with 1X protease and phosphatase inhibitor to each well (*see specific protocols for volume*)  
  
*Be sure to keep M-PER and samples on ice.*  
  
Note: HALT protease and phosphatase inhibitor *must* be diluted in lysis buffer *immediately* before use (cannot be prepared the day prior).
- 12 Scrape the cells from each well 1 at a time using a cell scraper. Pipette all of the M-PER containing the cell lysate from each well into a clean 1.5 mL micro-centrifuge tube and immediately place the tubes on ice.
- 13 Centrifuge the tubes for  00:10:00 at  4 °C at  10000 rcf .
- 14 Collect  100 µL aliquots of the supernatant and store them for future use (lysate that will be used within 1-2 weeks after harvesting can be stored at  -20 °C , otherwise store lysates in  -80 °C , per **Cell Signaling** recommendations).

## Protocol references

Promega FuGENE® 4K Transfection Reagent (#E5911) Technical Manual:

[https://www.promega.com/-/media/files/resources/protocols/technical-manuals/500/fugene-4k-transfection-reagent-protocol-tm694.pdf?rev=beb16dbcf2604534ace3cbe95b5080cb&sc\\_lang=en](https://www.promega.com/-/media/files/resources/protocols/technical-manuals/500/fugene-4k-transfection-reagent-protocol-tm694.pdf?rev=beb16dbcf2604534ace3cbe95b5080cb&sc_lang=en)

HALT Protease and Phosphatase Inhibitor:

[https://www.thermofisher.com/document-connect/document-connect.html?url=https://assets.thermofisher.com/TFS-Assets%2FLSG%2Fmanuals%2FMAN0011626\\_HaltProteasePhosphat\\_Inhibit\\_UG.pdf](https://www.thermofisher.com/document-connect/document-connect.html?url=https://assets.thermofisher.com/TFS-Assets%2FLSG%2Fmanuals%2FMAN0011626_HaltProteasePhosphat_Inhibit_UG.pdf)