




Sep 17, 2020

# PMN- 06 - Culture of Human PMN - TNF- $\alpha$ production

Marco Cosentino<sup>1</sup>, Elisa Storelli<sup>1</sup>, Alessandra Luini<sup>1</sup>, Massimiliano LM Legnaro<sup>1</sup>, Emanuela Rasini<sup>1</sup>, Marco Ferrari<sup>1</sup>, Franca Marino<sup>1</sup>

<sup>1</sup>Center for Research in Medical Pharmacology, University of Insubria (Varese, Italy)

1 Works for me dx.doi.org/10.17504/protocols.io.bkrakv2e

 Farmacologia Medica

## ABSTRACT

Published work using this protocol:

**- A Novel Standardized Cannabis sativa L. Extract and Its Constituent Cannabidiol Inhibit Human Polymorphonuclear Leukocyte Functions.** Mabou Tagne A, Marino F, Legnaro M, Luini A, Pacchetti B, Cosentino M. Int J Mol Sci. 2019 Apr 13;20(8):1833. doi: **10.3390/ijms20081833**

DOI

[dx.doi.org/10.17504/protocols.io.bkrakv2e](https://dx.doi.org/10.17504/protocols.io.bkrakv2e)

## PROTOCOL CITATION

Marco Cosentino, Elisa Storelli, Alessandra Luini, Massimiliano LM Legnaro, Emanuela Rasini, Marco Ferrari, Franca Marino 2020. PMN- 06 - Culture of Human PMN - TNF- $\alpha$  production. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.bkrakv2e>

## LICENSE

———— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Sep 03, 2020

## LAST MODIFIED

Sep 17, 2020

## PROTOCOL INTEGER ID

41474

## MATERIALS

NAME	CATALOG #	VENDOR
Fetal Bovine Serum (FBS)	ECS0180L-500 ml	EuroClone
RPMI 1640	ECM 0495L- 500 ml	EuroClone
Lipopolysaccharides from Escherichia coli O127:B8	L3137	Sigma
Penicillin/Streptomycin	ECB3001D - 100 ml	EuroClone
ELISA kit for TNF- $\alpha$	KHC0031	Thermo Fisher Scientific
Cannabidiol	not available	Linnea SA

## MATERIALS TEXT






Instrumentation needed:


Sterile and non sterile plastic disposables

ELISA plate reader

EQUIPMENT

NAME	CATALOG #	VENDOR
Allegra AVANTI 30	Beckman Italy	Beckman Coulter

- 1 Isolate PMN according to the protocol PMN-01a or PMN-01b.
  - 2 Resuspend PMN at  $10 \times 10^6$  cells/ml in RPMI medium supplemented with 10% FBS and 1% penicillin/streptomycin.
  - 3 Add  **300 µl** of cell suspension in a 5 ml tube for cell culture and incubate for  **21:00:00** at  **37 °C** in 5% CO<sub>2</sub>.  
Cells were incubated in resting condition or in presence of stimulus. Stimulus (**LPS 1 µg/ml** and/or **CBD 10<sup>-8</sup>-10<sup>-5</sup> M**) in 5 ml tubes for cell culture. **Other possible activating stimuli** can be selected and included.
  - 4 Centrifuge at  **1500 x g, 00:10:00**
- 

Allegra AVANTI 30  
Centrifuge  
Beckman Coulter   Beckman Italy
- 5 Harvest cell-free supernatant (  **150 µl** x2 aliquots) and cell pellet.
  - 6 Store the supernatant at -80°C until the TNF-α Enzyme Linked-Immuno-Sorbent Assay (ELISA) is performed.
  - 7 Store cell pellet at -80°C for gene expression analysis (RT PCR).