



JUN 21, 2023

🌐 LB Broth

sgoodwin^{1,2}

¹University of Sheffield; ²Sheffield iGEM



sgoodwin

ABSTRACT

Process for producing LB broth growth medium

OPEN ACCESS

DOI:

dx.doi.org/10.17504/protocols.io.yxmvm375nl3p/v1

Protocol Citation: sgoodwin 2023. LB Broth. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.yxmvm375nl3p/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development
We are still developing and optimizing this protocol



Created: Jun 21, 2023

Last Modified: Jun 21, 2023

PROTOCOL integer ID:
83812

Making the LB broth

3h 10m

1 Dilute  20 g of LB powder in  800 mL distilled water

5m

2 Stir until fully dissolved with magnetic stirrer

5m

3 Transfer to a Duran flask, and sterilise by autoclaving

3h