



Aug 14, 2020

Sandwich ELISA for investigating the binding of Protein-LG (SpLG) to avian immunoglobulins using anti-IgY-peroxidase as conjugate.

Angel A Justiz-Vaillant¹, Monica F. Smikle²

¹University of the West Indies St. Augustine; ²University of the West Indies. Mona Campus

1 Works for me dx.doi.org/10.17504/protocols.io.bjq4kmyw

University of the West Indies angel.vaillant@sta.uwi.edu

Angel Justiz-Vaillant
University of the West Indies St. Augustine

DOI

dx.doi.org/10.17504/protocols.io.bjq4kmyw

PROTOCOL CITATION

Angel A Justiz-Vaillant, Monica F. Smikle 2020. Sandwich ELISA for investigating the binding of Protein-LG (SpLG) to avian immunoglobulins using anti-IgY-peroxidase as conjugate.. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bjq4kmyw>

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 14, 2020

LAST MODIFIED

Aug 14, 2020

PROTOCOL INTEGER ID

40444

MATERIALS

NAME	CATALOG #	VENDOR
Anti-Chicken IgY, HRP Conjugate, 300ul	G1351	Promega
Nunc&trade; 96-Well Polystyrene Round Bottom Microwell Plates, V 96 well plate, Non-Treated, clear, without lid, Sterile	260210	Thermo Fisher
Protein-L from P. Magnus		
Streptococcal protein G by Sigma Aldrich		

- 1 This ELISA is used to study the interaction of protein-LG (SpLG) with diverse avian immunoglobulins.
- 2 The 96 well microtitre plate is coated overnight at 4°C with 2 µg/µl per well of recombinant SpLG or a mixture of SpL with SpG in carbonate-bicarbonate buffer pH 9.6.

- 3 Then plate is treated with bovine serum albumin solution and washed 4X with PBS-Tween.
- 4 50 µl of avian egg yolk or egg white (1 mg/ml) is added and incubated for 1.30h at RT and the microplate is then rewashed 4X with PBS-Tween.
- 5 Then 50 µl of peroxidase-labeled-anti-IgY conjugate diluted 1:15000 in PBS-non-fat milk is added to each well and incubated for 1.30h at RT. After that the plate is washed 4X with PBS-Tween.
- 6 Pipette 50 µl of 3,3',5,5' - tetramethylbenzidine (TMB; Sigma-Aldrich) to each well.
- 7 The reaction is stopped with 50 µl of 3M H₂SO₄ solution.
- 8 The plate is visually assessed for the development of colour and read in a microplate reader at 450 nm.
- 9 A cut-off point can be calculated as the mean of the optical density of negative controls x 3. The higher the OD value the higher will be the binding affinity of SpLG to avian immunoglobulins.