

Sep 30, 2020

in situ hybridization

Forked from [FISH and antibody staining](#)

1

1

1 Works for me

This protocol is published without a DOI.

Zhang Xue



1

PROTOCOL CITATION

2020. in situ hybridization. [protocols.io](https://protocols.io/view/in-situ-hybridization-bmvzk676)
<https://protocols.io/view/in-situ-hybridization-bmvzk676>

FORK FROM

Forked from [FISH and antibody staining](#), 1

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Sep 30, 2020

LAST MODIFIED

Sep 30, 2020

PROTOCOL INTEGER ID

42649

BEFORE STARTING

1. 1.5ml EP egg water 4% PFA 2. PBT 2 3. 1ml MetOH shake for 5min 3 4. -20°C 1

Day1

1 75%MetOH+25%PBT 1ml shake for 5min 00:05:00

2 50%MetOH+50%PBT 5min 00:05:00

3 25%MetOH+75%PBT 5min 00:05:00

4 100%PBT 5min 1/4 00:05:00



PK


5 100%PBT 5min 2/4 ⌚ 00:05:00

6 100%PBT 5min 3/4 ⌚ 00:05:00

7 100%PBT 5min 4/4 ⌚ 00:05:00

8 100%Proteinase K: PBT = 1 : 1000 (10µg/ml) 24h 5min 36h 10min 48h 15min 60h 20min 72h 30min
⌚ 00:00:00 ⌚ 00:10:00 ⌚ 00:15:00 ⌚ 00:20:00 ⌚ 00:30:00

9 4%PFA 30min

 PFA 4% 30min

⌚ 00:30:00

10 1xPBT 5min 1/4 ⌚ 00:05:00

11 1xPBT 5min 2/4 ⌚ 00:05:00

12 1xPBT 5min 3/4 ⌚ 00:05:00


13 1xPBT 5min 4/4 ⌚ 00:05:00

14 500µl Hybe buffer 68.5°C 2-5h 🔥 68.5 °C ⌚ 03:00:00

15 6µl probe+200µl Hybe buffer 68.5°C overnight 🔥 68.5 °C

Day2

16 -20°C 700-800µl 750 µl 🔥 68.5 °C

 68.5°C 29

17 Hybe buffer 10min ⌚ 00:10:00

18 75% Hybe buffer+25% 2×SSCT 10min ⌚ 00:10:00

19 50%Hybe buffer+50% 2×SSCT 10min ⌚ 00:10:00

20 25%Hybe buffer+75% 2×SSCT 10min ⌚ 00:10:00

21 2×SSCT 10min ⌚ 00:10:00

22 0.2×SSCT 15min 1/4 ⌚ 00:15:00

23 0.2×SSCT 15min 2/4 ⌚ 00:15:00

24 0.2×SSCT 15min 3/4 ⌚ 00:15:00

25 0.2×SSCT 15min 4/4 ⌚ 00:15:00



XXXXXXXXXXXXXXXX

26 75% 0.2× SSCT+25% MABt 5min ⌚ 00:05:00 ⚡ Room temperature 🧴 1 mL

27 50% 0.2× SSCT+50%MABt 5min ⌚ 00:05:00

28 25%0.2× SSCT+75%MABt 5min ⌚ 00:05:00

29 MABt 5min ⌚ 00:05:00

30 1×blocking 500μl 700μl 2h RT ⌚ 02:00:00

🔧 Room temperature 🧴 500 μl

31 Anti-Dig-AP 1:2000 200μl 4°C 4 °C 🧴 200 μl

 FISH Anti-Dig-AP

Day3

32 4°C

33 1×MABt 15min 1/8 ⌚ 00:15:00

34 1×MABt 15min 2/8 ⌚ 00:15:00

35 1×MABt 15min 3/8 ⌚ 00:15:00

36 1×MABt 15min 4/8 ⌚ 00:15:00

37 1×MABt 15min 5/8 ⌚ 00:15:00

38 1×MABt 15min 6/8 ⌚ 00:15:00

39 1×MABt 15min 7/8 ⌚ 00:15:00

40 1×MABt 15min 8/8 NTMT 3.5ml 00:15:00

20ml	NTMT
2ml	1M pH9.5 Tris-HCl
2ml	1M NaCl
0.5ml	10% Tween
15.5ml	ddH2O

41 Equilibrate with NTMT 5min 1/3 00:05:00

42 24 NTMT 5min 2/3 00:05:00

BCIP in NTMT 1:50 500µl

43 NTMT 5min 3/3 00:05:00

44 BCIP in NTMT 1:50 500µl 37°C 15min BF 500 µl 37 °C 00:30:00

45 1ml stop solution