



Dec 11, 2020

Antigen Presentation Protocol

Lilia Rodriguez¹¹Université de Montréal

1

Works for me

dx.doi.org/10.17504/protocols.io.bpk4mkyw

Lilia Rodriguez

ABSTRACT

This protocol details methods for 3-day Antigen Presentation Assay.

DOI

dx.doi.org/10.17504/protocols.io.bpk4mkyw

PROTOCOL CITATION

Lilia Rodriguez 2020. Antigen Presentation Protocol. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bpk4mkyw>



KEYWORDS

antigen, presentation, Ag

LICENSE

———— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Nov 11, 2020

LAST MODIFIED

Dec 11, 2020

OWNERSHIP HISTORY

Nov 11, 2020 Liz Brydon Protocols.io

Dec 04, 2020 Lilia Rodriguez

PROTOCOL INTEGER ID

44412

MATERIALS TEXT

REAGENTS

- LPS/IFN γ (including controls)
- RPMI media
- PFA
- PBS
- Wash buffer: DMEM with 0.1M glycine + 10% iFBS
- Stock peptide
- 1x lysis buffer (stock = 5x, with triton, pH = 7.8) in dH₂O
- 1M DTT
- CPRG buffer
- Water
- CPRG

CONSUMABLES

- 96-well plates
- 50ml conical tubes

EQUIPMENT

- Incubator: 37°C with 5% CO₂
- Centrifuge
- Cell counter
- Plate reader

SAFETY WARNINGS

Please refer to Safety Data Sheets (SDS) for health and environmental hazards.

ABSTRACT

This protocol details methods for 3-day Antigen Presentation Assay.

Day 1 1d

- 1 Prepare a 96 well plate:
 - Label appropriately (experimental condition – LPS/IFN γ controls, infection controls, etc.)
 - Duplicate labelled wells for the peptide control

2 

Scrape RAW cells and resuspend (pipette up and down 8x).

- 3 Count RAW cells and adjust concentration to **0.75 million cells/ml**.

4 

Add  **200 μ l RAW cells** to each well.

- 5 Add LPS/IFN γ accordingly.

6 

1d

Incubate for ⌚ **24:00:00** at 🌡 **37 °C** with **5% CO2** .

Day 2 10m

7 Prior to beginning, check the cell's media (red/orange = good, yellow = bad).

8 Prepare **1 % (v/v) PFA in PBS** at 🌡 **Room temperature** (**pH7.4**).

9 

Discard supernatant in the 96 well plate and add **50 µl 1% PFA** .

9.1 

10m

Incubate at 🌡 **Room temperature** for ⌚ **00:10:00** .

10 Prepare 2E2 cells while RAW cells are in 1% PFA:

10.1 

Pour flask into 50 ml conical and centrifuge at 🌀 **1500 rpm, 00:03:00** .

To keep more 2E2 cells growing, add **50 mL RPMI media** back into their flask and incubate at 🌡 **37 °C** , **5% CO2** .

10.2 Count 2E2 cells and adjust concentration to **0.4 million cells/ml**.

11 

Add **200 µl wash buffer** to each well in the 96 well plate and discard immediately.

Wash buffer = DMEM with 0.1M glycine, + 10% iFBS.

11.1

Repeat 2x more: Add **200 µl wash buffer** to each well in the 96 well plate and discard immediately. **(Wash 1/2)**

11.2

Repeat 1x more: Add **200 µl wash buffer** to each well in the 96 well plate and discard immediately. **(Wash 2/2)**

12

Add **250 µl 2E2 cells** per well in ½ of the samples (the non-peptide ones). Add an extra row of 2E2 cells alone for a negative control.

13 Dilute stock peptide (1 µg/ml) **5000 fold** in the 2E2 cells (final concentration: **0.2 nanogram per milliliter (ng/mL)**).

14

Add **250 µl peptide + 2E2 cell solution** to the remaining wells.

15

16h

Incubate at **37 °C** , **5% CO2** for **16:00:00 (NO LONGER)** .

Day 3 5m

16 Prepare 1x lysis buffer (stock = 5x, with triton, **pH7.8**) in dH₂O.

16.1

Add **30 µl 1M stock of DTT** in **10 mL lysis buffer** .

17 Prepare CPRG buffer (**pH7.8**).

18

Centrifuge 96 well plate at **2200 rpm, 00:01:00** .

19 

Add  **50 µl lysis buffer** to all wells.




19.1 

5m

Incubate at  **Room temperature** for  **00:05:00** (up to 20 minutes MAX).

20 


Prepare CPRG (recipe below = for each well)

-  **150 µl CPRG buffer**
-  **20.2 µl water**
-  **0.046 mg CPRG**

21 

When lysis is done, add  **170 µl CPRG solution** / well.

21.1 

Incubate either at  **Room temperature** or 37°C (37 speeds up reaction to about ~ 20 minutes for peptide samples).

22 

Transfer  **150 µl colored solution** to a new plate.

Take care not to transfer debris or make bubbles.

23 

Take reading at 595 nm or 570 nm.

* To stop reaction to leave overnight, incubate at  **4 °C** .