



Feb 07, 2021

# Action Spectra protocol for Opentrons OT-1 liquid handling robot

andrei.herdean <sup>1</sup><sup>1</sup>University of Technology Sydney

andrei.herdean : Climate Change Cluster

**1** Works for me [dx.doi.org/10.17504/protocols.io.br6vm9e6](https://dx.doi.org/10.17504/protocols.io.br6vm9e6)

Climate Change Cluster

 andrei.herdean

DOI

[dx.doi.org/10.17504/protocols.io.br6vm9e6](https://dx.doi.org/10.17504/protocols.io.br6vm9e6)

PROTOCOL CITATION

andrei.herdean 2021. Action Spectra protocol for Opentrons OT-1 liquid handling robot . **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.br6vm9e6>



LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Feb 07, 2021

LAST MODIFIED

Feb 07, 2021

PROTOCOL INTEGER ID

47029

## 1 OT-1 setup:

A3 - box with 1 ml tips

B3 - 48 well plate 1

D3 - 48 well plate 2

C1 - 48 well plate 3

48 plate 1 setup:

C8 - oxygen optode

D8 - temperature sensor

A7, B7, C7, D7, E7, F7 - blank media

F4, F5, F6 - ddH<sub>2</sub>O

A8, B8, E8, F8 - empty

All other wells have 1 ml of algae culture

48 plate 2 setup:

All wells empty, except:

F2 - temperature sensor

F1 - oxygen optode

48 well plate 3 setup:

All wells empty, except:

C8 - oxygen optodes\*

B8, D8 - temperature sensors

\* a fiber optic cable connected to the OL490 is positioned in front of well C8

## 2 OT-1 protocol:

```
from opentrons import robot, containers, instruments
robot.home()
#containers to load
plateC1 = containers.load('48-well-plate', 'C1')
change_light = containers.load('6-well-plate', 'A1')
plateD3 = containers.load('48-well-plate', 'D3')
plateB3 = containers.load('48-well-plate', 'B3')
tiprack_1 = containers.load('tiprack-1000ul', 'A3')
#optode well definition
dark_optode = 'F1'
light_optode = 'C8'

#pipette definition
p1000 = instruments.Pipette(
    axis="b",
    max_volume=1000,
    tip_racks=[tiprack_1])

#how much additional distance to raise the tip
robot.arc_height = 20

#starting tip position
p1000.reset()
p1000.start_at_tip(tiprack_1['A1'])

#Test trigger
p1000.touch_tip(change_light.wells('A1'),radius=1.7) # trigger light
p1000.touch_tip(change_light.wells('A1'),radius=1.7) # trigger light
p1000.touch_tip(change_light.wells('A1'),radius=1.7) # trigger light
p1000.touch_tip(change_light.wells('A1'),radius=1.7) # trigger light

#START
#transfer blank samples (water to optodes)
p1000.pick_up_tip()
p1000.aspirate(1000, plateB3.wells('A7'))
p1000.dispense(1000, plateB3.wells(light_optode))
p1000.blow_out()

p1000.aspirate(1000, plateB3.wells('B7'))
p1000.dispense(1000, plateC1.wells(light_optode))
p1000.blow_out()

p1000.aspirate(1000, plateB3.wells('C7'))
p1000.dispense(1000, plateD3.wells(dark_optode))
p1000.blow_out()
p1000.delay(seconds=600)

#pre-loop
#step 1 - empty All light optode
```

```

p1000.mix(3, 1000, plateB3.wells(light_optode))
p1000.aspirate(1000, plateB3.wells(light_optode))
p1000.dispense(1000, plateB3.wells('A7'))
p1000.blow_out()
#step 2 - sample1 in All light optode
p1000.mix(3, 1000, plateB3.wells('A1'))
p1000.aspirate(1000, plateB3.wells('A1'))
p1000.dispense(1000, plateB3.wells(light_optode))
p1000.blow_out()
p1000.delay(seconds=600)

#step 3 - empty OL490
p1000.mix(3, 1000, plateC1.wells(light_optode))
p1000.aspirate(1000, plateC1.wells(light_optode))
p1000.dispense(1000, plateB3.wells('B7'))
p1000.blow_out()
p1000.touch_tip(change_light.wells('A1'),radius=1.7) # trigger light

#step 4 - sample in OL 490
p1000.mix(3, 1000, plateB3.wells(light_optode))
p1000.aspirate(1000, plateB3.wells(light_optode))
p1000.dispense(1000, plateC1.wells(light_optode))
p1000.blow_out()
#Step 5 - sample 2 in All light optode
p1000.mix(3, 1000, plateB3.wells('A2'))
p1000.aspirate(1000, plateB3.wells('A2'))
p1000.dispense(1000, plateB3.wells(light_optode))
p1000.blow_out()
p1000.delay(seconds=600)

#Step 6 - empty Dark optode
p1000.mix(3, 1000, plateD3.wells(dark_optode))
p1000.aspirate(1000, plateD3.wells(dark_optode))
p1000.dispense(1000, plateB3.wells('C7'))
p1000.blow_out()
#Step 7 - Sample 1 in dark optode
p1000.mix(3, 1000, plateC1.wells(light_optode))
p1000.aspirate(1000, plateC1.wells(light_optode))
p1000.dispense(1000, plateD3.wells(dark_optode))
p1000.blow_out()
p1000.touch_tip(change_light.wells('A1'),radius=1.7) # trigger light

#Step 8 - Sample 2 in OL490
p1000.mix(3, 1000, plateB3.wells(light_optode))
p1000.aspirate(1000, plateB3.wells(light_optode))
p1000.dispense(1000, plateC1.wells(light_optode))
p1000.blow_out()
#Step 9 - Sample 3 in All light optode
p1000.mix(3, 1000, plateB3.wells('A3'))
p1000.aspirate(1000, plateB3.wells('A3'))
p1000.dispense(1000, plateB3.wells(light_optode))
p1000.blow_out()
p1000.delay(seconds=600)

#Loop
def run_loop(sample_well1 = 'A4', sample_well2 = 'A1'): #00:00
    #Step 1 - empty DARK optode
    p1000.mix(3, 1000, plateD3.wells(dark_optode))

```

```

p1000.aspirate(1000, plateD3.wells(dark_optode))
p1000.dispense(1000, plateB3.wells(sample_well2))
p1000.blow_out() #00:54
#Step 2 - OL490 to Dark
p1000.mix(3, 1000, plateC1.wells(light_optode)) #01:22
p1000.aspirate(1000, plateC1.wells(light_optode)) #01:24
p1000.dispense(1000, plateD3.wells(dark_optode)) #01:31
p1000.blow_out() #01:47
p1000.touch_tip(change_light.wells('A1'),radius=1.7) # trigger light

#Step 3 - ALL light to OL490
p1000.mix(3, 1000, plateB3.wells(light_optode))
p1000.aspirate(1000, plateB3.wells(light_optode))
p1000.dispense(1000, plateC1.wells(light_optode))
p1000.blow_out() #02:41
#Step 4 - new sample in ALL light
p1000.mix(3, 1000, plateB3.wells(sample_well1))
p1000.aspirate(1000, plateB3.wells(sample_well1))
p1000.dispense(1000, plateB3.wells(light_optode))
p1000.blow_out() #11:26
p1000.delay(seconds=600)

sample_wells = [name+number for name in 'ABCDEF' for number in '123456']

for i in range(3, len(sample_wells)):
    run_loop(sample_wells[i], sample_wells[i-3])

#Wash procedure
p1000.mix(3, 1000, plateB3.wells(light_optode))
p1000.aspirate(1000, plateB3.wells(light_optode))
p1000.dispense(1000, plateB3.wells('F4'))
p1000.blow_out()
p1000.return_tip()
p1000.pick_up_tip()

p1000.aspirate(1000, plateB3.wells('A7'))
p1000.dispense(1000, plateB3.wells(light_optode))
p1000.blow_out()
p1000.mix(3, 1000, plateB3.wells(light_optode))
p1000.aspirate(1000, plateB3.wells(light_optode))
p1000.dispense(1000, plateB3.wells('A7'))
p1000.blow_out()
p1000.return_tip()
p1000.pick_up_tip()

p1000.aspirate(1000, plateB3.wells('D7'))
p1000.dispense(1000, plateB3.wells(light_optode))
p1000.blow_out()

p1000.mix(3, 1000, plateC1.wells(light_optode))
p1000.aspirate(1000, plateC1.wells(light_optode))
p1000.dispense(1000, plateB3.wells('F5'))
p1000.blow_out()
p1000.return_tip()
p1000.pick_up_tip()

p1000.aspirate(1000, plateB3.wells('B7'))
p1000.dispense(1000, plateC1.wells(light_optode))

```

```

p1000.blow_out()
p1000.mix(3, 1000, plateC1.wells(light_optode))
p1000.aspirate(1000, plateC1.wells(light_optode))
p1000.dispense(1000, plateB3.wells('B7'))
p1000.blow_out()
p1000.return_tip()
p1000.pick_up_tip()

p1000.aspirate(1000, plateB3.wells('E7'))
p1000.dispense(1000, plateC1.wells(light_optode))
p1000.blow_out()
#
p1000.mix(3, 1000, plateD3.wells(dark_optode))
p1000.aspirate(1000, plateD3.wells(dark_optode))
p1000.dispense(1000, plateB3.wells('F6'))
p1000.blow_out()
p1000.return_tip()
p1000.pick_up_tip()

p1000.aspirate(1000, plateB3.wells('C7'))
p1000.dispense(1000, plateD3.wells(dark_optode))
p1000.blow_out()
p1000.mix(3, 1000, plateD3.wells(dark_optode))
p1000.aspirate(1000, plateD3.wells(dark_optode))
p1000.dispense(1000, plateB3.wells('C7'))
p1000.blow_out()
p1000.return_tip()
p1000.pick_up_tip()

p1000.aspirate(1000, plateB3.wells('F7'))
p1000.dispense(1000, plateD3.wells(dark_optode))
p1000.blow_out()
p1000.delay(seconds=600)

p1000.aspirate(1000, plateB3.wells(light_optode))
p1000.dispense(1000, plateB3.wells('D7'))
p1000.blow_out()

p1000.aspirate(1000, plateC1.wells(light_optode))
p1000.dispense(1000, plateB3.wells('E7'))
p1000.blow_out()

p1000.aspirate(1000, plateD3.wells(dark_optode))
p1000.dispense(1000, plateB3.wells('F7'))
p1000.blow_out()

p1000.return_tip()
robot.home()

```