



Preparation of a protein-AG conjugated to horseradish peroxidase by the periodate method.

Angel A Justiz-Vaillant¹

¹University of the West Indies St. Augustine

1 Works for me dx.doi.org/10.17504/protocols.io.bjk6kkze

University of the West Indies angel.vaillant@sta.uwi.edu

Angel Justiz-Vaillant
University of the West Indies St. Augustine

ABSTRACT

A recombinant protein that combines the IgG-binding domains of SpA and SpG was developed, and labelled to horseradish peroxidase. It was used as universal conjugate in ELISA for the assessment of antibodies against *Brucella* spp in cattle, sheep, dogs, goats and pigs. It was reported that similar results as the one shown using the chimeric protein AG were obtained when murine monoclonal antibody-enzyme conjugates were used [1,2].

References

1. Justiz-Vaillant AA, Akpaka PE, McFarlane-Anderson N, Smikle MF. Comparison of techniques of detecting immunoglobulin-binding protein reactivity to immunoglobulin produced by different avian and mammalian species. *West Indian Med J.* 2013;62(1):12-20.
2. Justiz-Vaillant AA, McFarlane-Anderson N, and Smikle M. "Bacterial Immunoglobulin (Ig)-Receptors: Past and Present Perspectives." *American Journal of Microbiological Research*, vol. 5, no. 2 (2017): 44-50. doi: 10.12691/ajmr-5-2-4.

DOI

dx.doi.org/10.17504/protocols.io.bjk6kkze

PROTOCOL CITATION

Angel A Justiz-Vaillant 2020. Preparation of a protein-AG conjugated to horseradish peroxidase by the periodate method.. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bjk6kkze>

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 11, 2020

LAST MODIFIED

Aug 11, 2020

PROTOCOL INTEGER ID

40318

GUIDELINES

All reagents but specially the enzyme and the sodium periodate solution have to be prepared freshly before mixing it with the enzyme.

MATERIALS

NAME	CATALOG #	VENDOR
Ammonium Sulfate		P212121
Sodium periodate	SB0875.SIZE.100g	Bio Basic Inc.
sodium borohydride	452882	Sigma Aldrich
Horseradish Peroxidase (HRP) type IV	P8375-25KU	Sigma Aldrich
Staphylococcal Protein-A		Sigma Aldrich
Streptococcal protein G by Sigma Aldrich		

MATERIALS TEXT

Pipettes
20ml to 1000 ml glass
Scale
Incubator
Refrigerator
Freezer
Centrifuges

- 1 Horseradish peroxidase (500 µg in 50 µl NaCO₃, pH 9.6) is mixed with freshly made sodium periodate solution (1.71 mg/ml) followed by incubation in the dark for 2 h.
- 2 Mix 500 µg of staphylococcal protein-A (SpA) with an equal amount (500 micrograms) of a mix of horseradish peroxidase-sodium periodate. On the other hand mix 500 µg of streptococcal protein-G (SpG) with an equal amount (500 micrograms) of the mix of horseradish peroxidase-sodium periodate.
- 3 The two mixtures are incubated separately for 3 hours at 4°C with gentle agitation.
- 4 Forty µl of freshly prepared NaBH₄ solution (5 mg NaBH₄ /ml 0.1 mM NaOH) is then added separately to the preparations, which are centrifuge (13,000rpm., 10 minutes at RT). Add to each preparation cold saturated ammonium sulphate solution and centrifuge again (10000rpm, 25 minutes at 4°C).
- 5 Now mix the SpA-HRP preparation with SpG-HRP and incubate the mixture for 90 min at 4°C in the dark with gentle agitation.
- 6 The mixture is then centrifuged for 25 min at 4°C and recover the pellet at the bottom of the tube.
- 7 The pellet (SpAG-HRP) is re-suspended in 500 µl of PBS pH=7.4 and dialysed against 1L of PBS for 24 h with 3 buffer changes.
- 8 An equal volume of glycerol is added to the dialysate followed by 200 µl of bovine serum albumin, BSA (20 mg/ ml).

- 9 The SpAG-HRP conjugate is then stored at -20°C until further used.