



Aug 17, 2021

VPS13D DNA plasmid generation

Andrés Guillén-Samander^{1,2}, Marianna Leonzino^{1,2}, Pietro De Camilli^{1,2}

¹Departments of Neuroscience and of Cell Biology, Howard Hughes Medical Institute, Program in Cellular Neuroscience, Neurodegeneration and Repair, Yale University School of Medicine, New Haven, Connecticut 06510, USA;

²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815

1 Works for me



Share

dx.doi.org/10.17504/protocols.io.bvgwn3xe

Andrés Guillén-Samander

ABSTRACT

This protocol describes the basic molecular cloning technique utilized for the generation of VPS13D constructs in <https://doi.org/10.1083/jcb.202010004>. This protocol and the enzymes included in it are commercialized by Takara Bio.

Due to low expression of big DNA constructs (>10kb), we decided to generate a codon-optimized cDNA encoding for human VPS13D, including an mScarlet fluorescent protein after aminoacid 1576 flanked by BamHI restriction enzyme sites. The construct was generated by and purchased from Genscript. From this initial construct, we generated several VPS13D constructs. The specific primers and enzymes used for each construct are included in Table 1 of our manuscript: <https://doi.org/10.1083/jcb.202010004>. For most of our cloning procedures, Infusion cloning (Takara) was used.

ATTACHMENTS

[dn3ybgtx.pdf](#)

DOI

dx.doi.org/10.17504/protocols.io.bvgwn3xe

PROTOCOL CITATION

Andrés Guillén-Samander, Marianna Leonzino, Pietro De Camilli 2021. VPS13D DNA plasmid generation.
protocols.io
<https://dx.doi.org/10.17504/protocols.io.bvgwn3xe>

KEYWORDS

VPS13D, DNA, Plasmid generation, ASAPCRN

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jun 02, 2021

LAST MODIFIED

Jul 02, 2021

OWNERSHIP HISTORY

Jun 02, 2021 Urmilas

Jun 08, 2021 Andrés Guillén-Samander

VPS13D DNA plasmid generation

- 1 Linearize the backbone with restriction enzymes or by PCR amplification.



Amplify the insert by PCR with primers including 15nt of overlapping sequence with the target backbone.

Suggestion: Takara has an online tool to help with primer design at: <https://go.shr.lc/3daJf0E>

- 3 Run both products in an agarose gel to confirm expected size and purify the DNA from gel using a NucleoSpin Gel and PCR Clean-up kit (Takara).
- 4 Ligate the linearized backbone and the insert with the Infusion enzyme mix (Takara).

Suggestion: We found that big DNA constructs accumulate mutations easily when amplified in bacteria, this can be minimized by growing bacteria at **30 °C**.