



Aug 13, 2020

Universal sandwich ELISA for investigating the binding of avian immunoglobulins to Staphylococcal protein-A (SpA) using anti-IgY-peroxidase conjugate.

Angel A Justiz-Vaillant¹, Monica F. Smikle²

¹University of the West Indies St. Augustine; ²University of the West Indies. Mona Campus

1 Works for me dx.doi.org/10.17504/protocols.io.bjqxkmxn

University of the West Indies angel.vaillant@sta.uwi.edu

Angel Justiz-Vaillant
University of the West Indies St. Augustine

DOI

dx.doi.org/10.17504/protocols.io.bjqxkmxn

PROTOCOL CITATION

Angel A Justiz-Vaillant, Monica F. Smikle 2020. Universal sandwich ELISA for investigating the binding of avian immunoglobulins to Staphylococcal protein-A (SpA) using anti-IgY-peroxidase conjugate..

protocols.io

<https://dx.doi.org/10.17504/protocols.io.bjqxkmxn>

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 13, 2020

LAST MODIFIED

Aug 13, 2020

PROTOCOL INTEGER ID

40439

MATERIALS

| NAME | CATALOG # | VENDOR |
|--|-----------|---------------|
| Anti-Chicken IgY, HRP Conjugate, 300ul | G1351 | Promega |
| Nunc® 96-Well Polystyrene Round Bottom Microwell Plates, V 96 well plate, Non-Treated, clear, without lid, Sterile | 260210 | Thermo Fisher |
| Staphylococcal Protein-A | | Sigma Aldrich |

- 1 This ELISA is used to study the interaction of Staphylococcal protein-A (SpA) with diverse avian immunoglobulins.
- 2 The 96 well microtitre plate is coated overnight at 4°C with 2 µg/µl per well of SpA in carbonate-bicarbonate buffer pH 9.6.

- 3 Then plate is treated with bovine serum albumin solution and washed 4X with PBS-Tween.
- 4 50 µl of avian egg yolk or egg white (1 mg/ml) is added and incubated for 1.30h at room temperature and the microplate is then rewashed 4X with PBS-Tween.
- 5 Then 50 µl of peroxidase-labeled-anti-IgY conjugate diluted 1:15000 in PBS-non-fat milk is added to each well and incubated for 1.30h at RT. The plate is washed 4X with PBS-Tween.
- 6 Pipette 50 µl of 3,3',5,5' - tetramethylbenzidine (TMB; Sigma-Aldrich) to each well.
- 7 The reaction is stopped with 50 µl of 3M H₂SO₄ solution.
- 8 The plate is visually assessed for the development of colour and read in a microplate reader at 450 nm.
- 9 A cut-off point should be calculated as the mean of the optical density of negative controls x 3. The higher the OD value the higher will be the affinity of SpA to avian immunoglobulins.