

6



Dec 21, 2020

# © nCoV-2019 McGill Artic PCR Protocol, 5 ul RT and V3 only at 63C

Forked from nCoV-2019 McGill Artic PCR Protocol, 5 ul RT and V3 only + LA1 at 63C

Sarah J Reiling<sup>1</sup>, Anne-Marie Roy<sup>1</sup>, Shu-Huang Chen<sup>1</sup>, Ioannis Ragoussis<sup>1</sup>

<sup>1</sup>McGill University

1 Works for me

dx.doi.org/10.17504/protocols.io.bj9hkr36

McGill Genome Centre

Anne-Marie Roy

**ABSTRACT** 

V3 only primers can be found here:

https://github.com/sarahreiling/artic-ncov2019/blob/master/primer\_schemes/nCoV-2019/V3/nCoV-2019\_V3only.scheme.bed

DOL

dx.doi.org/10.17504/protocols.io.bj9hkr36

PROTOCOL CITATION

Sarah J Reiling, Anne-Marie Roy, Shu-Huang Chen, Ioannis Ragoussis 2020. nCoV-2019 McGill Artic PCR Protocol, 5 ul RT and V3 only at 63C. **protocols.io** 

https://dx.doi.org/10.17504/protocols.io.bj9hkr36

FORK NOTE

FORK FROM

Forked from nCoV-2019 McGill Artic PCR Protocol, 5 ul RT and V3 only + LA1 at 63C, Sarah Reiling

LICENSE

This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 24, 2020

LAST MODIFIED

Dec 21, 2020

PROTOCOL INTEGER ID

40969

 MATERIALS TEXT
MATERIALS

Biolabs Catalog #M0492L

⋈ nuclease-free water Contributed by users

**⊠** Quant-iT™ PicoGreen™ dsDNA Assay Kit **Invitrogen - Thermo** 

Fisher Catalog #P11496

AmpureXP beads Beckman

Coulter Catalog #A63880

ABSTRACT

V3 only primers can be found here:

https://github.com/sarahreiling/artic-ncov2019/blob/master/primer\_schemes/nCoV-2019/V3/nCoV-2019\_V3only.scheme.bed

# Primer pool preparation

#### 1 PRIMER POOL PREPARATION

If required resuspend lyophilised primers at a concentration of 100  $\mu M$  each

V3 only primers for this protocol were designed using <u>Primal Scheme</u> and generate overlapping 400 nt amplicons. Primer names and dilutions are listed in the table below.

 $\frac{https://github.com/sarahreiling/artic-ncov2019/blob/master/primer\_schemes/nCoV-2019/V3/nCoV-2019\_V3only.scheme.bed}{2019\_V3only.scheme.bed}$ 

2 Generate primer pool stocks by adding **5 μl** of each primer pair to a **1.5 mL** Eppendorf labelled either "Pool 1 (100μM)" or "Pool 2 (100μM)". Total volume should be **490 μl** for Pool 1 (100μM) and **490 μl** for Pool 2 (100μM). These are your 100μM stocks of each primer pool.

Primers should be diluted and pooled in the **mastermix cabinet** which should be cleaned with decontamination wipes and UV sterilised before and after use.

3 Dilute this primer pool 1:10 in molecular grade water, to generate 10 µM primer stocks. It is recommend that multiple aliquots of each primer pool are made to in case of degradation or contamination.

Primers need to be used at a final concentration of  $0.015\mu M$  per primer. In this case both pools have 98 primers in so the requirement is  $3.65\mu L$  primer pools (10uM) per 25 $\mu L$  reaction. For other schemes, adjust the volume added

protocols.io
2
12/21/2020

Citation: Sarah J Reiling, Anne-Marie Roy, Shu-Huang Chen, Ioannis Ragoussis (12/21/2020). nCoV-2019 McGill Artic PCR Protocol, 5 ul RT and V3 only at 63C. https://dx.doi.org/10.17504/protocols.io.bj9hkr36

## Multiplex PCR

4



#### **MULTIPLEX PCR**

In the extraction and sample addition cabinet add 35 µl RT product to each tube and mix well by pipetting.

The **extraction and sample addition cabinet** should should be cleaned with decontamination wipes and UV sterilised before and after use.

5 In the mastermix hood set up the multiplex PCR reactions as follows in 0.2mL 8-strip PCR tubes:

Component	Pool 1 [10 uM primer]		Pool 2 [10 uM]
Q5 Hot Start High-Fidelity 2X Master Mix	<b>□</b> 12.5 μl	<b>□12.5</b> μ	ıl
Primer Pool 1 or 2 (10µM pool 1+2)	<b>⊒</b> 3.7 μl	<b>⊒</b> 3.7 μl	
Nuclease-free water	⊒3.8 μl	<b>⊒</b> 3.8 μl	
Total	<b>⊒20 μl</b>	<b>⊒20</b> µl	

Add 20ul of PCR mastermix to the 5 ul RT product of step 10.

A PCR mastermix for each pool should be made up in the **mastermix cabinet** and aliquoted into PCR strip tubes. Tubes should be wiped down when entering and leaving the mastermix cabinet.

- 6 Pulse centrifuge the tubes to collect the contents at the bottom of the tube.
- 7 Set-up the following program on the thermal cycler:

Step	Tempe	rature Time		Cycles
Heat Activation	8 98 °C	© 00:00:30	1	
Denaturation	8 98 °C	<b>© 00:00:15</b>	36	
Annealing	₿ 63 °C	<b>© 00:05:00</b>	36	
Hold	8 4 °C	Indefinite	1	

Cycle number should be 25 for Ct 18-21 up to a maximum of 36 cycles for Ct 35

PCR clean-up

12/21/2020

#### **PCR CLEANUP**

Combine the entire contents of "Pool 1" and "Pool 2" PCR reactions for each biological sample into to a single **1.5 mL** Eppendorf tube.

9 Clean-up the amplicons using the following protocol:

Add an equal volume (1:1) of AmpureXP beads to the sample tube and mix by pipetting.

Incubate for 5 min at room temperature.

Pellet on magnet for 5 min. Remove supernatant.

Add 200 ul of 80% ethanol to the pellet and wash twice.

Let the beads dry for 3 min.

Add 30 ul elution buffer and resuspend the beads. Incubate for 3 minutes.

Pellet on magnet for 5 min. Remove and keep eluate (30 ul).

Amplicon clean-up should be performed in the **post-PCR cabinet which** should should be cleaned with decontamination wipes and UV sterilised before and after use.

Amplicon Quantification and normalisation

10



### **AMPLICON QUANTIFICATION AND NORMALIZATION**

Quantify the amplicon pools using a fluorimetric dsDNA assay. (e.g. PicoGreen with a standard curve 0-200ng)

We expect following concentrations:

## Pool 1+2 combined:

100-150 ng/ul for Ct 14-24 30-80 ng/ul for Ct 25-29 10-30 ng/ul for Ct 30-36

11



# Nextera Flex Library Prep:

After quantification of Pool 1+2, take a new plate and add 150 ng of Pool 1+2 and add up with nuclease-free water to a total volume of **30 ul (= 5 ng/ul)**.

#### Nanopore Library Prep:

After quantification of Pool 1+2, take a new plate and add 200 ng of Pool 1+2 and add up with nuclease-free water to a total volume of **20 ul (= 10 ng/ul)**.