



Upload image

Version 3

May 22, 2020

Amplicon clean-up using SPRI beads for RAPID nanopore kit RBK004 V.3

Forked from [Amplicon clean-up using SPRI beads](#)

Muhammad Faisal¹, Olin Silander¹, Nikki Freed¹

¹Massey University

In Development

[dx.doi.org/10.17504/protocols.io.bgsrjwd6](https://doi.org/10.17504/protocols.io.bgsrjwd6)

Coronavirus Method Development Community

Nikki Freed

MATERIALS

NAME	CATALOG #	VENDOR
Agencourt AMPure XP beads		

STEPS MATERIALS

NAME	CATALOG #	VENDOR
Agencourt AMPure XP	A63880	Beckman Coulter

MATERIALS TEXT







Freshly prepared 80% ethanol
10 mM Tris-HCl pH 8.0 with 50 mM NaCl

- 1 Vortex SPRI beads thoroughly to ensure they are well resuspended, the solution should be a homogenous brown colour.

Agencourt AMPure XP
by Beckman Coulter
Catalog #: [A63880](#)

Ampure XP bead clean up

- 2 Add an equal volume (1:1) of SPRI beads to the sample tube and mix gently by either flicking or pipetting. For example add **50 µl** room temperature SPRI beads to a **50 µl** reaction.
- 3 Pulse centrifuge to collect all liquid at the bottom of the tube.
- 4 Incubate for **00:05:00** at room temperature.
- 5 Place on magnetic rack and incubate for **00:02:00** or until the beads have pelleted and the supernatant is completely clear.

- 6 Carefully remove and discard the supernatant, being careful not to touch the bead pellet.
- 7 Add  **200 µl** of freshly prepared room-temperature  **80 % volume** ethanol to the pellet.
- 8 Keeping the magnetic rack on the benchtop, rotate the bead-containing tube by 180°. Wait for the beads to migrate towards the magnet and re-form a pellet. Remove the ethanol using a pipette and discard.
- 9  and repeat ethanol wash.
- 10 Pulse centrifuge to collect all liquid at the bottom of the tube and carefully remove as much residual ethanol as possible using a P10 pipette.
- 11 With the tube lid open incubate for  **00:01:00** or until the pellet loses its shine (if the pellet dries completely it will crack and become difficult to resuspend).
- 12 Remove the tube from the magnetic rack. Resuspend pellet in  **10 µl** 10 mM Tris-HCl pH 8.0 with 50 mM NaCl, mix gently by flicking and incubate at room temperature for  **00:02:00**.
- 13 Place on magnet and transfer sample to a clean 1.5mL Eppendorf tube ensuring no beads are transferred into this tube.