



Version 2

Feb 02, 2021

# An X-HTDC method for estimating particulate phosphorus from microalgae V.2

Yingyu YY Hu<sup>1</sup>, Zoe V. Finkel<sup>1</sup><sup>1</sup>Dalhousie University

1

Works for me

This protocol is published without a DOI.

Marine Microbial Macroecology Lab

Tech. support email: [ruby.hu@dal.ca](mailto:ruby.hu@dal.ca)Yingyu Hu  
Dalhousie University

## ABSTRACT

Here we describe a protocol to estimate particulate phosphorus associated with microalgae using an extra-high temperature dry combustion method (X-HTDC). An 800 °C-combustion can completely decompose intracellular macromolecules, improving estimates of particulate phosphorus relative to methods that use lower combustion temperatures. Based on tests on *Prorocentrum triestinum* in exponential growth in f/2 medium, estimates of intracellular particulate phosphorus are increased by 10% compared with the traditional 500 °C combustion method.

We recommend microalgae be collected on polycarbonate filters and then combusted at 800 °C before digestion by 0.2 M HCl for 30 minutes at 90 °C. The resulting orthophosphate is detected by mixing the digested sample with a mixture of molybdate and ascorbic acid to produce molybdenum blue (Chen 1956).

The working range of this assay is 0 to 200 uM orthophosphate; and the recommended maximum phosphorus per sample filter is about 30 ug.

In order to assess the intracellular phosphorus in microalgae, we recommend an oxalate reagent (Tovar-Sanchez 2003) to wash the microalgae collected on the filter to remove surface adsorbed phosphorus.

P.S. Chen, T.Y. Toribara and Huber Warner. Microdetermination of Phosphorus. Anal. Chem..

<https://doi.org/10.1021/ac60119a033>

AntonioTovar-Sanchez, Sergio A Sañudo-Wilhelmy, Manuel Garcia-Vargas, Richard S Weaver, Linda C Popels, David A Hutchins. A trace metal clean reagent to remove surface-bound iron from marine phytoplankton. Marine Chemistry.

[https://doi.org/10.1016/S0304-4203\(03\)00054-9](https://doi.org/10.1016/S0304-4203(03)00054-9)

## PROTOCOL CITATION

Yingyu YY Hu, Zoe V. Finkel 2021. An X-HTDC method for estimating particulate phosphorus from microalgae. **protocols.io**

<https://protocols.io/view/an-x-htdc-method-for-estimating-particulate-phosph-brhtm36n>

Version created by Yingyu Hu

## WHAT'S NEW

(1) How to prepare crucibles for the combustion (2) How to collect samples

#### KEYWORDS

particulate phosphorus, intracellular phosphorus, phosphomolybdenum-ascorbic reduction, orthophosphate, oxalate reagent, adsorbed phosphorus, X-HTDC, High temperature dry combustion

#### LICENSE

\_\_\_\_\_ This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

#### CREATED

Jan 15, 2021

#### LAST MODIFIED

Feb 02, 2021

#### PROTOCOL INTEGER ID

46355

#### SAFETY WARNINGS

Polycarbonate filter can release toxic gas and smoke during combustion. An exhaust system is required for muffle furnace while using the X-HTDC method.

#### BEFORE STARTING

Acid wash crucibles and covers as other glassware. Pre-combust crucibles and covers at 500 °C for 6 hours (ramp-rate: 150 °C/hour). Some crucibles or covers might shatter. The remainders maintain in good shape when are used to process samples at 800 °C.

### Sampling

## 1 Sampling microalgae for total particulate phosphorus (i.e. intracellular phosphorus and adsorbed phosphorus)

### 1.1 Filter microalgae in liquid media onto polycarbonate filters, using gentle vacuum pressure (5 inches Hg).

Filter forceps  
blunt end, stainless steel  
Millipore      XX6200006P

### 1.2 Rinse samples with filtered seawater

### 1.3 Place sample filters in 2 mL Cryogenic Vials.

Cryogenic Vials with Closures  
Polypropylene, 2 mL  
Corning®      66021-974

1.4 Filter blank media (without cells) through polycarbonate filter as blank.

1.5 Flash freeze filters and store at  $-20^{\circ}\text{C}$ .

## 2 Sampling microalgae for intracellular particulate phosphorus

2.1 Filter microalgae in liquid media onto polycarbonate filters, using gentle vacuum pressure (5 inches Hg).

Filter forceps  
blunt end, stainless steel  
Millipore      XX6200006P

2.2 Add  $5\text{ mL}$  oxalate reagent onto the filter, and let oxalate reagent sit in the filter funnel for 5m  
⌚ 00:05:00



Preparation of oxalate reagent  
by Yingyu Hu,  
Dalhousie University

PREVIEW

RUN

2.2.1 Add  $50\text{ mL}$  MilliQ water in a 250 mL beaker.

2.2.2 Weigh  $40\text{ g}$  NaOH and slowly pour into the beaker.

2.2.3 Use squeeze bottle to rinse the weighing boat and transfer rinse water into the same beaker.

2.2.4 Use glass rod to gently stir and fully dissolve NaOH.

The solution is very hot and corrosive. It can cause skin burns and eye damage.

Carefully transfer NaOH solution into 100 mL volumetric flask by using glass rod.

2.2.5

2.2.6 Rinse beaker with small amount of MilliQ water three times, transfer rinse water into the flask.

2.2.7 Mix the solution by gently shaking the capped volumetric flask and top to 100 mL with MilliQ water.

2.2.8 Transfer the prepared reagent into a 250 mL PP bottle.

2.2.9 Label the bottle with SDS pictogram.



2.10 In a 1000 mL beaker with stir bar, add **600 mL** MilliQ water.

2.11 Add **18.6 g** EDTA, **14.7 g** sodium citrate, **0.74 g** KCl and **5 g** NaCl into the beaker, stir until all ingredients are dissolved. **pH5.7**

2.12 **[M]10 Molarity (M)** NaOH is added dropwise to bring pH in between 6 to 7 by using a transfer pipet

2.13 Add **12.6 g** oxalic acid to the solution, stir the mixture while heating.

2.14 After oxalic acid is completely dissolved, stop heating and let it cool to room temperature. A water bath filled with tap water can be used to speed up cooling. **pH3.3**

2.15 Add **[M]10 Molarity (M)** NaOH dropwise to bring pH to **pH8**

2.16 Top to 1 L in volumetric flask with MilliQ water.

2.17 Filter oxalate reagent by rapid flow to a 1 L PP bottle.

Sterile Disposable Filter Units with PES  
Membrane  
Thermo Scientific™ Nalgene™ Rapid- 5964  
Flow™ 520

2.18 Label the bottle and keep it at 🧊 **Room temperature** .

2.3 Drain and then rinse the sample with filtered seawater once

2.4 Place sample filters in 2 mL Cryogenic Vials.

Cryogenic Vials with Closures  
Polypropylene, 2 mL  
Corning® 66021-974

2.5 Filter blank media (without cells) through polycarbonate filter as blank.

2.6 Flash freeze filters and store at 🧊 **-20 °C** .

#### X-HTDC-ing

- 3 Mark number at the bottom of each crucible with pencil, log the following information:
- (1) The number of crucible
  - (2) The code of sample in the crucible

Porcelain crucibles  
40 mL  
VWR 89037-996

Crucible cover

VWR 71000-146

4 Transfer sample to crucible with clean filter forceps.

5 [M]0.17 M  $\text{MgSO}_4$  reagent:

Dissolve 1.023 g  $\text{MgSO}_4$  in 50 mL MilliQ water

⊗ Magnesium sulfate anhydrous Fisher

Scientific Catalog #M65500

6 Add 200  $\mu\text{l}$  [M]0.17 M  $\text{MgSO}_4$  to each sample and blank filter.

7 Partly cover the crucibles and place in the oven at 90 °C until samples are completely dry.

Forced air oven

VWR 89511-410

8 Combust dried samples at 800 °C for 09:00:00

9h

Muffle furnace

F30428C

Thermo 10-505-13

Map the location of crucibles in the oven, in case pencil mark disappears under 800°C.

Ramp rate should be controlled at < 200 °C /hour or follow the instruction provided by manufacture, otherwise the crucibles might shatter.

SP.RAT: 150/PAMPU: hour

Or

SP.RAT: 2.5/PAMPU: minute

- 9 Allow samples to gradually cool down in the muffle furnace.
- 10 Pencil mark on crucibles should be still visible, however, it can be easily removed by water. Therefore, when removing samples out of the furnace, label the lid and crucible with sharpie immediately.

#### Digesting

- 11 **[M]0.2 M** HCl reagent:

In a reagent bottle, dissolve one part of **[M]12 N** HCl in 59 parts of MilliQ water

**[M]12 N Hydrochloric acid Contributed by users**

Volume of HCl\_0.2M\_mL = (5\_mL) X (#Sample + #Blank)

- 12 Preheat oven to **90 °C**

- 13 Add 2.5 mL 0.2 M HCl to each crucible.

- 14 Gently swirl the crucible.

- 15 Cover the crucibles and place crucibles in the oven for **00:30:00**

30m

- 16 Cool samples down to **Room temperature**

- 17 Gently swirl the crucible and then transfer 500 ul solution to 2 mL microtube. Duplicate each sample and blank.

Maxymum Recovery® Snaplock  
Microcentrifuge Tube  
2.0 mL, Polypropylene, Clear, Nonsterile,  
Axygen® MCT-200-L-C

#### Preparing standard working solutions

2h

- 18 Standard working solutions and reagents can be prepared during sample digestion.

## 19 KH<sub>2</sub>PO<sub>4</sub> primary standard stock solution ( $\approx$ 1 mM)

 Potassium dihydrogen orthophosphate ACP

Chemicals Catalog #P-4550

19.1 Transfer about 1 g KH<sub>2</sub>PO<sub>4</sub> into a beaker, cover the beaker with foil

19.2 Place the beaker into an oven, dry KH<sub>2</sub>PO<sub>4</sub> at  $\uparrow$  110 °C for at least  02:00:00

2h

19.3 Move KH<sub>2</sub>PO<sub>4</sub> into a vacuum desiccator, allow KH<sub>2</sub>PO<sub>4</sub> to cool to room temperature

19.4 Dissolve around 0.136 g dried KH<sub>2</sub>PO<sub>4</sub> in 1 L milliQ water.

- Use 1 L volumetric flask
- Take notes of the actual weight of KH<sub>2</sub>PO<sub>4</sub> for final concentration of standard stock solution

19.5 Transfer standard stock solution into a 1 L bottle and store in the fridge.

This stock solution lasts quite a long time, unless there is evidence for growth of algae or other extraneous biotic material.

## 20 Standard working solution

KH <sub>2</sub> PO <sub>4</sub>	Primary (ul)	MilliQ (ul)
S1	0	1000
S2	5	995
S3	10	990
S4	20	980
S5	50	950
S6	100	900
S7	150	850
S8	200	800

21 Transfer 500 ul of each standard working solution to 2 mL microtube.

Preparing working reagents 2h

## 22

All reagents are freshly prepared before colorimetric measurement.



## 23 [M]6 N (3 M) sulfuric acid reagent:

Carefully add 1 part [M]18 M concentrated sulfuric acid into 5 part MilliQ water

[18M sulfuric acid](#) **Contributed by users**

## 24 [M]2.5 % ammonium molybdate reagent:

Weigh [0.25 g](#) ammonium molybdate in a Falcon tube and top to [10 g](#) with MilliQ water.

Cap and shake until totally dissolved.

[Ammonium molybdate](#) **Sigma**

**Aldrich Catalog #09878-100G**

## 25 [M]10 % ascorbic acid reagent:

Weigh [1 g](#) ascorbic acid in a Falcon tube and top to [10 g](#) with MilliQ water;

Cap and shake until all dissolved.

[Ascorbic acid](#) **Sigma**

**Aldrich Catalog #A5960-100G**

Wrap the tube with foil if the reagent is not used right after prepared.

## 26 Calculate the volume of molybdate-ascorbic reagent:

Total volume of reagent\_mL = (0.5 mL) X (#standard working solution + #samples + #blanks)

## 27 Mix the reagents into Falcon tube:

Reagent	Part(s) as in volume
MilliQ	2
6N sulphuric acid	1
2.5% ammonium molybdate	1
10% ascorbic acid	1

Colorimetric measurement

2h


## 28 Preheat incubator/shaker to [37 °C](#)

SHAKING INCUBATOR

71L

Corning® LSE™

6753

29 Add  500 µl reagent to each standard, sample and blank

Finntip Stepper Tips

5 mL

Thermo Scientific

9404200

Before dispensing the reagent, wipe or dab the liquid drop on the outside of the tip, avoid wiping the open tip.

30 Vortex each tube.

31 Incubate at  37 °C for  03:00:00 while shaking at 200 rpm

3h

32 Load microplate with 250 µl reactant from each tube, duplicate.

	<u>1</u>	<u>2</u>	<u>3</u>	<u>4</u>	<u>5</u>	<u>6</u>	<u>7</u>	<u>8</u>	<u>9</u>	<u>10</u>	<u>11</u>	<u>12</u>
<u>A</u>	S1	S1	Samples and sample blanks: 40 with duplicate									
<u>B</u>	S2	S2										
<u>C</u>	S3	S3										
<u>D</u>	S4	S4										
<u>E</u>	S5	S5										
<u>F</u>	S6	S6										
<u>G</u>	S7	S7										
<u>H</u>	S8	S8										

96-Well Microplates, Polystyrene, Clear,

Greiner Bio-One

655101

33 Read plate in microplate reader

A	B
Shake duration	00:00:05
Shaking type	Continuous
Shaking force	High
Shaking speed [rpm]	600
Wavelength [nm]	820
Use transmittance	No
Pathlength correction	No
Measurement Time [ms]	100

Varioskan LUX Multimode Microplate  
Reader

Thermo Fisher      VL0L00D0

#### Calculating

- 34 Subtract the average absorbance at 820 nm of the blank standard replicates from the absorbance at 820 nm of all other standard working solutions.
- 35 Subtract the average absorbance at 820 nm of the blank sample (i.e. blank filter) replicates from the absorbance at 820 nm of all other individual samples.
- 36 Prepare a standard curve by plotting the average blank-corrected 820 nm absorbance for each standard working solution versus its concentration in  $\mu\text{M}$ .
- 37 Use the standard curve to determine the orthophosphate concentration of each unknown sample by using its blank-corrected 820 nm absorbance.
- 38  $(\text{P per sample})_{\mu\text{g}} = (\text{orthophosphate})_{\mu\text{M}} \times (V_{\text{HCl}})_{\text{mL}} \times (0.001) \times (30.97)$