


May 13, 2021

# Artisanal Production of Prefusion-Stabilized SARS-CoV-2 Spikes

Jean-Francois Rivest<sup>1</sup>, Claudia Goupil<sup>1</sup>, Yannick Doyon<sup>1</sup><sup>1</sup>Centre Hospitalier Universitaire de Québec Research Center and Faculty of Medicine, Laval University, Quebec City, QC G1V 4G2, Canada.

1 Works for me



Share

[dx.doi.org/10.17504/protocols.io.buyfnxtn](https://dx.doi.org/10.17504/protocols.io.buyfnxtn)

Coronavirus Method Development Community Doyon Lab

Jean-Francois Rivest

## ABSTRACT

The two-part protocol presented here describes the establishment of a stable pool of 293-F cells expressing the HexaPro variant of the severe acute respiratory syndrome coronavirus 2 (SARS-CoV2) spike (S) protein and its production and purification from the culture medium. CRISPR-Cas9-driven targeted integration of the expression cassette at the AAVS1 safe harbor locus streamlines cell line production. Engineered 293-F cells grown as suspension cultures in animal origin-free, chemically defined, protein-free medium allow purification to be performed directly from the culture medium with minimal cell processing steps. Yields of native and purified Spikes typically average 30-50 mg per liter of culture medium following tandem affinity purification using nickel-coupled and Strep-Tactin resins. Purified proteins bind recombinant ACE2 and mAb CR3022 in ELISA-based assays. Establishing an engineered cell line from cryopreserved 293-F cells takes approximately 4 weeks, and subsequent production and purification take between 9 and 14 days. As shown for the D614G mutation, this system can be readily adapted to study Spike variants.

## ATTACHMENTS

Artisanal Production of  
Prefusion-Stabilized SARS-  
CoV-2  
Spikes\_manuscript.pdf

## DOI

[dx.doi.org/10.17504/protocols.io.buyfnxtn](https://dx.doi.org/10.17504/protocols.io.buyfnxtn)

## PROTOCOL CITATION

Jean-Francois Rivest, Claudia Goupil, Yannick Doyon 2021. Artisanal Production of Prefusion-Stabilized SARS-CoV-2 Spikes. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.buyfnxtn>

## KEYWORDS

SARS-CoV2, spike protein, 293-F, stable cell line, genomic safe harbor, CRISPR-Cas9, AAVS1 targeting

## LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

May 12, 2021

LAST MODIFIED

May 13, 2021

PROTOCOL INTEGER ID

49895

MATERIALS TEXT

### Equipment

Note: All vessels used for cell culture should be sterilized if necessary. All glassware used for protein purification should be rinsed before use with ultrapure water.

| A   |
|---|
| 0,45 µm sterile filters   |
| 1.5 × 10 cm glass chromatography column (Bio-Rad #7375011)                          |
| 15 and 50 ml sterile polypropylene conical tubes                                    |
| 250 ml Erlenmeyer culture flasks with vented cap                                    |
| 2,8 L Fernbach baffled culture flask with vented cap                                |
| 2.5 × 20 cm glass chromatography column (Bio-Rad #7374252)                          |
| 37°C tissue culture incubator with shaking capacity up to 120 RPM and adhesive mats |
| 37°C water bath   |
| 2.5 × 20 cm glass chromatography column (Bio-Rad #7374252)                          |
| 37°C tissue culture incubator with shaking capacity up to 120 RPM and adhesive mats |
| 37°C water bath   |
| 50 mL Erlenmeyer flask with phenolic screw cap                                      |
| 500 ml Nalgene PPCO Centrifuge Bottles  |
| 12-wells cell culture plates  |
| Beckman Avanti series centrifuge with JLA-10.500 fixed angle rotor                  |
| Benchtop centrifuge   |
| Benchtop rocking platform   |
| Hemocytometer or automated cell counter   |
| Magnetic stirrer  |
| Stain-free enabled imager (optional)  |

### List of reagents

| A  | B                  | C              |
|--|--------------------|----------------|
| REAGENTS   | Supplier           | Catalog number |
|  |                    |                |
| <b>For cell culture and protein production</b>                         |                    |                |
|  |                    |                |
| 293-F Cells  | Gibco              | 11625019       |
| AAVS1_puro_CAG_HexaPro   | Addgene            | 164077         |
| eSpCas9(1.1)_No_FLAG_AAVS1_T2  | Addgene            | 79888          |
| FreeStyle 293 Expression Medium  | Gibco              | 12338018       |
| Lipofectamine 3000   | Invitrogen         | L3000015       |
| Opti-MEM I Reduced Serum Medium  | Gibco              | 31985062       |
| Puromycin Dihydrochloride  | Gibco              | A1113803       |
| Trypan Blue Solution, 0.4%   | Gibco              | 15250061       |
| <b>For protein purification</b>  |                    |                |
| 10x Buffer BXT   | IBA Lifesciences   | 2-1042-025     |
| 10x Buffer W   | IBA Lifesciences   | 2-1003-100     |
| 2,2,2-Trichloroethanol (optional)                                      | Sigma              | T54801         |
| 4– 20% Mini-PROTEAN TGX Stain-Free Gels (or hand-cast acrylamide gels) | Bio-Rad            | 4568095        |
| Imidazole  | Fisher BioReagents | BP30550        |
| Magnesium chloride, anhydrous  | Sigma              | M8266          |
| Ni Sepharose 6 Fast Flow   | Cytiva             | 17531801       |
| Potassium chloride   | Sigma              | P9541          |
| Potassium phosphate monobasic  | Sigma              | P5379          |
| QC Colloidal Coomassie Stain   | Bio-Rad            | 1610803        |
| Sodium chloride  | Sigma              | S7653          |
| Sodium hydroxide solution, 1N  | Fisher             | SS266-1        |
| Sodium phosphate dibasic   | Sigma              | S0876          |
| Sodium phosphate monobasic   | Sigma              | S0751          |
| Strep-Tactin HRP   | IBA Lifesciences   | 2-1502-001     |
| Strep-Tactin XT 4Flow high-capacity resin                              | IBA Lifesciences   | 2-5030-010     |

#### REAGENTS SETUP

##### PBS 10X (1 L):

Add the powders to 600 ml of ultrapure water in a 1 L glass bottle while stirring, adjust pH to 7,4 using NaOH and fill to 1 L with ultrapure water. Dilute to 1X with ultrapure water before use.

|                                  |       |
|----------------------------------|-------|
| NaCl                             | 80g   |
| KCl                              | 2g    |
| Na <sub>2</sub> HPO <sub>4</sub> | 14,4g |
| KH <sub>2</sub> PO <sub>4</sub>  | 2,4g  |

##### His wash buffer (1 L):

Add the powders to 600 ml of ultrapure water in a 1 L glass bottle while stirring, adjust pH to 8,0 using NaOH and fill to 1 L with ultrapure water.

|  |         |
|--|---------|
| NaH <sub>2</sub> PO <sub>4</sub> (50 mM) | 6 g     |
| NaCl (300 mM)                            | 17,54 g |
| Imidazole (20 mM)                        | 1,36 g  |

#### His elution buffer (1 L):

Add the powders to 600 ml of ultrapure water in a 1 L glass bottle while stirring, adjust pH to 8,0 using NaOH and fill to 1 L with ultrapure water.

|  |         |
|--|---------|
| NaH <sub>2</sub> PO <sub>4</sub> (50 mM) | 6 g     |
| NaCl (300 mM)                            | 17,54 g |
| Imidazole (250 mM)                       | 17 g    |

#### Buffer W and QXT 1X (1 L):

These buffers are provided as 10X stocks. Dilute to 1X with ultrapure water before use.

#### Puromycin diluted solution (1 ml):

The puromycin used here is provided at a stock concentration of 10mg/ml. We recommend diluting this stock with sterile PBS to a concentration of 0,5mg/ml. The stock solution should be kept at -20°C for long term storage. We recommend keeping the diluted solution at 4°C up to a year. To make 1ml of diluted puromycin:

|                 |       |
|-----------------|-------|
| Puromycin stock | 50µl  |
| Sterile PBS     | 950µl |

### Thawing and general handling of 293-F cells before transfection (Timing: 11 days)

1w 4d

- 1 Before thawing the cells, prewarm 30 ml of FreeStyle 293 Expression Medium in a 250 ml tissue culture shaker flask at 37°C, preferably in a humidified atmosphere with 8% CO<sub>2</sub>. Also, prewarm a 9 ml aliquot of culture medium at 37°C in a 50 ml sterile polypropylene conical tube.

Note: While we used commercial 293-F cells from Gibco, it is likely that this protocol can be adapted for other 293-derived cell lines adapted for suspension and serum-free growth. We recommend culturing 293-F cells in medium without antibiotic, as we have seen reduced growth fitness when using Pen-Strep. Cells should be tested for the absence of mycoplasma contamination before use.

- 2 Rapidly thaw the cells in a 37°C water bath, disinfect with 70% ethanol and resuspend in the 9 ml medium aliquot. Centrifuge the cells at 340g for 5 minutes and discard the supernatant. Resuspend the cell pellet in the 30 ml aliquot of prewarmed medium.
- 3 Assess post-thaw cell count and viability using trypan blue and a hemocytometer or an automated cell counter. Incubate the cells in a 37°C incubator in a humidified atmosphere with 8% CO<sub>2</sub> and orbital shaking at 120 RPM.
- 4 Monitor cell growth and viability everyday, until cells reach a confluency of about  $2-3 \times 10^6$  cells/ml. At this point, cells should be split at  $0,5 \times 10^6$  cells/ml. <sup>4d</sup>
- 5 We recommend passaging the cells at least twice after thawing to ensure proper fitness before transfection. We also recommend transfecting cells at most at 10 passages post-thaw to create a stable cell line, hence the importance of creating a stock of low-passage frozen cells after the initial thaw. <sup>1w</sup>

**Lipofectamine transfection of 293-F cells (Timing: 3 days)**

3d 0h 10m

- 6 Prepare the two plasmid vectors ordered from Addgene (eSpCas9(1.1)\_No\_FLAG\_AAVS1\_T2 and AAVS1\_puro\_CAG\_HexaPro) as high-quality DNA preps. We recommend transfecting a control vector such as a mammalian GFP expression plasmid (750 ng) to monitor transfection efficiency and cell death during the subsequent puromycin selection step.

Note: Lipofection conditions, such as the ratios of Lipofectamine 3000 and P3000 reagents, may benefit from further optimisation. However, we found that the conditions presented here were acceptable for generating stable cell lines using 293-F cells. The transfection conditions used here can be scaled up or down.

- 7 Count and plate  $2 \times 10^5$  cells in a volume of 2ml per well in a 12-wells cell culture plate for each transfection condition.
- 8 Using Opti-MEM medium prewarmed at 37°C, prepare the following mixes in two tubes:

| Tube A             |                        |
|--------------------|------------------------|
| Reagent            | Amount                 |
| Lipofectamine 3000 | 3 $\mu$ l              |
| Opti-MEM medium    | Complete to 50 $\mu$ l |

| A                             | B                      |
|-------------------------------|------------------------|
| Tube B                        |                        |
| Reagent                       | Amount                 |
| eSpCas9(1.1)_No_FLAG_AAVS1_T2 | 250 ng                 |
| AAVS1_puro_CAG_HexaPro        | 500 ng                 |
| P3000 reagent                 | 2 $\mu$ l              |
| Opti-MEM medium               | Complete to 50 $\mu$ l |

- 9 Transfer the content of tube **A** to tube **B** and mix thoroughly by pipetting. Incubate this mix 10 minutes at room temperature. 10m
- 10 Add the content of tube **B** to the corresponding well of the 12-wells plate.
- 11 Incubate the cells in a 37°C incubator in a humidified atmosphere with 8% CO<sub>2</sub> without shaking for 72h. 3d

**Puromycin selection of the transfected 293-F cells (Timing: 8-11 days)**

1w 4d

- 12 Prewarm 8 ml FreeStyle Expression Medium in a sterile 50 ml Erlenmeyer flask. Transfer the transfected cells to the flask. Cells may clump and attach to the bottom of the wells, in which case gentle up and down pipetting is generally sufficient to resuspend them. Incubate the cells overnight in a 37°C incubator in a humidified atmosphere with 8% CO<sub>2</sub> and orbital shaking at 125 RPM. We recommend lightly unscrewing the cap from the flasks during incubation, as narrow flask caps such as these are usually not vented. 1d

- 13 The next day, add 10 µl of the diluted puromycin solution to the medium so that the final concentration is 0,5 µg/ml in a 10 ml total volume.
- 14 Monitor cell death and growth every 2-3 days until a confluency of  $2-2,5 \times 10^6$  cells/ml, generally 7-10 days after adding the antibiotic. At this point, the puromycin selection is generally complete. If cells reach confluency too early, an additional 1:2 to 1:3 split may be performed. 1w 3d
- 15 Once the pool of puromycin-selected cells reaches confluency, split the cells in culture medium without puromycin. At this point, we highly recommend taking an aliquot of medium from a confluent culture to confirm the secretion of the spike protein by western blot with a Strep-Tactin HRP conjugate.

#### Production and harvest from a stable spike-expressing cell line (Timing: 7-12 days)

1w 5d

16

To start a production, split the cells in a 2,8L Fernbach baffled culture flask to obtain  $0,5 \times 10^6$  cells/ml in a total volume of 500ml fresh 293 FreeStyle Expression medium.

Note: Volumes of medium and reagents are provided for protein production from 500ml of culture medium and may be scaled up or down accordingly. As with other chromatography columns, ensure that the volume of liquid inside the column is always sufficient to keep the resin submerged.

- 17 Incubate the cells 7-12 days in a 37°C incubator in a humidified atmosphere with 8% CO<sub>2</sub> and orbital shaking at 55 RPM. **Figure 2A** illustrates spike accumulation in the culture medium at several time points. 1w 5d
- 18 Once the desired time point has been reached, harvest the cultures, and divide equally between two 500 ml PPCO Nalgene bottles. Centrifuge at 1850g for 20 minutes. Transfer the supernatant to two new Nalgene bottles and repeat centrifugation. 20m
- 19 Harvest the supernatant from both bottles after the second centrifugation step and pass through a 0,45 µm sterile filter. Use filtered supernatant for subsequent purification.

#### Histidine tag purification of the spike protein from filtered culture supernatant (Timing: 1 day)

1d

- 20 Transfer the filtered supernatant to a clean 1 L container and add 30 ml of Ni Sepharose 6 Fast Flow resin pre-rinsed with PBS 1X. Incubate overnight on a rocking platform at room temperature. **Optional:** save an aliquot of the initial filtered supernatant, referred to here as "Hisinput". 12h
- 21 Transfer the supernatant mixed with the resin to a clean 2,5 cm x 20 cm chromatography column and wait for the beads to settle down. Next, collect the flowthrough. **Optional:** save an aliquot of this fraction, referred to here as "His flowthrough".
- 22 Add 150 ml of 1X PBS to wash the column and collect the flowthrough. **Optional:** save an aliquot of this fraction, referred to here as "His wash 1".
- 23 Add 150 ml of His wash buffer and collect the flowthrough. **Optional:** save an aliquot of this fraction, referred to here as "His wash 2".

- 24 Add 150 ml of His elution buffer and collect the following fractions:  
  
E1: 20ml (E1 His)  
  
E2: 40ml (E2 His)  
  
E3: 30ml (E3 His)  
  
E4: 30ml (E4 His)  
  
E5: 30 ml (E5 His)
- 25 Wash the column with 240 ml PBS and store at 4°C.
- 26 Assess protein recovery and purification efficiency by loading the different fractions onto an 8% polyacrylamide gel and performing polyacrylamide gel electrophoresis (PAGE). Total proteins can be visualized by staining the gel with colloidal Coomassie Blue. A typical gel showing the fractions obtained at the different steps of purification is shown in **Figure 2**. Total proteins can also be rapidly visualized by employing Stain-Free gels and a stain-free enabled imager. A less expensive alternative to commercial Stain-Free gels can be achieved by adding 0,5% (v/v) 2,2,2-Trichloroethanol (TCE) to a handcast polyacrylamide gel and visualizing total proteins on a stain-free enabled imager or a UV gel imager<sup>26</sup>. In our experience, the second elution (E2 His) contains most of the purified protein and is selected for subsequent Twin-Strep tag purification.

#### Twin-Strep tag purification of the spike protein from a histidine purified sample (Timing: 1 day)

1d

- 27 Add 10ml of Strep-Tactin XT 4Flow high-capacity resin to the selected His eluate fraction(s) in a clean 50-100 ml container. Incubate overnight on a rocking platform at room temperature. 12h
- 28 Transfer the His eluate mixed with the resin on a 1,5 cm x 10 cm chromatography column and wait for the beads to settle down. Collect the flowthrough. **Optional:** save an aliquot of this fraction, referred to here as "Strep-Tactin flowthrough".
- 29 Wash the column with 5 x 10 ml of 1X buffer W and collect the flowthrough. **Optional:** save an aliquot of each of these five fractions, referred to here as "Strep-Tactin wash 1-5".
- 30 Elute with 30 ml of 1X buffer BXT and collect the following fractions:  
  
E1: 5 ml (E1 ST)  
  
E2: 20 ml (E2 ST)  
  
E3: 5 ml (E3 ST)
- 31 Pass 150 ml of 3M MgCl<sub>2</sub> to regenerate the column, and re-equilibrate using 80 ml of 1X buffer W. Store the column at 4°C.
- 32 Assess protein recovery and purification efficiency by loading the different fractions onto an acrylamide gel as described in **step 26**. Protein yield may be assessed by methods such as Bradford assay and absorbance at 280nm using a UV spectrophotometer. We found that both resins performed well when reused to purify up to 3-4 batches. The


resulting purified HexaPro spike protein can be used in downstream applications. We observed that the protein seemed to remain stable at 4°C on the short term but advise flash-freezing followed by -80°C keeping for long term storage.

## Anticipated results

- 33 By following this protocol, we consistently obtain yields of 30-50 mg of purified protein per liter of culture supernatant as determined by Bradford assay and A280 absorbance. Most of the purified protein is eluted in the second (E2 ST) and third (E3 ST) fractions of the final elution step, at concentrations of up to 2 mg/ml. **Figure 1C** shows the consistency in purity obtained from different production lots using the protocol described here. **Figure 3** shows that the purified protein binds both to recombinant ACE2 and mAb CR3022 and that this protocol can be adapted to produce a D614G HexaPro spike variant. The vectors used in this protocol are available from Addgene (#79888 and #164077).


## Figures

34




**Figure 1** | Production and purification of the HexaPro variant of the SARS-CoV2 spike protein from AAVS1 targeted 293-F cells. **A)** Representation of the HexaPro variant of the spike protein, adapted from *Wrapp et al.* SS, signal sequence; NTD, N-terminal domain; RBD, receptor-binding domain; SD1/SD2, subdomains 1 and 2; S2', protease cleavage site; FP, fusion peptide; HR1, heptad repeat 1; CH, central helix; CD, connector domain; HR2, heptad repeat 2. The six stabilizing proline substitutions are indicated. **B)** Schematic of the AAVS1 safe harbor locus before and after CRISPR-driven targeting of the HexaPro expression construct. SA, splice acceptor site; 2A, 2A self-cleaving peptide sequence; Puro, puromycin resistance gene; pA, polyadenylation sequence; CAG, CMV early enhancer/chicken  $\beta$ -actin promoter; T4 Foldon, T4 fibrin trimerization motif; 8xHis, octa-histidine affinity tag; Twin-Strep, Twin-Strep affinity tag. The Addgene numbers of the vectors used in this protocol are shown. **C)** Batch-to-batch consistency in purity between different lots of protein produced and purified using this protocol. Purified proteins (10  $\mu$ g per sample) were loaded on a 4–20% acrylamide gel for PAGE followed by colloidal Coomassie staining. **D)** Illustration of the general workflow of this protocol. The first part of the procedure is the establishment of a stable pool of 293-F cells expressing the HexaPro spike protein. The second part is the production and purification of this protein from the culture medium.





**Figure 2** | Spike secretion in the culture medium and acrylamide gel visualization of typical purification fractions. **A)** Aliquots of culture medium from a culture of spike-expressing 293-F cells were taken at the indicated time points. 20µl per sample of culture medium mixed with Laemmli buffer were loaded on a handcast acrylamide gel containing 0,5% (v/v) 2,2,2-Trichloroethanol (TCE) for PAGE. The gels were visualized with a stain free-enabled imager (ChemiDoc MP, Bio-Rad). **B)** The different fractions from 8xHis tag purification with Ni Sepharose 6 Fast Flow that are described in the protocol were loaded on a gel and imaged as in **A)**. **C)** Same as in **B)** for Twin-Strep tag purification with Strep-Tactin XT 4Flow high-capacity resin. FT: flowthrough, E2 His: second eluted fraction from the histidine tag purification.



**Figure 3** | Functional binding of the purified spike in ELISA assays. ELISA assays were conducted with four different lots of recombinant spike produced using this protocol against either **A**) recombinant ACE2-Fc or **B**) the mAb CR3022 antibody. **C**) Same as in **A**) but with protein from lot 4 (HexaPro D614) and protein harboring the D614G mutation (HexaPro G614). Each data point is the average of two technical duplicates. ELISA assays were performed using a modified version of the protocol developed by *Anamat et al.* Briefly, wells from a 96-wells plates were coated overnight at 4°C with 100µl per well of 2µg/ml of purified spike diluted in PBS. Washes were performed between each subsequent step with 200µl of PBS 1X with 0,1% (vol/vol) Tween-20 (PBST). Blocking was performed by incubating 1h at room temperature with 200µl per well of 2% (wt/vol) bovine serum albumin (BSA, Sigma) prepared in PBST. Next, 100µl per well of blocking solution (PBST with 2% (wt/vol) BSA) containing either recombinant ACE2-Fc (a gift from Manuel Carruso) or mAb CR3022 were added to each well before incubation 1h at room temperature. After blocking, 100µl of a 1:60,000 (vol/vol) dilution of goat anti-human IgG-horseradish peroxidase (HRP) conjugated secondary antibody (AP112, Emd Millipore) prepared in the blocking buffer was added per well, followed by incubation 1h at room temperature. Revelation was performed by adding 100µl of Pierce TMB (3,3',5,5' tetramethylbenzidine) substrate (Thermo Scientific) per well. The plates were incubated 20 min at room temperature in the dark before the reaction was stopped by adding 100µl of 2M sulfuric acid per well. Optical density at 450nm was measured using a Tecan Spark 10M microplate reader.

Amanat, F., D. Stadlbauer, S. Strohmeier, T. H. O. Nguyen, V. Chromikova, M. McMahon, K. Jiang, G. A. Arunkumar, D. Jurczyszak, J. Polanco, M. Bermudez-Gonzalez, G. Kleiner, T. Aydiello, L. Miorin, D. S. Fierer, L. A. Lugo, E. M. Kojic, J. Stoeve, S. T. H. Liu, C. Cunningham-Rundles, P. L. Felgner, T. Moran, A. García-Sastre, D. Caplivski, A. C. Cheng, K. Kedzierska, O. Vapalahti, J. M. Hepojoki, V. Simon and F. Krammer (2020). A serological assay to detect SARS-CoV-2 seroconversion in humans. *Nature Medicine*.  
<https://doi.org/10.1038/s41591-020-0913-5>

Yuan, M., N. C. Wu, X. Zhu, C.-C. D. Lee, R. T. Y. So, H. Lv, C. K. P. Mok and I. A. Wilson (2020). A highly conserved cryptic epitope in the receptor binding domains of SARS-CoV-2 and SARS-CoV. *Science*.  
<https://doi.org/10.1126/science.abb7269>