




May 17, 2020

DH5 α bacteria transformation protocol

Yoan Coudert¹¹Ecole Normale Supérieure de Lyon**1** Works for me dx.doi.org/10.17504/protocols.io.bghijt4e

CoudertLab

 Yoan Coudert
Ecole Normale Supérieure de Lyon

- 1 Thaw bacteria on ice
- 2 Dispense 50 μ l in a microtube
- 3 Add 1 to 5 μ l of DNA, and carefully flick the tube to mix cells and DNA
- 4 Place the mixture on ice for 10 minutes
- 5 Heat shock at 42°C for 30 seconds
- 6 Place on ice for 5 minutes
- 7 Pipette 450 μ l of room temperature SOC (or LB) into the mixture
- 8 Incubate at 37°C for 1 hour. Shake vigorously at 250 rpm.
- 9 Warm selection plates at 37°C
- 10 Spread desired volume (150-350 μ l) onto selection plate
- 11 Incubate overnight at 37°C

12 Store at 4°C after incubation of needed