




Jun 28, 2021

Immunofluorescence for free-floating brain sections

Jill R. Crittenden¹, Tomoko Yoshida¹, Margaret I. Davis², Ann M. Graybiel¹¹McGovern Institute for Brain Research, Massachusetts Institute of Technology, Cambridge, MA, 02138, USA;²Laboratory for Integrative Neuroscience, National Institute on Alcohol Abuse and Alcoholism, National Institutes of Health, USA

1 Works for me



Share

dx.doi.org/10.17504/protocols.io.kracy2e

Jill Crittenden

ABSTRACT

This protocol describes how to prepare fixed brain sections and process them for immunofluorescence.

DOI

dx.doi.org/10.17504/protocols.io.kracy2e

PROTOCOL CITATION

Jill R. Crittenden, Tomoko Yoshida, Margaret I. Davis, Ann M. Graybiel 2021. Immunofluorescence for free-floating brain sections. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.kracy2e>



KEYWORDS

fluorescence, immunohistology, free-floating

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Nov 09, 2017

LAST MODIFIED

Jun 28, 2021

PROTOCOL INTEGER ID

8706

SAFETY WARNINGS

Paraformaldehyde (PFA) is toxic, a strong irritant, and flammable as a solid. Prepare and use solutions only under adequate ventilation and with eye and skin protection.

Solution preparation

1 0.4 M NaKPO₄ Phosphate Buffer Stock (dilute this to make 0.1M PB for PFA, mounting solution, and rinse buffer)

112 g K₂HPO₄, (potassium phosphate dibasic; MilliporeSigma cat. # PX1570)

21.2 g Na₂H₂PO₄H₂O (sodium phosphate monobasic; VWR cat. #BDH9298)

2 L Distilled H₂O

0.1 M NaPO₄ Phosphate Buffer Stock (use this to make PBST)

23.1 g Na₂HPO₄ (sodium phosphate dibasic; VWR cat. # 0404)

5.04 g NaH₂PO₄H₂O (sodium phosphate monobasic; VWR cat. # BDH9298)

PBST working solution

200 ml	0.1M NaPO ₄ Phosphate Buffer Stock for PBST
1800 ml	Distilled H ₂ O
17.5 g	NaCl sodium chloride, Granular (Macron Fine Chemicals cat. #7532)
0.4 g	KCl potassium chloride (Macron Fine Chemicals cat. #6858)
4 ml	Triton X-100 (Sigma cat. # T9284)

4% PFA (freshly depolymerized from paraformaldehyde) for perfusion

Prepare in a fume hood, formaldehyde is a biohazard.

Add 4 % w/v paraformaldehyde (Macron Chemicals cat. #2621-59) to 0.1M PB.

Heat mixture to 80-85°C (do not let it go beyond this temperature).

Remove from heat and stir with a stir-bar until all solids have dissolved and the solution is clear.

Allow to cool to room temperature, then filter the solution through filter paper (VWR cat. #28331-081).

If using the same day, keep at room temperature, otherwise, store at 4 °C.

Glycerol cryoprotectant solution

25 ml glycerol (VWR cat. #0854)

75 ml 0.1M PB with 0.1% sodium azide (MP Biomedicals cat. #0210289190)

TSA blocking buffer

Heat 500ml PBST to 50°C

Add 2.5g Blocking Reagent (0.5% final; Perkin Elmer cat. #FP1012)

Stir until dissolved and cool to room temperature.

Use a funnel to filter through paper (VWR cat. #28331-081).

Mounting Solution

Heat 500ml 0.1M PB with 0.1% sodium azide to 55°C

Add 0.4g gelatin (0.08% final; Sigma cat. #G1890)

Cool to room temperature.

Use a funnel to filter through paper (VWR cat. #28331-081).



Paraformaldehyde (PFA) is toxic, a strong irritant, and flammable as a solid. Prepare and use solutions only under adequate ventilation and with eye and skin protection.

<https://www.sigmaaldrich.com/MSDS/MSDS/DisplayMSDSPage.do?country=US&language=en&productNumber=158127&brand=SIAL&PageToGoToURL=https%3A%2F%2Fwww.sigmaaldrich.com%2Fcatalog%2Fproduct%2Fsial%2F158127%3Flang%3Den>

Brain tissue preparation

- 2 Brains are dissected from deeply-anesthetized mice that were trans-cardially perfused with 10 ml of room-temperature 0.9% saline followed by 60 ml of 4% PFA. PFA must be used with appropriate protection and ventilation. Brains are then post-fixed in 4% PFA in glass vials for 2 hours with gentle rocking at 4 °C. Replace PFA solution with glycerol cryoprotectant solution until the brain sinks (overnight), at 4 °C. Store the brains in the same solution at 4 °C until sectioning. For sectioning, remove the brains from storage cryoprotectant solution, freeze in powdered dry ice on a chuck, and cut at 30 µm intervals on a sliding microtome. Put sequential sections into adjacent wells of a multi-well plastic tray containing 0.1% sodium azide in 0.1 M PB for storage at 4 °C until immunostaining.



Paraformaldehyde (PFA) is toxic, a strong irritant, and flammable as a solid. Prepare and use solutions only under adequate ventilation and with eye and skin protection.

<https://www.sigmaaldrich.com/MSDS/MSDS/DisplayMSDSPage.do?country=US&language=en&productNumber=158127&brand=SIAL&PageToGoToURL=https%3A%2F%2Fwww.sigmaaldrich.com%2Fcatalog%2Fproduct%2Fsial%2F158127%3Flang%3Den>

Immunofluorescence reactions

- 3 Rinse sections in PBST, 3 x 2 minutes
- 4 Incubate in blocking solution, 20 minutes
- 5 Incubate in primary antibodies diluted in blocking solution, 1 overnight at 4°C
- 6 Incubate in secondary antibodies diluted in blocking solution, 2 hours at room temp. in dark
[We use Alexa fluorophore conjugated secondary antibodies (Thermo Fisher Scientific), diluted 1:300]
- 7 Rinse in PBST, 3 x 2 min
- 8 Rinse in 0.1M PB, 3 x 2 minutes
- 9 Submerge sections in mounting solution in a large glass Petri dish and use a paintbrush to move them onto glass microslides (VWR #48311-703)
- 10 Let sections sit flat, covered loosely by foil, until they appear almost dry, 15 - 60 minutes, depending on the room humidity
- 11 Coverslip with ProLong™ Gold Antifade Mountant with DAPI (Thermo Fisher Scientific cat. #P36931)
- 12 Store slides in the dark at 4 °C