

3

Jan 04, 2022

Pilot Exercise: Generating the Illumina SampleSheet and sharing data via BaseSpace V.3

Ruth Timme¹, Maria Balkey², Tunc Kayikcioglu², Candace Bias², Cameron Boerner², James Pettengill²

¹US Food and Drug Administration; ²US Food And Drug Administration

1

protocol .

GenomeTrakr
Tech. support email: genomeTrakr@fda.hhs.gov

Maria Balkey

This protocol describes the required steps to generate Illumina SampleSheet and provides guidelines for sharing data via BaseSpace with FDA for the 2021 Pilot Exercise

V2:
- added more details regarding data transfer through BaseSpace, including steps for sharing a project with FDA-CFSAN.
- replaced samplesheet example.

V3:
- replaced samplesheet example.

Ruth Timme, Maria Balkey, Tunc Kayikcioglu, Candace Bias, Cameron Boerner, James Pettengill 2022.
Pilot Exercise: Generating the Illumina SampleSheet and sharing data via BaseSpace. [protocols.io](https://protocols.io/view/pilot-exercise-generating-the-illumina-samplesheet-b3gqqjvw)
<https://protocols.io/view/pilot-exercise-generating-the-illumina-samplesheet-b3gqqjvw>
Maria Balkey

protocol ,

Jan 04, 2022

Jan 04, 2022

56560

Generating Sequencing SampleSheet

- 1
- Create your SampleSheet using Excel or a text editor. Name your sample sheet according your internal protocols, and use a *.csv extension. The sample sheet is organized in sections titled Header, Reads, Settings and Data. Section headings are case-sensitive and shown in brackets []. [SampleSheet.csv](#) Details for each section are described below.

1.1 Header Section.

A	B
Parameter	Description
Investigator Name	First Name and Last Name for laboratory scientist
Experiment Name	Pilot Exercise SARS-Cov-2
Date	YYYY-MM-DD
Workflow	GenerateFASTQ
Application	FASTQ Only
Instrument Type	MiSeq
Library Prep Kit	Swift Normalase Amplicon NGS Panels (SNAP)
Index Kit	Swift Normalase Amplicon NGS Panels (SNAP)
Description	FDA Wastewater Protocol Pilot Exercise
Chemistry	Amplicon

Table 1: Header Section in Sequencing SampleSheet

Reads Section

1.2

**Enter the length of reads you are generating (could be different from 151)

A	B
Parameter	Description
Number of cycles in Read 1	151
Number of cycles in Read 2	151

Table 2: Reader Section in Sequencing SampleSheet

1.3 Settings Section

A	B
ReverseComplement	0

Table 3: Settings Section in Sequencing SampleSheet

1.4 Data Section: Include the machine ID in the Sample_Name column. e.g. sample_B01_M01234 where M01234 corresponds to the machine ID. This is the ID that will get populated to the fastq file names.

A	B	C	D	E	F	G	H	I	J
Sample_ID	Sample_Name	Sample_Plate	Sample_Well	I7_Index_ID	index	I5_Index_ID	index2	Sample_Project	Description
WPP_sample_B.01	WPP_sample_B.01_M01234	1	A1	SU001	TTGTTCTTG	I5_1	CCGAACAACA	LabAbbreviation-WPP	
WPP_sample_C.01	WPP_sample_C.01_M01234	1	B1	SU002	AGGTTGTGT	I5_2	AAGAATCGGC	LabAbbreviation-WPP	
WPP_sample_SA-1.01	WPP_sample_SA-1.01_M01234	1	C1	SU003	TATTGGTTGG	I5_3	CACTTCGCTT	LabAbbreviation-WPP	
WPP_sample_SA-2.01	WPP_sample_SA-2.01_M01234	1	D1	SU004	CACTTCGCTT	I5_4	TATTGGTTGG	LabAbbreviation-WPP	

Table 4: Data Section in Sequencing SampleSheet (dashes converted to underscores for Illumina requirements)

1.5 If a new library or sequencing protocol is utilized for processing the pilot exercise samples, make the changes in the corresponding parameters for the sequencing sample sheet.

2 Starting your MiSeq run and connecting to BaseSpace:

When you set up the run on the MiSeq, select the option to log in to BaseSpace.

Transfer sequencing RUN through BaseSpace

3 Click the **Runs** tab in the Illumina BaseSpace website .

3.1 Select the **run** that you would like to share with the FDA LFFM wastewater project team.

3.2 Go to the summary tab and click at the share button.

3.3 Enter the email address for the FDA team (**covidtrkr@fda.hhs.gov**), and then click Add Collaborator.

3.4 Click Save Settings. Your run will be automatically shared with the FDA LFFM wastewater project team.

Share this run

210720 M04822 0084 000000000-JRT3L

Invite a collaborator
covidtrkr@fda.hhs.gov|

Optional Message

Add a personal message (optional)

Remaining characters : 140/140

☒ Share the associated Project(s) as well (**Recommended**)
Unchecking this box will limit what the receiptient will have access to.

ADD COLLABORATOR

Collaborators

GenomeTrakr@fda.h... Pending Read Only

CANCEL SAVE SETTINGS

Figure 1: Sharing data via BaseSpace by Email Invitation

Transfer sequencing PROJECT through BaseSpace

- 4 Click the **Projects** tab in the Illumina BaseSpace website .
 - 4.1 Select the **project** that you would like to share with the FDA LFFM wastewater project team.
 - 4.2 Click at the share project button.
 - 4.3 Enter the email address for the FDA team (**covidtrkr@fda.hhs.gov**), and then click Add Collaborator.
 - 4.4 Click Save Settings. Your project will be automatically shared with the FDA LFFM wastewater project team.