



SEP 28, 2023

OPEN ACCESS



**DOI:**  
[dx.doi.org/10.17504/protocols.io.ewov1q212gr2/v1](https://dx.doi.org/10.17504/protocols.io.ewov1q212gr2/v1)

**Protocol Citation:** NUS iGEM 2023. PCR Clean-up.  
**protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.ewov1q212gr2/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working  
 We use this protocol and it's working

**Created:** Sep 28, 2023

**Last Modified:** Sep 28, 2023

## PCR Clean-up

NUS iGEM<sup>1</sup>

<sup>1</sup>National University of Singapore

NUS iGEM 2023



NUS iGEM

National University of Singapore

### ABSTRACT

2023 NUS-Singapore iGEM team followed this protocol to isolate the known DNA fragments from the PCR product without running the gel electrophoresis.

### PROTOCOL MATERIALS



PB  
 buffer Qiagen Catalog #19066

Step 1



Buffer PE Qiagen Catalog #19065





In [3 steps](#)







### SAFETY WARNINGS



Proper lab PPE must be worn at all times.

**Keywords:** PCR Clean-up,  
PCR, DNA, DNA Isolation

- 1 Add 5 times the sample's volume of  PB buffer Qiagen Catalog #19066 into the PCR tube with the PCR product.
- 2 Transfer the whole solution into a QIAquick Spin Column (purple tube) and ensure that the sample is dripped onto the white membrane in the column.
- 3 Centrifuge the tube at  13 rpm, 00:01:00 . 1m
- 4 Discard the flow-through and place the QIAquick column back into the same tube.
- 5 Add  700 µL of  Buffer PE Qiagen Catalog #19065 into the QIAquick column.
- 6 Centrifuge it at  13 rpm, 00:01:00 .
- 7 Discard the flow-through and place the QIAquick column back into the same tube.
- 8 Add  700 µL of  Buffer PE Qiagen Catalog #19065 again into the QIAquick column.

- 9 Centrifuge it at  13 rpm, 00:01:00 .
- 10 Discard the flow-through and place the QIAquick column back into the same tube.
- 11 Centrifuge the emptied QIAquick column at  13 rpm, 00:01:00 to remove residual  Buffer PE Qiagen Catalog #19065 .
- 12 Transfer the QIAquick column into the newly labelled Eppendorf tube.
- 13 Add  30 µL of DI water into the QIAquick column.
- 14 Centrifuge it at  13 rpm, 00:01:00 , ensuring that the direction of the Eppendorf tube's cap is the same as the direction of spinning to avoid breaking. 
- 15 Discard the QIAquick column, the solution left in the Eppendorf tube contains the DNA fragments.
- 16 Use the Nanodrop to measure and record the purity and concentration of the DNA fragments.

Equipment	
<b>NanoDrop™ One/OneC Microvolume UV-Vis Spectrophotometer</b>	NAME
UV-Vis Spectrophotometer	TYPE
Thermo Scientific	BRAND
ND-ONE-W	SKU