



Version 3

Dec 09, 2020

SDS-PAGE gel electrophoresis V.3

Steven J Burgess¹¹University of Illinois at Urbana-Champaign

In Development

This protocol is published without a DOI.

Ag SynBio Lab UIUC



Steven Burgess

University of Illinois at Urbana-Champaign

ABSTRACT

SDS-PAGE gel electrophoresis protocol for analyzing samples from plant leaf tissue via immunofluorescence. In this protocol no Coomassie blue is added to samples, the reason is that this interferes with the fluorescent signal during immunoblot. Instead, samples have already been prepared in Laemmli buffer (minus coomassie, protein extraction procedure), the leading edge of samples can be visualized due to the presence of chlorophyll.

Note

- When using 15 well, 0.75 mm comb, try to limit the volume loaded to 10 µL to minimize the risk of spillover of protein between wells.
- Ensure to wipe the tip on the rim of the sample tube to remove the sample stuck to the outside of the tip.

Literature:

http://www.bio-rad.com/webroot/web/pdf/lsr/literature/Bulletin_6040.pdf
<https://www.bio-rad.com/webroot/web/pdf/lsr/literature/10026447.pdf>

PROTOCOL CITATION

Steven J Burgess 2020. SDS-PAGE gel electrophoresis. **protocols.io**
<https://protocols.io/view/sds-page-gel-electrophoresis-bqiemube>
Version created by Steven Burgess

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Dec 08, 2020

LAST MODIFIED

Dec 09, 2020

PROTOCOL INTEGER ID

45350

MATERIALS TEXT

- 4–20% Mini-PROTEAN[®] TGX[™] Precast Protein Gels, 15-well, 15 µL (Bio-Rad Laboratories; [4561096](#))
- Opening lever (Bio-Rad Laboratories; 456-0000)
- Chameleon[™] Duo Pre-stained Protein Ladder (LI-COR Biosciences; [NC0738562](#))
- Mini-PROTEAN[®] Tetra Vertical Electrophoresis Cell, 4-gel (Bio-Rad Laboratories; [1658004](#))
- 10x Tris/Glycine/SDS (Bio-Rad Laboratories; [1610732](#))
- Fisherbrand[™] Gel-Loading Tips, 1-200 µL (Fisher Scientific; [02-707-181](#))
- Imperial[™] Protein Stain (Thermo Scientific; [24615](#))s
- Bio-Rad Gel-Doc Imager (optional)

ABSTRACT

SDS-PAGE gel electrophoresis protocol for analyzing samples from plant leaf tissue via immunofluorescence. In this protocol no Coomassie blue is added to samples, the reason is that this interferes with the fluorescent signal during immunoblot. Instead, samples have already been prepared in Laemmli buffer (minus coomassie, protein extraction procedure), the leading edge of samples can be visualized due to the presence of chlorophyll.


Note

- When using 15 well, 0.75 mm comb, try to limit the volume loaded to 10 µL to minimize the risk of spillover of protein between wells.
- Ensure to wipe the tip on the rim of the sample tube to remove the sample stuck to the outside of the tip.

Literature:

http://www.bio-rad.com/webroot/web/pdf/lsr/literature/Bulletin_6040.pdf
<https://www.bio-rad.com/webroot/web/pdf/lsr/literature/10026447.pdf>

Prepare gel tank and buffers

- 1 Create a working dilution of Tris-Glycine running buffer (~  1 L is required per gel tank) by diluting 1:10 with d H₂O.
- 2 Carefully remove the comb from the precast gel and the tape across the bottom.
- 3 Assemble the Mini-PROTEAN electrophoresis cell and fill the inner chamber with buffer and the outer chamber up to the recommended mark

the volume varies depending on whether running 2 or 4 gels, the level is marked on the tank

- 4 Wash the wells with running buffer by pipetting up and down





This is done to remove residual acrylamide that may have collected in wells

Prepare Samples


10m

- 5 In fresh centrifuge tubes, create a dilution of each sample using 1x PEB, such that each sample is set at a concentration of 3 µg /mL of total soluble protein.

Recommended final volume ~  100 µl (this will allow for 10 samples) but will depend on the application

- 6 Heat samples at  95 °C for  00:05:00 5m
- 7 Spin down samples at  10000 x g for  00:05:00 to pellet debris 5m

- 8 Load  3 µl of Chameleon™ Duo Pre-stained Protein Ladder to the first well

- 9 Load  **10 µl** of each sample (30 µg of total soluble protein) per lane

Running Gel 10m

- 10 Run precast gels at 200 V for ~  **00:30:00** .

30m

Mini-PROTEAN Tetra Cell
Gel Electrophoresis Tank

Bio-rad Laboratories 1658005EDU 

Or until the samples have reached the end of the gel. For self-made gels, run at 80-120 V

- 11 Carefully open precast gel case using an opening lever, by inserting where the black arrows indicate on the gel case.

- 12 Remove stacking gel with a blade

- 13 Proceed either directly to protein transfer



Protein Transfer

by Steven Burgess,

University of Illinois at Urbana-Champaign

PREVIEW

RUN

or if the gel is as a loading control, place in a container and cover with

 **Imperial**  **Protein Stain** **Thermo**

Fisher Catalog #24615

- 13.1 Place the membrane and bottom stack in the middle of the cassette base

- 13.2 Place pre-cast gel in the middle of the membrane and roll to remove air-bubbles

- 13.3 Place the top stack on top of the gel, gently roll

13.4 Close the cassette lid, taking care not to disturb the gel

13.5 Place the cassette in the Trans-Blot Turbo and follow the instructions on the machine (Fast protocol TGX gel).

Trans-Blot® Turbo™
Protein transfer apparatus

Bio-rad Laboratories 1704150EDU [↗](#)

13.6 Either dry membrane and store at [4 °C](#) for later use or proceed immediately to fluorescent western protocol

Visualizing Gel (optional) 10m

30m

14 Gently agitate the gel on a rocking platform for [00:30:00](#).

15 Pour off the Imperial Protein Stain (collect as hazardous waste)

[Imperial Protein Stain Thermo](#)

Fisher Catalog #24615

16 Rinse the gel with dH₂O.

17 Cover the stained gel with dH₂O and gently agitate in a rocking platform for [00:30:00](#)

30m

18 Repeat steps 15-17 until the background of the gel is clear and blue protein bands can be clearly visualized

19 Image the gel on the Bio-Rad Gel Doc XR system trans-white illumination

Gel Doc XR+ Gel Documentation System
Gel Documentation System

Bio-rad Laboratories 1708195 [↗](#)