

May 14, 2024

Enzymatic treatment of free floating sections

DOI

dx.doi.org/10.17504/protocols.io.81wgbzyyqgpk/v1

Bryan Killinger¹

¹Rush University

Killinger



Bryan Killinger

Rush University

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.81wgbzyyqgpk/v1

Protocol Citation: Bryan Killinger 2024. Enzymatic treatment of free floating sections. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.81wgbzyyqgpk/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 14, 2024

Last Modified: May 14, 2024

Protocol Integer ID: 99781

Funders Acknowledgement:

Michael J. Fox Foundation

Grant ID: ASAP-021030

NINDS

Grant ID: 1R01NS128467

Abstract

General protocol for enzymatic treatment of free floating brain sections.



Prepare 40 micron thick free floating brain sections for enzymatic treatments.

30m

1 Wash in DM

10m

2 Wash in DM

10m

3 Wash in DM

10m

Equilibrate tissue sections in reaction buffer.

10m

4 For GluC treatment, incubate in 1X GluC reaction buffer (New England Biolabs)

5m

5 For Trypsin treatment, incubate in TBS (50mM Tris-HCl, pH 7.6, 150mM NaCl)

5m

Enzymatic Treatment

1d 8h

6 For Gluc, combine 100ul 2X GluC reaction buffer and 100ul GluC (New England BioLabs, 100 ng/ μ l dissolved in high-purity water). Incubate tissues with mixture at 37C for 16h.

16h

7 For Trypsin, combine 4ul of trypsin stock (Sigma-Aldrich, 1 mg/ml dissolved in 1 mM HCl, pH 3) and 196ul of TBS. Incubate tissues with the mixture at 37C for 16h.

16h

Wash

30m

8 wash tissue in DM

10m

9 wash tissue in DM

10m

10 wash tissue in DM

10m