



Aug 29, 2020

# Reverse transcription of complementary DNA

PLOS One

Madeline M Glennon<sup>1</sup>, Austin Skinner<sup>1</sup>, Mara Krutsinger<sup>1</sup>, [Marino Resendiz<sup>1</sup>](#)

<sup>1</sup>University of Colorado Denver

**1** *Works for me* [dx.doi.org/10.17504/protocols.io.birskd6e](https://dx.doi.org/10.17504/protocols.io.birskd6e)

[Marino Resendiz](#)

## EXTERNAL LINK

<https://doi.org/10.1371/journal.pone.0235102>

## THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Glennon MM, Skinner A, Krutsinger M, Resendiz MJE (2020) Translesion synthesis by AMV, HIV, and MMLV reverse transcriptases using RNA templates containing inosine, guanosine, and their 8-oxo-7,8-dihydropurine derivatives. PLoS ONE 15(8): e0235102. doi: [10.1371/journal.pone.0235102](https://doi.org/10.1371/journal.pone.0235102)

## DOI

[dx.doi.org/10.17504/protocols.io.birskd6e](https://dx.doi.org/10.17504/protocols.io.birskd6e)

## PROTOCOL CITATION

Madeline M Glennon, Austin Skinner, Mara Krutsinger, Marino Resendiz 2020. Reverse transcription of complementary DNA. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.birskd6e>

## MANUSCRIPT CITATION please remember to cite the following publication along with this protocol

Glennon MM, Skinner A, Krutsinger M, Resendiz MJE (2020) Translesion synthesis by AMV, HIV, and MMLV reverse transcriptases using RNA templates containing inosine, guanosine, and their 8-oxo-7,8-dihydropurine derivatives. PLoS ONE 15(8): e0235102. doi: [10.1371/journal.pone.0235102](https://doi.org/10.1371/journal.pone.0235102)

## EXTERNAL LINK

<https://doi.org/10.1371/journal.pone.0235102>

## LICENSE

———— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Jul 20, 2020

## LAST MODIFIED

Aug 29, 2020












## PROTOCOL INTEGER ID


39442



## STEPS MATERIALS

NAME	CATALOG #	VENDOR
AMV Reverse Transcriptase - 1,000 units	M0277L	New England Biolabs

## Reverse Transcription using AMV RT

- 1 Obtain an RNase-free microcentrifuge tube  **400 µl**  
Add the following components in the described order:
  - Water  **15 µl**
  - 10 X Buffer  
 **2 µl** **50 mM Tris-acetate, 75 mM Potassium Acetate, 8 mM Magnesium Acetate, 10 mM DTT (pH 8.3)**
  - 5'-P32 Radiolabeled DNA  **1 µl** < 1 pmol
  - RNA  **2 µl** app. 2 pmol
- 2 Place the tube in a heat block at  **90 °C** and turn off the heat block to allow for slow cooling. ⌚ **02:00:00**
- 3 Remove the tube upon reaching room temperature and centrifuge for ⌚ **00:00:05** using a microcentrifuge
- 4 The solution is then transferred to a separate microcentrifuge tube  **3 µl**
- 5 Add the desired dNTP solution  **1 µl** [**10 Millimolar (mM)**]
- 6 Add the reverse transcriptase solution  **2 µl** . The reverse transcriptase solution consists of RT  **2 µl** dissolved in water  **58 µl** . This enzyme was purchased from NEB.



**AMV Reverse Transcriptase - 1,000 units**  
by New England Biolabs  
Catalog #: **M0277L**
- 7 Place the solution in an incubator at  **37 °C** for ⌚ **00:45:00**
- 8 Remove from water bath and add urea  **7 µl** [**6 Molarity (M)**]

- 9 Centrifuge for 🕒 00:00:05 using a microcentrifuge
- 10 Place the tube in a heat block at 🔥 90 °C for 🕒 00:10:00
- 11 Centrifuge for 🕒 00:00:05 using a microcentrifuge while allowing to cool to rt
- 12 Withdraw the solution 📏 8 µl **This typically accounts for app 4000 cpm** and load onto a denaturing 20 % PAGE
- 13 On a separate well, add a marker composed of xylene cyanol and bromophenol blue in 90% formamide buffer
- 14 Turn on the power supply and allow to run until the bromophenol blue marker is 3/4 down the gel. Settings will vary depending on the gel size and thicknes.
- 15 Remove the gel from the gel stand and remove one glass slide.
- 16 Develop using an autoradiography cassette 🕒 Overnight
- 17 Visualize and quantify the gel using a Molecular Dynamics Phosphorimager 840

#### Reverse Transcription using MMLV RT

- 18 Repeat all steps above, except with the amount of reverse transcriptase added. Add the reverse transcriptase solution 📏 2 µl . The reverse transcriptase solution consists of RT 📏 1 µl dissolved in water 📏 10 µl . This enzyme was purchased from NEB M0253L at 200,000 U/mL

#### Reverse Transcription using HIV RT

- 19 Repeat all steps above, except with the amount of reverse transcriptase added. Add the reverse transcriptase solution 📏 2 µl . The reverse transcriptase solution consists of RT 📏 0.75 µl dissolved in water 📏 1.5 mL ; followed by a second 1:1 dilution in water. This enzyme was purchased from Worthington LS05003 at 27.5 units/uL