


Sep 23, 2022

Embryo stage *C. elegans* dissociation for FACS isolation and RNA-seq analysis of intestine-specific cells

 In 1 collectionRobert TP Williams¹, Erin Osborne Nishimura¹¹Colorado State University, Fort Collins

1 Works for me

 Sharedx.doi.org/10.17504/protocols.io.dm6gpbw9plzp/v1

Osborne Nishimura Lab

Robert TP Williams
Colorado State University, Fort Collins

ABSTRACT

This protocol is for generating a single cell suspension suitable for isolation of intestine-specific cells through Fluorescence Activated Cell Sorting (FACS) from embryo stage *C. elegans*. This protocol utilizes treatment with Chitinase and Pronase E to disrupt the cuticle. Embryos are mechanically homogenized with 21G syringe needle.

DOI

dx.doi.org/10.17504/protocols.io.dm6gpbw9plzp/v1

PROTOCOL CITATION

Robert TP Williams, Erin Osborne Nishimura 2022. Embryo stage *C. elegans* dissociation for FACS isolation and RNA-seq analysis of intestine-specific cells.
protocols.io
<https://protocols.io/view/embryo-stage-c-elegans-dissociation-for-facs-isolation-b35vqq66>



FUNDERS ACKNOWLEDGEMENT

National Institutes of Health

Grant ID: R35GM124877

Bridge to Doctorate at Colorado State University

Grant ID: 1612513

National Science Foundation

Grant ID: 2143849

Webb-Waring Biomedical Research Award

Grant ID: E.O.N.

COLLECTIONS ⓘ



Protocol collection: Dissociation and FACS isolation of embryonic and post-embryonic *C. elegans* intestine cells for RNA-seq analysis

KEYWORDS

C. elegans, FACS, single cell suspension, cell dissociation, embryo stage, intestine

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jan 21, 2022

LAST MODIFIED

Sep 23, 2022

PROTOCOL INTEGER ID

57237

PARENT PROTOCOLS

Part of collection

[Protocol collection: Dissociation and FACS isolation of embryonic and post-embryonic *C. elegans* intestine cells for RNA-seq analysis](#)

MATERIALS TEXT

Strains:

- FACS control *C. elegans* strain, i.e. N2
- FACS sorting *C. elegans* strain, i.e. JM149 *cal-51*[*elt-2p::GFP::HIS-2B::unc-54* 3'UTR + *rol-6*(su1006)]

Reagents:

L15-10 solution

- Leibovitz's L-15 Medium (Thermo 21083027)
- Fetal Bovine Serum (heat inactivated) (Thermo 10438026)
- 100X Penicillin Streptomycin solution (Thermo 15140148)
- Sucrose powder

Stock solutions for egg buffer

- 2M NaCl
- 2M KCl
- 1M CaCl₂
- 1M MgCl₂
- 1M HEPES pH 7.2

Enzymes

- Chitinase from *Streptomyces griseus* (Sigma C6137-5UN)
- Pronase E, Protease from *Streptomyces griseus* (Sigma P8811-1G)

Consumables:

- standard 1.5 ml tubes
- Stericup 0.2 micron filter (Fisher S2GPU05RE)
- 21 gauge 1 inch needle (fisher 14-826C)
- 1 ml syringe (fisher 14-823-30)
- 35-micron nylon mesh filter caps (Stellar Scientific FSC-FLTCP)
- 5 ml sterile polypropylene round-bottom tube (STEMCELL Technologies 38057)
- Bio-Rad TC20 automated cell counting slide (Bio-rad 1450011)

Equipment:

- Fixed angle rotor centrifuge (Eppendorf 5424)
- Swinging bucket rotor refrigerated centrifuge (eppendorf 5810R)
- 15 ml tube and 1.5 ml tube adapter (eppendorf 022638704, eppendorf 022638742)
- Fluorescent microscope
- Nutating mixer
- Bio-Rad TC20 automated cell counter

Before beginning

- 1 Prepare reagents in advance

L15-10 Buffer: Mix 500 ml Leibovitz's L-15 Medium, 50 ml Fetal Bovine Serum (heat inactivated), 50 ul of 100x Penicillin-Streptomycin solution and 7.7 g sucrose. Filter with 0.2 micron pore filter. Store at 4°C.

Egg Buffer: Mix 29.5 ml of 2M NaCl, 12 ml of 2M KCl, 1 ml of 1M CaCl₂, 1 ml of MgCl₂, 12.5 ml

of 1M HEPES-NaOH pH 7.2 and 435 ml molecular grade water. Filter with 0.2 micron pore filter. Store at 4°C.

Chitinase solution (1 U/ml): Dissolve 5 units of Chitinase from *Streptomyces griseus* (Sigma C6137-5UN) in 5 ml of Egg Buffer. Nutate the solution for approximately 10 minutes until dissolved. Prepare 1 ml aliquots in 1.5 ml tubes. Store aliquots at -20°C.

Pronase E solution (15 mg/ml): Weigh 150 mg of Protease from *Streptomyces griseus* (Sigma P8811-1G) into a 15 ml tube. Dissolve the enzyme in 10 ml of Egg Buffer. Nutate the solution for approximately 10 minutes until dissolved. Prepare 1 ml aliquots in 1.5 ml tubes. Store aliquots at -20°C.

2 On day of protocol:

Cool swinging bucket centrifuge to 4°C

Thaw pronase and chitinase aliquots at room temperature

Place L15-10 and egg buffer on ice

3 Starting material:

Worm suspension in 15 ml tube (material generated from [this protocol](#))

Strains: N2, fluorescent sorting strain

Perform this protocol on both strains in parallel

Note: The volumes for enzymatic treatments in this protocol require an embryo pellet less than 200 uL. If embryo pellet exceeds 200 ul, utilize 2x the embryo pellet volume.

Chitinase Treatment

4 Centrifuge embryo suspension at 2000 rcf for 1 minute in swinging bucket centrifuge

5 Resuspend the embryo pellet in 1 ml of M9 and transfer to a 1.5 ml tube.

6 Pellet the embryos at 2000 rcf for 1 minute in a centrifuge

7 Aspirate and discard the supernatant

8





Transfer 10 ul of embryo pellet to 1ml of Qiazol and store at -80°C for downstream RNA analysis

- 9 Resuspend the embryo pellet from Step 7 in 0.5 ml egg buffer and 1 ml chitinase (1 U/ml)
- 10 Incubate for 20 min rotating/nutating at room temperature



Verify eggshell digestion by visualizing the chitinase treated embryos under a microscope. Early embryos should change shape, and pretzel stage embryos should release from their eggshell.

- 12 Pellet the embryos at 200 rcf for 5 min in fixed angle rotor centrifuge
- 13 Aspirate and discard the supernatant

Pronase treatment and dissociation

- 14 Resuspend the chitinase treated embryo pellet in 200 ul pronase (15 mg/ml) and 500 ul egg buffer
- 15 Attach a 21 guage 1¼ inch needle to a sterile 1 ml syringe
- 16 Disrupt the embryo vitelline membrane and release cells by passing the embryo suspension through the needle 100 times, generating a worm slurry
- 17 Visually confirm embryo dissociation by viewing a 2 ul sample of worm slurry on a fluorescent microscope
- 18 Quench the pronase treatment by adding 800 ul of L15-10 media to worm slurry

19 Store the sample on ice until all strains are completed

Wash and harvest single cells

20 Wash away excess pronase from the worm slurry

20.1 Pellet the worm slurry at 500 rcf for 5 mins in swinging bucket centrifuge cooled to 4°C

20.2 Aspirate and discard the supernatant

20.3 Resuspend the worm slurry in 1 ml of L15-10 media.

20.4 Pellet the worm slurry at 500 rcf for 5 mins in swinging bucket centrifuge cooled to 4°C

20.5 Aspirate and discard the supernatant

20.6 Resuspend the worm slurry in 1 ml of L15-10 media.

21 Harvest the cells

21.1 

Pellet undissociated embryos at 100 rcf for 1 minute in swinging bucket centrifuge cooled to 4°C

NOTES:

- This step will separate the dissociated cells from intact embryos
- Cells will remain in the supernatant
- Ensure your cell type of interest is not lost during this step.
- Visually confirm fluorescent cells **are present** in the supernatant.
- Visually confirm fluorescent cells **are not present** in the pellet.
- You may need to reduce the centrifuge speed and/or time if fluorescent cells are in the pellet of this step.

21.2 Aspirate 1 ml of the cell-containing supernatant. Keep the pipette away from the pelleted worm debris.

21.3 Dispense the cell suspension through a 35-micron nylon mesh filter into a 5 ml flow cytometry tube

21.4 Pellet undissociated embryos at 100 rcf for 1 minute in swinging bucket centrifuge cooled to 4°C

22 Perform an additional round of cell harvest for the sorting strain only (Step 21)

Total cell suspension volumes:

- Control strain = 1ml
- Sorting strain = 2ml

23 - Transfer 70 ul of cells to 1 ml of Qiazol and store at -80°C for downstream RNA analysis

- Continue to step 24
- Retain the remaining ~2ml of cells for [FACS](#)

Measure approximate cell concentration

24 Load 10 ul of cell suspension to a Bio-Rad TC20 automated cell counting slide

25 This protocol should yield between 2×10^6 to 4×10^6 total cells

26 Dilute the sample to $1 \times 10^6 \frac{\text{cells}}{\text{ml}}$ if above this concentration with L15-10

27 Microscopically confirm fluorescent cells are present in the cell suspension

28 Move on to [FACS protocol](#)