

Aug 29, 2024

## Purification of BCL2L13-GFP

DOI

**[dx.doi.org/10.17504/protocols.io.36wgqno2ogk5/v1](https://dx.doi.org/10.17504/protocols.io.36wgqno2ogk5/v1)**

Elias Adriaenssens<sup>1</sup>

<sup>1</sup>Sascha Martens lab, University of Vienna, Max Perutz Labs - Vienna



Elias Adriaenssens

Sascha Martens lab, University of Vienna, Max Perutz Labs - ...

OPEN  ACCESS



DOI: **[dx.doi.org/10.17504/protocols.io.36wgqno2ogk5/v1](https://dx.doi.org/10.17504/protocols.io.36wgqno2ogk5/v1)**

**Protocol Citation:** Elias Adriaenssens 2024. Purification of BCL2L13-GFP. **protocols.io**

**<https://dx.doi.org/10.17504/protocols.io.36wgqno2ogk5/v1>**

**License:** This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** May 24, 2024

**Last Modified:** August 29, 2024

**Protocol Integer ID:** 103090

**Keywords:** ASAPCRN

**Funders Acknowledgement:**

**Aligning Science Across  
Parkinson's (ASAP)**

**Grant ID:** ASAP-000350

**Marie Skłodowska-Curie  
MSCA Postdoctoral  
fellowship**

**Grant ID:** 101062916

## Abstract

This protocol details the purification of BCL2L13-GFP.

## Materials

### Lysis buffer:

A	B
Tris-HCl	50 mM
pH	7.4
NaCl	300 mM
Triton X-100	1%
glycerol	5%
MgCl <sub>2</sub>	2 mM
DTT	1 mM
β-mercaptoethanol	2mM
cOmplete EDTA-free protease inhibitors (Roche)	
CIP protease inhibitor (Sigma)	
DNase (Sigma)	

### Wash buffer:

Tris-HCl	50 mM
pH	7.4
NaCl	300 mM
DTT	1 mM

### High salt wash buffer:

Tris-HCl	50 mM
pH	7.4
NaCl	700 mM
DTT	1 mM

### SEC buffer:

Tris-HCl	25 mM
pH	7.4
NaCl	300 mM
DTT	1 mM




- pET-DUET1 vector (available on Addgene)  pETDuet-1 TIM9,10 **addgene Catalog #170280**

- BCL2L13 W275A/I278A ( $\Delta$ LIR1)(available on Addgene)
- BCL2L13 Y213A/I216A/W275A/I278A ( $\Delta$ LIR1+2) (available on Addgene)
- BCL2L13 I224A/L227A/W275A/I278A ( $\Delta$ LIR1+3) (available on Addgene)
- BCL2L13 W275A/I278A/I307A/V310A ( $\Delta$ LIR1+4) (available on Addgene)
- BCL2L13 I224A/L227A/W275A/I278A/I307A/V310A ( $\Delta$ LIR1+3+4) (available on Addgene)
- Rosetta pLysS cells (Novagen Cat# 70956-4)

 Rosetta™(DE3)pLysS Competent Cells - Novagen **Merck Catalog #70956-4**

- SORVAL RC6+ centrifuge with an F21S-8x50Y rotor (Thermo Scientific)
- Glutathione Sepharose 4B beads (GE Healthcare)
- 10 kDa cut-off Amicon filter (Merck Millipore)

 Amicon® Ultra Centrifugal Filter, 10 kDa MWCO **Merck MilliporeSigma (Sigma-Aldrich) Catalog #UFC801008**

- GFP-tagged BCL2L13-GFP (available from Addgene),
- BCL2L13(W275A/I278A)-GFP ( $\Delta$ LIR1) (available from Addgene),
- BCL2L13(Y213A/I216A/W275A/I278A)-GFP ( $\Delta$ LIR1+2) (available from Addgene),
- BCL2L13(I224A/L227A/W275A/I278A)-GFP ( $\Delta$ LIR1+3)(available from Addgene),
- BCL2L13(W275A/I278A/I307A/V310A)-GFP ( $\Delta$ LIR1+4) (available from Addgene),
- BCL2L13(I224A/L227A/W275A/I278A/I307A/V310A)-GFP ( $\Delta$ LIR1+3+4)(available from Addgene)

## Purification - BCL2L13-GFP

20h 46m

### 1 To purify GFP-tagged

- BCL2L13-GFP (available from Addgene),
- BCL2L13(W276A/I279A)-GFP ( $\Delta$ LIR1) (available from Addgene),
- BCL2L13(Y213A/I216A/W276A/I279A)-GFP ( $\Delta$ LIR1+2) (available from Addgene),
- BCL2L13(I224A/L227A/W276A/I279A)-GFP ( $\Delta$ LIR1+3)(available from Addgene),
- BCL2L13(W276A/I279A/I307A/V310A)-GFP ( $\Delta$ LIR1+4) (available from Addgene),
- BCL2L13(I224A/L227A/W276A/I279A/I307A/V310A)-GFP ( $\Delta$ LIR1+3+4)(available from Addgene),

fuse the cytosol-exposed domain of BCL2L13 (1-463aa) to a C-terminal GFP-tag through cloning into a pET-DUET1 vector (available on Addgene).

### 2 Introduce the point mutants by in vitro mutagenesis to generate

- BCL2L13 W276A/I279A ( $\Delta$ LIR1)(available on Addgene),
- BCL2L13 Y213A/I216A/W276A/I279A ( $\Delta$ LIR1+2) (available on Addgene),
- BCL2L13 I224A/L227A/W276A/I279A ( $\Delta$ LIR1+3) (available on Addgene),
- BCL2L13 W276A/I279A/I307A/V310A ( $\Delta$ LIR1+4) (available on Addgene),
- BCL2L13 I224A/L227A/W276A/I279A/I307A/V310A ( $\Delta$ LIR1+3+4) (available on Addgene).

### 3 After the transformation of the pET-DUET1 vector encoding BCL2L13-GFP wild-type or mutants in E. coli Rosetta pLysS cells (Novagen Cat# 70956-4), grow the cells in 2x Tryptone Yeast extract (TY) medium at 37 °C until an OD<sub>600</sub> of 0.4 and then continued at 18 °C .



### 4 Once the cells reaches an OD<sub>600</sub> of 0.8, induce the protein expression with 100 micromolar ( $\mu$ M) isopropyl $\beta$ -D-1-thiogalactopyranoside (IPTG) for 16:00:00 at 18 °C . Collect the cells by centrifugation and resuspend in lysis buffer.

16h




Lysis buffer:

A	B
Tris-HCl pH 7.4	50 mM
NaCl	300 mM
Triton X-100	1%
Glycerol	5%
MgCl <sub>2</sub>	2 mM
DTT	1 mM




A	B
β-mercaptoethanol	2mM
cOmplete EDTA-free protease inhibitors (Roche)	
CIP protease inhibitor (Sigma)	
DNase (Sigma)	


- 5 Sonicate the cell lysates twice for 30 s and clear by centrifugation at  18000 rpm, 4°C, 00:45:00 in a SORVAL RC6+ centrifuge with an F21S-8x50Y rotor (Thermo Scientific).

45m





- 5.1 Sonicate the cell lysates for  00:00:30 (1/2).

30s

- 5.2 Sonicate the cell lysates for  00:00:30 (2/2).

30s

- 6 Collect the supernatant and incubate with pre-equilibrated Glutathione Sepharose 4B beads (GE Healthcare) for  02:00:00 at  4 °C with gentle shaking to bind BCL2L13-GFP.

2h



- 7 Centrifuge the samples to pellet the beads and remove the unbound lysate.



- 8 Wash the beads twice with wash buffer, once with high salt wash buffer, and two more times with wash buffer.




Wash buffer:

A	B
Tris-HCl pH 7.4	50 mM
NaCl	300 mM
DTT	1 mM

High salt wash buffer:

A	B
Tris-HCl pH 7.4	50 mM
NaCl	700 mM
DTT	1 mM




9 To cleave off the GST-tag  Overnight , elute the GFP-tagged cargo receptor from the GSH beads by the addition of TEV protease.



10 To collect the supernatant, collect the beads by centrifugation.



11 Wash the beads twice with  4 mL of wash buffer, and collect the supernatant.



12 Pool the supernatant fractions, filter through a 0.45 µm syringe filter, concentrate with 50 kDa cut-off Amicon filter (Merck Millipore), and load onto a pre-equilibrated Superdex 200 Increase 10/300 GL column (Cytiva).

13 Elute the proteins with SEC buffer.

SEC buffer:


A	B
Tris-HCl pH 7.4	25 mM
NaCl	300 mM
DTT	1 mM

14 Analyze fractions by SDS-PAGE and Coomassie staining. Pool fractions containing purified BCL2L13-GFP.

15 After concentrating the purified protein, aliquote the protein and snap-frozen in liquid nitrogen.



#### Note

Store the proteins at  -80 °C .