

MAR 10, 2023

Stern REHAB

Abby Myers¹

¹HHMI



Abby Myers

ABSTRACT

Instructions for taking care of Stern's REHAB vials.



Protocol Citation: Abby Myers 2023. Stern REHAB. **protocols.io**

https://protocols.io/view/stern-rehab-cqxbvxin

License: This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development
We are still developing and optimizing this protocol

Created: Mar 10, 2023

Last Modified: Mar 10, 2023

PROTOCOL integer ID:

78531

Adding stocks to Stern REHAB

- When a vial looks sick, first check all other copies (including drobot copy) and if EVERY vial looks sick, consolidate all of the flies in a fresh food vial with moist cotton roll
- 1.1 If not all vials look sick, top off your sick vial with the other healthy copy

1.2	Cotton rolls are located on the back of the REHAB cart underneath the purple REHAB vials
2	Add vial with consolidated flies and all older vials to designated box (organized by species) on REHAB cart next to Stern HF carts
	Flipping REHAB (every Friday)
3	Flip flies from vial in front to a fresh vial with a moist cotton roll
4	Check vials behind the first one. If flies have emerged, flip those into the new vial as well If vial is looking healthy, separate into two vials and wait another week. If both are healthy by the
Ü	end of the next week, ready to add back to Drobot/HF box.
6	Discard oldest vial if it looks like nothing else will emerge
	Adding stocks back after REHAB
7	If vial is healthy and ready to go back, add one vial to its spot in the HF boxes (These are organized by species and in ascending order of stock ID)
	These are organized by species and in ascending order by stock ID

8 Add other vial back to Drobot the next time its species is being flipped (Drobot 1)