

Aug 29, 2024

Purification of BCL2L13-GST

DOI

dx.doi.org/10.17504/protocols.io.rm7vzjj12lx1/v1

Elias Adriaenssens¹

¹Sascha Martens lab, University of Vienna, Max Perutz Labs - Vienna



Elias Adriaenssens

Sascha Martens lab, University of Vienna, Max Perutz Labs - ...

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.rm7vzjj12lx1/v1

Protocol Citation: Elias Adriaenssens 2024. Purification of BCL2L13-GST. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.rm7vzjj12lx1/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 23, 2024

Last Modified: August 29, 2024

Protocol Integer ID: 101126

Keywords: ASAPCRN

Funders Acknowledgement:
Aligning Science Across
Parkinson's (ASAP)
Grant ID: ASAP-000350
Marie Skłodowska-Curie
MSCA Postdoctoral
fellowship
Grant ID: 101062916

Abstract

This protocol details the purification of BCL2L13-GST.

Materials

Lysis buffer:

| A | B |
|--|--------|
| Tris-HCl | 50 mM |
| pH | 7.4 |
| NaCl | 300 mM |
| Triton X-100 | 1% |
| glycerol | 5% |
| MgCl ₂ | 2 mM |
| DTT | 1 mM |
| β-mercaptoethanol | 2mM |
| cOmplete EDTA-free protease inhibitors (Roche) | |
| CIP protease inhibitor (Sigma) | |
| DNase (Sigma) | |

Wash buffer:

| | |
|----------|--------|
| Tris-HCl | 50 mM |
| pH | 7.4 |
| NaCl | 300 mM |
| DTT | 1 mM |




High salt wash buffer:

| | |
|----------|--------|
| Tris-HCl | 50 mM |
| pH | 7.4 |
| NaCl | 700 mM |
| DTT | 1 mM |

SEC buffer:

| | |
|----------|--------|
| Tris-HCl | 25 mM |
| pH | 7.4 |
| NaCl | 300 mM |
| DTT | 1 mM |








- pET-DUET1 vector (available on Addgene)  pETDuet-1 TIM9,10 **addgene Catalog #170280**
- BCL2L13 W275A/I278A (Δ LIR1)(available on Addgene)
- BCL2L13 Y213A/I216A/W275A/I278A (Δ LIR1+2) (available on Addgene)
- BCL2L13 I224A/L227A/W275A/I278A (Δ LIR1+3) (available on Addgene)
- BCL2L13 W275A/I278A/I307A/V310A (Δ LIR1+4) (available on Addgene)
- BCL2L13 I224A/L227A/W275A/I278A/I307A/V310A (Δ LIR1+3+4) (available on Addgene)
- Rosetta pLysS cells (Novagen Cat# 70956-4)
 Rosetta™(DE3)pLysS Competent Cells - Novagen **Merck Catalog #70956-4**
- SORVAL RC6+ centrifuge with an F21S-8x50Y rotor (Thermo Scientific)
- Glutathione Sepharose 4B beads (GE Healthcare)
- 10 kDa cut-off Amicon filter (Merck Millipore)
 Amicon® Ultra Centrifugal Filter, 10 kDa MWCO **Merck MilliporeSigma (Sigma-Aldrich) Catalog #UFC801008**



Purification - BCL2L13-GST


16h

- 1 To purify BCL2L13-GST, fuse the cytosol-exposed domain of BCL2L13 (1-463aa) to a C-terminal GST-tag through cloning into a pET-DUET1 vector (available on Addgene).
- 2 Introduce the point mutants by in vitro mutagenesis to generate
 - BCL2L13 W275A/I278A (Δ LIR1)(available on Addgene),
 - BCL2L13 Y213A/I216A/W275A/I278A (Δ LIR1+2) (available on Addgene),
 - BCL2L13 I224A/L227A/W275A/I278A (Δ LIR1+3) (available on Addgene),
 - BCL2L13 W275A/I278A/I307A/V310A (Δ LIR1+4) (available on Addgene),
 - BCL2L13 I224A/L227A/W275A/I278A/I307A/V310A (Δ LIR1+3+4) (available on Addgene).
- 3 After the transformation of the pET-DUET1 vector encoding BCL2L13-GST wild-type or mutants in *E. coli* Rosetta pLysS cells (Novagen Cat# 70956-4), grow the cells in 2x Tryptone Yeast extract (TY) medium at  37 °C until an OD₆₀₀ of 0.4 and then continued at  18 °C .
- 4 Once the cells reaches an OD₆₀₀ of 0.8, induce the protein expression with  100 micromolar (μ M) isopropyl β -D-1-thiogalactopyranoside (IPTG) for  16:00:00 at  18 °C . Collect the cells by centrifugation and resuspend in lysis buffer.


16h

Lysis buffer:


| A | B |
|--|--------|
| Tris-HCl pH 7.4 | 50 mM |
| NaCl | 300 mM |
| Triton X-100 | 1% |
| Glycerol | 5% |
| MgCl ₂ | 2 mM |
| DTT | 1 mM |
| β -mercaptoethanol | 2mM |
| cOmplete EDTA-free protease inhibitors (Roche) | |
| CIP protease inhibitor (Sigma) | |
| DNase (Sigma) | |

- 5 Sonicate the cell lysates twice for 30 s and clear by centrifugation at  18000 rpm, 4°C, 00:45:00 in a SORVAL RC6+ centrifuge with an F21S-8x50Y rotor (Thermo Scientific).



45m

5.1 Sonicate the cell lysates for  00:00:30 (1/2).

30s

5.2 Sonicate the cell lysates for  00:00:30 (2/2).

30s

6 Collect the supernatant and incubate with pre-equilibrated Glutathione Sepharose 4B beads (GE Healthcare) for  02:00:00 at  4 °C with gentle shaking to bind BCL2L13-GST.

2h



7 Centrifuge the samples to pellet the beads and remove the unbound lysate.



8 Wash the beads twice with wash buffer, once with high salt wash buffer, and two more times with wash buffer.






Wash buffer:

| A | B |
|-----------------|--------|
| Tris-HCl pH 7.4 | 50 mM |
| NaCl | 300 mM |
| DTT | 1 mM |

High salt wash buffer:

| A | B |
|-----------------|--------|
| Tris-HCl pH 7.4 | 50 mM |
| NaCl | 700 mM |
| DTT | 1 mM |


9 Incubate the beads  Overnight with  4 mL of [M] 50 millimolar (mM) reduced glutathione dissolves in wash buffer at  4 °C , to elute BCL2L13-GST from the beads.

2h



10 To collect the supernatant, collect the beads by centrifugation.



11 Wash the beads twice with  4 mL of wash buffer, and collect the supernatant.





- 12 Pool the supernatant fractions, filter through a 0.45 µm syringe filter, concentrate with 10 kDa cut-off Amicon filter (Merck Millipore), and load onto a pre-equilibrated Superdex 200 Increase 10/300 GL column (Cytiva).
- 13 Elute the proteins with SEC buffer.

SEC buffer:

| A | B |
|-----------------|--------|
| Tris-HCl pH 7.4 | 25 mM |
| NaCl | 300 mM |
| DTT | 1 mM |

- 14 Analyze fractions by SDS-PAGE and Coomassie staining. Pool fractions containing purified BCL2L13-GST.
- 15 After concentrating the purified protein, aliquote the protein and snap-frozen in liquid nitrogen.



Note

Store the proteins at  -80 °C .