



Aug 16, 2022

Washing Protocol for Intact Proteoform MALDI on Human Pancreas

Kevin J. Zemaitis¹, Dusan Velickovic¹, Ljiljana.PasaTolic¹

¹Pacific Northwest National Laboratory



dx.doi.org/10.17504/protocols.io.5qpvobe17l4o/v1

Human BioMolecular Atlas Program (HuBMAP) Method Development Community PNNL-TTD

Kevin J. Zemaitis

ABSTRACT

Scope:

A detailed protocol entailing the sample washing protocols developed for human pancreas, this includes several timed washing steps for fixing the tissue, removing small molecule interferents, and desalting the tissue surface optimized for use with a Spectroglyph EP-MALDI-2 source mounted on a custom Thermo Scientific UHMR Q Exactive HF Orbitrap.

Expected Outcomes:

Human pancreas slides ready for pre-extraction and/or matrix application.

DOI

dx.doi.org/10.17504/protocols.io.5qpvobe17l4o/v1

PROTOCOL CITATION

Kevin J. Zemaitis, Dusan Velickovic, Ljiljana.PasaTolic 2022. Washing Protocol for Intact Proteoform MALDI on Human Pancreas. **protocols.io** https://protocols.io/view/washing-protocol-for-intact-proteoform-maldi-on-huce2ctgaw

FUNDERS ACKNOWLEDGEMENT

National Institutes of Health (NIH) Common Fund, Human Biomolecular Atlas Program (HuBMAP)

Grant ID: UG3CA256959-01



1

LICENSE

This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 08, 2022

LAST MODIFIED

Aug 16, 2022

PROTOCOL INTEGER ID

68388

PARENT PROTOCOLS

In steps of

Tissue Preparation for Intact Proteoform MALDI-MSI on Human Tissue

MATERIALS TEXT

Chemicals:

- Ethanol (200 Proof)
- Nanopure water (LC-MS Grade)
- Chloroform (LC-MS Grade)
- Acetic acid (Glacial, LC-MS Grade)
- Trifluoroacetic acid (LC-MS Grade)

Equipment:

- Regulated nitrogen gas supply with nozzle
- Coplin jars

SAFETY WARNINGS

All steps of this protocol working with solvents should be performed within a fume hood as to minimize exposure to fumes from volatile organic solvents.

BEFORE STARTING

Prepare 100 mL of all necessary solvents for the tissue washing steps and clean Coplin jars prior to use.

Preparation

While the tissue is within the vacuum desiccator, pour all necessary solvents and solutions into the Coplin jar. These have a volume of approximately **75 mL** and to reduce the variability in extraction fill the jars consistently.

protocols.io

2 Prior to submerging the tissue section in solvent and solutions, ensure that the timers are set beforehand. While three timers are not necessary, it vastly improves the process.

3 🐧

If any metadata is present on the slide written in permanent marker remove it at this time. Take a photo and note within a notebook prior to removal, this is necessary as this will contaminate the Coplin jars within all washes.

Tissue fixation wash 3m

- 4 Submerge the tissue section within a Coplin jar filled with a solution of **70% ethanol** for © **00:00:30**
- 5 Immediately move the tissue section with tweezers to then next coplin jar containing 100% ethanol for © 00:00:30

Lipid depletion wash 3m

- After completion of the tissue fixation steps, immediately move the tissue section with tweezers to then next Coplin jar containing **Carnoy's solution** for **© 00:00:20**
 - **Carnoy's solution** is a 6:3:1 volumetric solution of ethanol: chloroform: glacial acetic acid, while other solutions could be put within plastic staining jars a glass Coplin jar is highly recommended for this and other polar aprotic solvents.
 - 6.2

Note: dipping, agitation, and any movement of the slide within the Coplin jar will dramatically change outcomes of all washes. Tissue is prone to detach from the surface within this step.

7 Following exposure to Carnoy's solution, submerge the slide with **100% ethanol** for **© 00:00:30**

30s

Desalting wash 3m

8 After completion of the lipid depletion washes, submerge the tissue sections within **nanopure** water with 0.2% trifluoroacetic (TFA) acid for 5 dips.

8.1

Note: dipping, agitation, and any movement of the slide within the Coplin jar will dramatically change outcomes of all washes. Tissue is prone to detach from the surface within this step.

30s

- After exposure to water, submerge the slide within 100% ethanol for © 00:00:30
 - 9.1 Ensure this is a new solution of **100% ethanol** not used within previous steps.

Drying tissue 3m

- After the completion of all steps, dry the tissue for **© 00:00:30** under a diffuse stream of nitrogen, take care to ensure that liquid does not pool within the surface. The tissue section is now prepared for the next steps within the tissue preparation protocol.
- 11 Discard all solvents and solutions properly, and best practice is to use fresh solvents and solutions for subsequent washes.