



Version 2

Apr 23, 2021

Making electrocompetent *Agrobacterium tumefaciens*

V.2

Johannes Wolfram JWD Debler¹¹Curtin University

1

Works for me

This protocol is published without a DOI.

Johannes Wolfram Debler
Curtin University

ABSTRACT

This protocol yields about 5 ml of electrocompetent *Agrobacterium tumefaciens* cells, aliquotted into 80 ul you get 50-60 tubes out of it.

Rifampicin is dark red and the colour will change over time, it does however NOT interfere with OD600 measurements, so you can use LB to blank.

PROTOCOL CITATION

Johannes Wolfram JWD Debler 2021. Making electrocompetent *Agrobacterium tumefaciens*. **protocols.io**
<https://protocols.io/view/making-electrocompetent-agrobacterium-tumefaciens-bufpntmn>
Version created by Johannes Wolfram Debler

KEYWORDS

Agrobacterium tumefaciens, electrocompetent

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Apr 23, 2021

LAST MODIFIED

Apr 23, 2021

PROTOCOL INTEGER ID

49359

Things you need to prepare

1

- LB medium (for 1 litre combine 10 g Tryptone, 5 g Yeast Extrad and 10 g NaCl, adjust pH to 7.0 and autoclave)
- 10% glycerol solution (ice cold)
- H₂O (ice cold)
- liquid nitrogen
- 50 - 60 x 1.5 ml tubes (pre-chill in fridge)
- 10 ml pipette tips
- 8 x 50 ml tubes (pre-chill in fridge)
- Your favourite strain of *Agrobacterium tumefaciens* grown on a plate

Depending on your strain of *Agrobacterium* you'll need the following antibiotics:

- **AGL1**: rifampicin 100 ng/ul
- **LBA 4404**: rifampicin 100 ng/ul, streptomycin 100 ng/ul --> grows slow and clumps
- **GV3101::pMP90**: rifampicin 100 ng/ul, gentamicin 25 ng/ul --> grows fastest of the three

Day 1

- Inoculate **10 ml LB** (containing the appropriate antibiotics) with **1 colony of *Agrobacterium tumefaciens*** in the morning and incubate at **28°C at 250 rpm** for the rest of the day.

🌡 28 °C

🧴 3 mL LB

🕒 05:00:00

- In the evening inoculate **3 x 150 ml LB** (containing the appropriate antibiotics) with **3 ml each of the culture** grown throughout the day.

🌡 28 °C

🧴 450 mL LB

- Incubate at **28°C at 250 rpm** until **OD600 reaches 0.8 - 1.0** (This usually takes overnight, depending on Agro strain im might take longer).

🌡 28 °C

🕒 30:00:00

Day 2

- Check OD600, and combine all 3 cultures once **OD600 reaches 0.8-1.0**

- Chill culture by **putting it on ice for 15 minutes and swirling it regularly.**

From here on everything needs to happen on ice or at 4°C!

🕒 00:15:00

🌡 4 °C

- Distribute culture into **8 x 50 ml tubes** and pellet by centrifugation at **3000 - 4000 x g for 20-30 minutes at 4°C.**

🕒 00:30:00

🌡 4 °C

- Discard supernatant and resuspend each pellet in **10 ml ice cold H2O.**

🧴 10 mL H2O

- Adjust volume to **50 ml each with ice cold H2O.**

- Pellet cells by centrifugation at **3000 - 4000 x g for 20-30 minutes at 4°C.**

🕒 00:30:00

🌡 4 °C

- Discard supernatant and resuspend each pellet in **10 ml ice cold H2O.**

 **10 mL H2O**

12 Combine solutions into **4 x 50 ml** tubes and adjust to **50 ml with ice cold H2O**.

13 Pellet cells by centrifugation at **3000 - 4000 x g for 20-30 minutes at 4°C**.

 **00:30:00**


 **4 °C**

14 Discard supernatant and resuspend each pellet in **10 ml ice cold H2O**.

 **10 mL H2O**

15 Combine solutions into **2 x 50 ml** tubes and adjust to **50 ml with ice cold H2O**.

16 Pellet cells by centrifugation at **3000 - 4000 x g for 20-30 minutes at 4°C**.

 **00:30:00**

 **4 °C**

17 Discard supernatant and resuspend each pellet in **10 ml ice cold H2O**.

 **10 mL H2O**

18 Combine solutions into **1 x 50 ml** tubes and adjust to **50 ml with ice cold H2O**.

19 Pellet cells by centrifugation at **3000 - 4000 x g for 20-30 minutes at 4°C**. (Use a counterbalance tube)

 **00:30:00**

 **4 °C**

20 Discard supernatant and resuspend pellet in **4.5 ml ice-cold 10% glycerol**.

 **4.5 mL 10% glycerol**

21 Dispense **80 ul aliquots** into pre-chilled 1.5 ml tubes and flash freeze in liquid nitrogen.

22 Store aliquots at **-80°C** until use.