



Oct 09, 2022

## SDS-PAGE

Ana Belem García González<sup>1</sup>, Georgina Diego<sup>1</sup>,  
Irán Alessandra Chaparro Rodríguez<sup>1</sup>, Jair Alexis Gardea Sáenz<sup>1</sup>

<sup>1</sup>Tecnologico de Monterrey Campus Chihuahua

1 Works for me

Share

This protocol is published without a DOI.

**D** Denisse Chacón Ramírez

### ABSTRACT

This protocol shows the steps carried out by team Tec-Chihuahua to perform SDS-PAGE

### PROTOCOL CITATION

Ana Belem García González, Georgina Diego, Irán Alessandra Chaparro Rodríguez, Jair Alexis Gardea Sáenz 2022. SDS-PAGE. **protocols.io**  
<https://protocols.io/view/sds-page-cezqtf5w>



### LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

### CREATED

Aug 08, 2022

### LAST MODIFIED

Oct 09, 2022

### PROTOCOL INTEGER ID

68368

### MATERIALS TEXT

#### SOLUTIONS:

- **Acrylamide/Bis 30%**
  - 29.2 g of Acrylamide.
  - 0.8 g of Bis Acrylamide
  - Dissolve in 50 mL of distilled water in constant agitation, gauge to a volume of 100 mL.
  - Store at 4 °C degrees while protected from light.
- **Tris-HCl/ SDS 4X pH 8.8 (1.5 M Tris-HCl, 0.4% SDS)**

- 18.17 g of TRIZMA® base
  - 0.4 g of SDS
  - Dissolve in 80 mL of distilled water
  - Adjust to a pH de 8.8 with HCl and gauge to a 100 mL
  - Filtrate the solution
  - Store at 4 °C degrees
- 
- **Tris-HCl/ SDS 4X pH 6.8 (0.5 M Tris-HCl, 0.4% SDS)**
  - 6.06 g de TRIZMA® base
  - 0.4 g de SDS
  - Dissolve en 80 mL of distilled water.
  - Adjust to a pH de 6.8 with HCl and gauge to a 100 mL
  - Filtrate the solution
  - Store at 4 °C degrees
- 
- **10% Ammonium Persulfate**
  - 500 mg of APS (Ammonium Persulfate)
  - Dissolve en 5 mL of distilled ultra pure water
  - Store at -20 °C for up to two weeks
  - Store at 4 °C degrees for one use only
- 
- **2X Loading Buffer**
  - 2 mL of glycerol
  - 400 uL of Mercapto
  - 0.02 g of bromophenol blue
  - 0.4 g of SDS
  - 2.5 mL of 8.6 pH buffer
- 
- **1X Running Buffer**
  - 14 g of glycine.
  - 3 g of Trizma base.
  - 1 g of SDS.
  - Gauge to a 1 L.
- 
- **Staining Solution**
  - 1.25 g of brilliant blue R in 250 mL of methanol
  - 200 mL of distilled water.
  - 50 mL of Glacial Acetic Acid
  - Store in room temperature while protected from light
- 
- **Destaining Solution**
  - 250 mL of methanol
  - 62.5 mL of Glacial Acetic Acid
  - 312.5 mL of distilled water.
  - Store in room temperature

**Usage:** Pour solution into a container and allow to stir until the SDS gel is clear or electrophoresis bands are visible.

- 1
- Clean the components of the electrophoresis camera with 70% ethanol and gauze 's.
- 2
- Assemble the chamber and check that there are no leaks by pouring distilled water between the glasses.
- 3
- Prepare polyacrylamide gels:








A	B
Reagents	2 minigels
Distilled water	3.4 mL
Acrylamide/Bis 30%	4 mL
Tris-HCl/SDS 4X pH 8.8 (1.5 M Tris-HCl, 0.4% SDS)	2.5 mL
ASP 10%	100 µL
TEMED	4 µL

Separation gel (12%)

- 4
- Pour the solution between the glasses with a 1 mL micropipette, leaving a space of 1.5 cm<sup>30m</sup> for the concentrating gel. To level, distilled water is added and allowed to settle for ⌚00:30:00 or until a line is seen between the gel and the water.

A	B
Reagents	2 minigels
Distilled water	2.7 mL
Acrylamide/ Bis 30%	1 mL
Tris-HCl/ SDS 4X pH 8.8 (1.5 M Tris-HCl, 0.4% SDS)	1.3 mL
ASP 10%	50 µL
TEMED	4 µL

■ Concentrating gel (6%)

- 5 Pour solution onto separating gel using a 1 mL micropipette.
- 6  00:30:00  00:30:00 Insert the comb (carefully avoiding the formation of bubbles) and leave to solidify for  00:30:00 1h 30m
- 7 When the gels are polymerized, prepare an electrophoresis chamber with 1X running buffer until it covers the gels and it reaches the line of two gels.
- 8 Sample PTake the pellets contained in Eppendorf tubes
  - Take the pellets contained in Eppendorf tubes
  - Add  300 µL of 1X loading buffer and resuspend the pellet.reparation:
- 9 To denature proteins, heat samples in boiling water for 5 min.
- 10 Load gels with the hot sample.
- 11 In the first well add  7 µL of the molecular weight marker.
- 12 Once the samples are loaded, run the gel at 80 V for  00:20:00 and then at 180 V for  00:45:00 1h 5m
- 13 Turn off the camera and disarm it.
- 14 Remove the gels from the glasses and place in a container with staining solution. Leave stirring for one hour.

- 15 Remove staining solution after one hour, add destaining solution and leave stirring for 40m  
🕒 00:20:00 Change the destaining solution and leave again for 🕒 00:20:00
- 16 Leave stirring until the gel is transparent.
- 17 Analyze the gels.