



Dec 02, 2021

Primer TE (10mM Tris, 0.1mM EDTA, pH8.0)

Jacquelina.Woods¹¹FDA

1

dx.doi.org/10.17504/protocols.io.byccpssw**GenomeTrakr**Tech. support email: genomeTrakr@fda.hhs.gov

Jessica Jones

US Food and Drug Administration

This method was developed at the FDA's Center for Food Safety and Applied Nutrition for GenomeTrakr's pandemic response project, monitoring SARS-CoV-2 variants in wastewater. Protocols developed for this project cover wastewater collection, concentration, RNA extraction, RT-qPCR detection, library prep, genome sequencing, quality control checks, and data submission to NCBI. This method provides a reagent formula required in the extraction of RNA from viral concentrates using the RNeasy and Zymo kits.

DOI

dx.doi.org/10.17504/protocols.io.byccpssw

Jacquelina.Woods 2021. Primer TE (10mM Tris, 0.1mM EDTA, pH8.0).

protocols.io<https://dx.doi.org/10.17504/protocols.io.byccpssw>

protocol ,

Sep 17, 2021

Dec 02, 2021

53348

In steps of

[RNA Extraction from Wastewater Concentrates Using RNeasy and Zymo Kits](#)[RNA Extraction from Wastewater Concentrates Using RNeasy and Zymo Kits](#)[Reagent Mixes for RT-qPCR Detection of SARS-CoV-2 from Wastewater](#)

 [Tris \(1 M\) pH 8.0 RNase-free](#) **Thermo Fisher**

Scientific Catalog #AM9856 Step 1.1

 [EDTA \(0.5 M\) pH 8.0 RNase-free](#) **Thermo Fisher**

Scientific Catalog #AM9261 Step 1.2

 [Nuclease-Free Water](#) **Thermo Fisher**

Scientific Catalog #AM9937 Step 1.3

1 Mix components together.

1.1 100 µL

 [Tris \(1 M\) pH 8.0 RNase-free](#) **Thermo Fisher**

Scientific Catalog #AM9856

, or equivalent.

1.2 20 µL

 [EDTA \(0.5 M\) pH 8.0 RNase-free](#) **Thermo Fisher**

Scientific Catalog #AM9261

, or equivalent.

1.3 9.88 mL

 [Nuclease-Free Water](#) **Thermo Fisher**

Scientific Catalog #AM9937

, or equivalent.

2 Store at **Room temperature** .