




VERSION 1

NOV 08, 2022



SHARE

WORKS FOR ME 1

Extraction and ONT MinLibrary Preparation of uHMW gDNA V.1

This protocol is published without a DOI.

Kaylee S. Herzog¹, jfauver¹

¹University of Nebraska Medical Center



Kaylee S. Herzog

COMMENTS 0

ABSTRACT

This custom protocol optimizes extraction, purification, and Oxford Nanopore Technologies (ONT) MinION library preparation for ultra-high molecular weight genomic DNA (uHMW gDNA) from parasitic nematodes. It can be used effectively with both low-input samples (e.g., a single adult hookworm) and high-input samples (e.g., a large chunk of tissue from an *Ascaris* sp. adult).

Protocols on which this workflow is based:

- [Zymo® Quick-DNA™ Magbead Plus Kit protocol](#)
- [Oxford Nanopore Technologies® SQK-LSK-109 gDNA Ligation Sequencing protocol](#)
- Zymo® DNA Clean & Concentrator™ Magbead Kit protocol (best-testing phase only)

PROTOCOL CITATION

Kaylee S. Herzog, jfauver 2022. Extraction and ONT MinLibrary Preparation of uHMW gDNA.
protocols.io
<https://protocols.io/view/extraction-and-ont-minlibrary-preparation-of-uhmw-ciylhft6>

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Nov 07, 2022

LAST MODIFIED

Nov 08, 2022

PROTOCOL INTEGER ID

72425

BEFORE STARTING

- Add 1,040 µL Zymo Proteinase K Storage Buffer to each tube of Zymo Proteinase K (20 mg) prior to use. The final concentration of Proteinase K is ~20 mg/ml. Store resuspended Proteinase K at -20°C after mixing.
- For best results, allow AMPure XP beads (stored at 4°C) to come to RT prior to use.

Part 1: Ultra-HWM gDNA extraction | Zymo Quick-DNA HWM MagBeads

1 Set dry bath to 55 °C

2 For each sample, add the following to a clean 1.5 mL microcentrifuge tube to create a master mix:

95 µL	Zymo DNA Elution Buffer Zymo Research Catalog #D3004-4-1
95 µL	Zymo Biofluid & Solid Tissue Buffer Zymo Research Catalog #D4081-3-25
10 µL	Zymo Proteinase K Zymo Research Catalog #D3001-2-20

2.1 Vortex the master mix gently to mix, then spin down and keep on ice

3 Using a new pipette tip or sterilized forceps, add one whole worm (or a piece of tissue) directly from tissue preservative to the bottom of a clean 1.5 mL microcentrifuge tube,

Note

Transfer as little tissue preservative liquid as possible to the new tube during this process

4 Use a new **NEB Monarch Pestle NEB Catalog #T3002-1** to grind and crush the tissue in the tube. Keep the pestle in the tube

- 5 Add 200 μ L master mix (prepared in Part 1 Step 1) to each tube containing tissue and pestle
- 6 Continue using the pestle to grind the tissue within the master mix until fully homogenized. Remove the pestle, being careful to keep any tissue in the tube by wiping the pestle on the tube edges as it is removed
- 7 Close the tube and mix by inverting and flicking gently, then spin down briefly to recollect tissue and liquids
- 8 Incubate sample in dry bath at 55°C for 02:30:00 or until tissue solubilizes. During incubation, flick tube every 00:20:00 to agitate tissues, then briefly spin down to recollect liquids and replace tube in dry bath

2h 50m



Note

If a very large amount of input tissue was used: It is likely there will still be visible tissue even after hours of lysis. If so, centrifuge the sample for 00:01:00 at 10000 $\times g$ or greater to pellet debris, then pipette all liquids into a clean 1.5 μ L microcentrifuge tube. (The majority of gDNA will be contained in the layer of liquid just above the pellet, so pipette carefully to get as much liquid as possible without disturbing the debris.) Discard the tube contain the pelleted debris and use the supernatant to move forward for Part 2.

2h 20m

Part 2: Ultra-HWM gDNA purification | Zymo Quick-DNA HWM MagBea

- 9 Set dry bath to 37°C
- 10 Add 400 μ L Zymo Quick-DNA™ MagBinding Buffer Zymo Research Catalog #D4077-1-150 to each sample
- 11 Flick tubes to mix, then spin down briefly to recollect liquids

- 12 Add  33 µL  Zymo MagBinding Beads **Zymo Research Catalog #D4100-2-6** to each sample

Note

MagBinding Beads settle quickly, so ensure beads are kept in suspension while dispensing by vortexing each time before beads are added to a sample



- 13 To ensure DNA binds to beads, mix on a rotator mixer at a low speed for 40 min at room temperature. Spin down briefly before proceeding with the next step



- 14 Set sample tubes on a magnetic stand until beads have separated from solution, then remove and discard the supernatant. Remove sample tubes from the magnetic stand.

Note

Some beads may adhere to the sides of the tube. When removing supernatant, aspirate slowly to allow these beads to be pulled to the magnet as the liquid level is lowered.

- 15 Add  500 µL  Zymo Quick-DNA™ MagBinding Buffer **Zymo Research Catalog #D4077-1-150** to each sample



- 16 Flick to mix initially, then mix on a rotator mixer at a low speed for 20 min. Spin down briefly before proceeding with the next step



- 17 Set sample tubes on a magnetic stand until beads have separated from the solution, then remove and discard the supernatant. Remove sample tubes from the magnetic stand

Note

Some beads may adhere to the sides of the tube. When removing supernatant, aspirate slowly to allow these beads to be pulled to the magnet as the liquid level is lowered.

18 Add  500 µL  Zymo DNA Pre-Wash Buffer **Zymo Research Catalog #D3004-5-250** to each sample

19 Flick to mix, then spin down briefly

20 Set sample tubes on a magnetic stand until beads have separated from solution, then remove and discard the supernatant. Remove sample tubes from the magnetic stand

Note

Some beads may adhere to the sides of the tube. When removing supernatant, aspirate slowly to allow these beads to be pulled to the magnet as the liquid level is lowered.

21 Add  900 µL  Zymo g-DNA Wash Buffer **Zymo Research Catalog #D3004-2-200** to each sample

22 Flick to mix, then spin down briefly

23 Transfer the entire sample (all liquid and beads) to a new 1.5 mL microcentrifuge tube



24 Set samples (now in new tubes) on a magnetic stand until beads have separated from the solution, then remove and discard the supernatant. Remove sample tubes from the magnetic stand

Note

Some beads may adhere to the sides of the tube. When removing supernatant, aspirate slowly to allow these beads to be pulled to the magnet as the liquid level is lowered

25 Add  900 µL  Zymo g-DNA Wash Buffer **Zymo Research Catalog #D3004-2-200** to each sample

26 Flick to mix, then spin down briefly

27 Transfer the entire sample (all liquid and beads) to a new 1.5 mL microcentrifuge tube




28 Set samples (now in new tubes) on a magnetic stand until beads have separated from the solution, then remove and discard the supernatant. Remove sample tubes from the magnetic stand

Note

Some beads may adhere to the sides of the tube. When removing supernatant, aspirate slowly to allow these beads to be pulled to the magnet as the liquid level is lowered



28.1 Use a P10 pipette to remove any residual liquid from the bottom of the tube

29 Air dry the beads for up to  00:20:00 and proceed to next step once beads are dry, but not over-dry

20m

Note

It may take less time for the beads to dry, so check them often during this process! Beads will change in appearance from glossy black when still wet to a matte black/brown when fully dry. Over drying the beads may result in lower DNA recovery.

30 Add  50 µL  Zymo DNA Elution Buffer **Zymo Research Catalog #D3004-4-50** to each sample and flick gently several times to mix. Spin down briefly

31 Incubate in dry bath at 37°C for 02:00:00 with occasional gentle flicking and spinning down throughout incubation

2h

32 Incubate on bench top at Room temperature overnight.

33 After overnight incubation, set tubes on a magnetic stand until beads have separated from solution, then move the supernatant (now containing eluted DNA) to a new tube

Note

The eluted DNA can be used immediately or stored at 4°C or -20°C for future use

33.1 Re-suspend beads in 20 μL of Nuclease-free Water Contributed by users in case there is no (or not enough) gDNA in final elution

34 Use 1 μL of final elution to quantify extraction via Qubit analysis

35 Use 1 μL of final elution to assess fragment size distribution via TapeStation

*

1h

Part 3: DNA repair and end-prep | Zymo Clean & Concentrator, ONT L

36 Set dry bath to 65°C

37 Defrost the needed NEB DNA and End Repair reagents on ice (see Part 3 Step 38)

38 For each sample, add the following to a clean 0.2 mL PCR tube to create a master mix, pipetting 10–20 times between each addition to mix:

3.5 µL	NEBNext® FFPE DNA Repair Buffer	New England Biolabs Catalog #E7180S
2 µL	NEBNext FFPE DNA Repair Mix - 96 rxns	New England Biolabs Catalog #M6630L
3.5 µL	NEBNext Ultra II End Prep Reaction Buffer	New England Biolabs Catalog #E7647
3 µL	NEBNext Ultra II End Prep Enzyme Mix	New England Biolabs Catalog #E7646

38.1 Keep master mix on ice

39 Add 12 µL of master mix (prepared in Part 3 Step 38) from the PCR tube directly into each 1.5 mL microcentrifuge tube containing extracted & purified uHWM gDNA (from Part 2). Mix all components by gently flicking, and spin tubes down to recollect liquids

40 Incubate samples at Room temperature for 00:10:00

10m



41 Incubate samples at 65 °C for 00:10:00

10m



42 Add 4 volumes of Zymo DNA MagBinding Buffer Zymo Research Catalog #D4012-1-50 to each sample and mix well by flicking and inverting

Note



Example for calculating 4 volumes: If input is 49 µL gDNA, add 196 µL DNA MagBinding Buffer

43 Spin samples down briefly and add 20 µL

Zymo MagBinding Beads Zymo Research Catalog #D4100-5-2

Note



MagBinding Beads settle quickly, so ensure beads are kept in suspension while dispensing by vortexing each time before beads are added to a sample

44 Mix samples on rotating mixer at a low speed at  Room temperature for  00:20:00


20m



45 Briefly spin down samples and pellet on a magnetic stand (1–2 min) until the supernatant is clear and colorless. With the tubes still on the magnet, pipette off and discard the supernatant

46 Add  500 μ L  Zymo DNA Wash Buffer **Zymo Research Catalog #D4003-2-24** and then remove from magnetic stand, and mix well by flicking and inverting

47 Briefly spin samples down briefly and transfer to magnetic stand to allow beads to pellet until solution is clear (1–2 min). With the tubes still on the magnet, pipette off and discard the supernatant


48 Air dry the beads for  00:10:00

10m

Note

MagBinding Beads utilize a different chemistry than SPRI beads (e.g., AMPure XP beads) so there is not the same risk of over-drying. It is important for optimal elution that the residual buffer is completely removed/evaporated from the beads

49 Add  50 μ L  Zymo DNA Elution Buffer **Zymo Research Catalog #D3004-4-10**


50 Manually agitate samples by gently flicking (and occasionally spinning down) for  00:10:00

10m















Note

This volume is too small to be able to use the rotator mixer effectively at this step. Manual agitation is necessary







- 51 Briefly spin samples down and pellet the beads on a magnet until the eluate is clear and colorless (1–2 min)
- 52 Remove and retain the 50 µL of eluate (containing repaired & end-prepped DNA) to a new clean 1.5 mL microcentrifuge tube
- 53 Use  1 µL of final elution to quantify via Qubit assay

Part 4: Adaptor ligation and clean up | ONT Ligation Sequencing & NEB

- 54 Set dry bath to  37 °C
- 55 Remove  AMPure XP Beads **Beckman Coulter Catalog #A63880** from storage at  4 °C and allow them to come to  Room temperature
- 56 Spin down  ONT Adaptor Mix (AMX) **Oxford Nanopore Technologies** and  Quick T4 DNA Ligase **New England Biolabs Catalog #E7180S** and place on ice
- 57 Thaw  ONT Ligation Buffer (LNB) **Oxford Nanopore Technologies** at  Room temperature, spin down, and mix by pipetting. Place on ice immediately after thawing and mixing
- 58 Thaw  Elution Buffer (EB) **Oxford Nanopore Technologies** at  Room temperature, vortex to mix, spin down, and place on ice





59 Thaw one tube each of  Short Fragment Buffer (SFB) **Oxford Nanopore Technologies** and  Long Fragment Buffer (LFB) **Oxford Nanopore Technologies** at room temperature, vortex to mix, spin down, and place on ice

60 For each sample, add the following, in order, to a clean 1.5 mL microcentrifuge tube pipetting 10–20 times between each addition to mix:


 25 µL	 ONT Ligation Buffer (LNB) Oxford Nanopore Technologies
 10 µL	 Quick T4 DNA Ligase New England Biolabs Catalog #E7180S
 5 µL	 ONT Adaptor Mix (AMX) Oxford Nanopore Technologies

60.1 Keep master mix on ice after mixing

61 For each sample, prepare 1:3 SFB:LFB titrated wash mix by adding the following to a 1.5 mL microcentrifuge tube and vortex to mix:

 125 µL	 Short Fragment Buffer (SFB) Oxford Nanopore Technologies
 375 µL	 Long Fragment Buffer (LFB) Oxford Nanopore Technologies

61.1 Keep titrated wash mix on ice after vortexing


62 Pipette  40 µL of master mix (prepared in Part 4 Step 60) directly into entire volume of repaired and end-prepped gDNA from Part 3. Mix all components by gently flicking and spin tube down to recollect liquids

63 Incubate the reaction  00:15:00 at  Room temperature

15m



Note

If you have omitted the bead-based purification steps from the second half of Part 3, do not incubate the reaction for longer than  00:10:00



64 Resuspend  AMPure XP Beads **Beckman Coulter Catalog #A63880** by vortexing and add 0.4X volume resuspended beads to each sample, then flick to mix

Note

AMPure XP Beads settle quickly, so ensure beads are kept in suspension while dispensing by vortexing each time before beads are added to a sample

Note

Example for calculating 0.4X volume: If input is 89 μL (after adding master mix), add 35.6 μL AMPure XP Beads



65 Mix on a rotator mixer at a low speed for  01:00:00 at  Room temperature

1h

66 Spin down the sample and pellet on a magnetic stand. Keeping the tube on the stand, pipette off and discard the supernatant

67 Wash the beads by adding 250 μL 1:3 SFB:LFB titrated wash mix (prepared in Part 4 Step 61). Flick the beads to resuspend, spin down, then return to the magnetic rack and allow the beads to pellet. Remove the supernatant using a pipette and discard


68 Wash the beads by adding 250 μL 1:3 SFB:LFB titrated wash mix (prepared in Part 4 Step 61). Flick the beads to resuspend, spin down, then return to the magnetic rack and allow the beads to pellet. Remove the supernatant using a pipette and discard

69 Spin down the beads and place them back on the magnetic rack. Use a P10 pipette to pipette off any residual liquid and allow beads to air-dry for  00:00:30 to  00:02:00

2m 30s



Note

Do not allow the pellet to dry to the point of cracking! Over-drying beads will result in reduced yields

70 Remove the tube from the magnetic stand and resuspend the pellet in  15 μL

71



Briefly spin down and incubate in dry bath at  37 °C for  02:00:00 with occasional gentle flicking and spinning down throughout incubation

2h

Note


For HMW & uHMW gDNA, incubation at 37°C for longer times can improve the recovery of long fragments

72




Incubate on the bench top at  Room temperature overnight

73


1. After overnight incubation, pellet the beads on a magnet until the eluate is clear and colorless (at least  00:01:00)

1m


74

Remove and retain the  15 µL of eluate (containing the prepared library) to a new clean 1.5 mL microcentrifuge tube

75

Use  1 µL of final elution to quantify library via Qubit analysis

Note

Note: For same-day or near-future sequencing, store the prepared library on ice (or at  4 °C) until ready to be loaded onto a flow cell. Otherwise, store libraries at -20°C