

Aug 10, 2020

Total RNA and DNA in microalgae: The extraction and quantification

Yingyu YY Hu¹, Zoe V Finkel¹¹Dalhousie University**1** Works for me This protocol is published without a DOI.

Marine Microbial Macroecology Lab

Tech. support email: ruby.hu@dal.caYingyu Hu
Dalhousie University

ABSTRACT

This assay protocol for total RNA and DNA is adapted from Berdalet E. et al. and Liefer J.D. et al.

RNA and DNA are extracted from microalgae samples and then quantified by fluorochrome SYBR Green II.



Berdalet E, Roldán C, Olivar MP, Lysnes K. Quantifying RNA and DNA in planktonic organisms with SYBR Green II and nucleases. Part A. Optimisation of the assay. *Scientia Marina*.
<https://doi.org/10.3989/scimar.2005.69n11>



Berdalet E, Roldán C, Olivar MP. Quantifying RNA and DNA in planktonic organisms with SYBR Green II and nucleases. Part B. Quantification in natural samples. *Scientia Marina*.
<https://doi.org/10.3989/scimar.2005.69n117>



Liefer JD, Garg A, Fyfe AH, Irwin AJ, Benner I, Brown CM, Follows MJ, Omta AW and Finkel ZV. The Macromolecular Basis of Phytoplankton C:N:P Under Nitrogen Starvation. *Front. Microbiol.*
<https://doi.org/10.3389/fmicb.2019.00763>

PROTOCOL CITATION

Yingyu YY Hu, Zoe V Finkel 2020. Total RNA and DNA in microalgae: The extraction and quantification.
protocols.io
<https://protocols.io/view/total-rna-and-dna-in-microalgae-the-extraction-and-bbnwimfe>

KEYWORDS

RNA, DNA, SYBR Green II, DNase, RNase, microalgae, fluorescence

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jan 22, 2020

LAST MODIFIED

Aug 10, 2020

PROTOCOL INTEGER ID

32182

GUIDELINES

The requirement of microalgae sample is of no higher than 50 ug protein (or 100 ug carbon) per sample, considering the optimized linear range of standard curve for DNA is 0~100 ng/mL and RNA is 0~300 ng/mL.


Microalgae samples are collected on polycarbonate filter and kept frozen under -80 °C.


MATERIALS

NAME	CATALOG #	VENDOR
Nuclease decontamination solution	11-05-03-01	IDT
SYBR™ Green II RNA Gel Stain (10000X concentrate in DMSO stored at -20 °C)	S7564	ThermoFisher
Tris solution (1M pH 8 RNase-free)	AM9851	ThermoFisher
EDTA buffer (0.5M DNase/RNase free)	EDT333.100	Bioshop
N-Lauroylsarcosine sodium salt solution (20% RNase/DNase free)	L744-50mL	Sigma Aldrich
DEOXYRIBONUCLEASE1 RNase and Protease Free	DRB002.10	Bioshop
Ribonuclease A from bovine pancreas	R6513-50MG	Sigma Aldrich
Deoxyribonucleic acid from calf thymus	D4522-1MG	Sigma Aldrich
E. coli Total RNA	AM7940	ThermoFisher
Magnesium chloride solution	63069-100ML	Sigma Aldrich
Calcium chloride solution	21115-100ML	Sigma Aldrich
UltraPure™ DNase/RNase-Free Distilled Water	10977023	ThermoFisher

MATERIALS TEXT

 FreeZone® 2.5 L Benchtop Freeze Dryers
Labconco® 700202000

 A2 Biosafety Cabinet
4' REDISHIP Purifier Logic+
Labconco 302419000

 BT Barrier Pipet Tips
Pre-Sterile
Neptune® BT1250, BT100, BT10



Filter forceps
blunt end, stainless steel
Millipore XX6200006P



LYSING TUBES
MATRIX D 2 mL/15 mL
MP BIOMEDICALS 116913500/116933050



Microcentrifuge Tubes
1.7 mL/0.6 mL
Axygen Scientific MCT-175-C/MCT-060-L-C



Screw-Cap Centrifuge Tube
5 mL
VWR 10002-738



Falcon® Centrifuge Tubes
Polypropylene, Sterile, 15 mL
Corning® 352096



Falcon® Centrifuge Tubes
Polypropylene, Sterile, 50 mL
Corning® 352070



Fastprep-24 5G™ Sample Preparation
Instrument
MP Biomedicals 116005500



CENTRIFUGE 5430 R
Eppendorf MP2231000510



VWR ANALOG VORTEX MIXER
VWR 10153-838
With tube insert



SHAKING INCUBATOR
71L
Corning® LSE™ 6753



Finnpipette Stepper Pipette
Thermo Scientific™ 4540000



Finntip™ Stepper Pipette Tips
500 uL
Thermo Scientific™ 9404170



96-Well Black Microplates
Polystyrene
Greiner Bio-One 655076



Black Vinyl Films for Fluorescence and
Photoprotection
VWR 89087-692



Varioskan LUX Multimode Microplate Reader
Thermo Fisher VL0L00D0

EQUIPMENT

NAME	CATALOG #	VENDOR
BT Barrier Pipet Tips	BT1250, BT100, BT10	VWR international Ltd
Falcon® Centrifuge Tubes	352096	VWR international Ltd
A2 Biosafety Cabinet	302419000	
CENTRIFUGE 5430 R	MP2231000510	
VWR ANALOG VORTEX MIXER	10153-838	VWR international Ltd
Finnpipette Stepper Pipette	4540000	Fisher Scientific
Microcentrifuge Tubes	MCT-175-C/MCT-060-L-C	VWR international Ltd
Screw-Cap Centrifuge Tube	10002-738	VWR international Ltd
Falcon® Centrifuge Tubes	352070	VWR international Ltd
SHAKING INCUBATOR	6753	
Varioskan LUX Multimode Microplate Reader	VL0L00D0	
Filter forceps	XX6200006P	Emdmillipore
LYSING TUBES	116913500/116933050	VWR international Ltd
Fastprep-24 5G™ Sample Preparation Instrument	116005500	VWR international Ltd
FreeZone® 2.5 L Benchtop Freeze Dryers	700202000	VWR international Ltd
Finnpipette Stepper Pipette Tips	9404170	Fisher Scientific
96-Well Black Microplates	655076	
Black Vinyl Films for Fluorescence and Photoprotection	89087-692	

RNA and DNA extraction

2h

- 1 Freeze dry samples.
Freeze-drying should be as short as possible to reduce sample degradation.
- 2 Expose biosafety cabinet to UV light for 00:15:00 and clean working surface with decontamination solution.
- 3 Prepare falcon tubes and tube rack in biosafety cabinet

Number	Tube type	Usage
2	5 mL falcon tubes	EDTA20% sarcosine
3	15 mL falcon tubes	1 M Tris, 5 mM Tris, 1% STEB

- 4 Prepare Tris buffer 5 mM 8.0

- 4.1 Pour 1 M 8.0 Tris into an RNase free 15 mL Falcon tube

4.2 Add Tris into 500 mL RNase free water in the original package.

5 Prepare STEB ()

Require 500 uL for each sample and 500 uL for preparing 0.05% STEB in assay.

5.1 Pour sarcosine () into an RNase free 5 mL falcon tube.

5.2 Pour EDTA () into an RNase free 5 mL falcon tube.

5.3 Pour Tris buffer (,) into an RNase free 15 mL falcon tube.

5.4 Mix sarcosine (), EDTA () and + Tris buffer (,) to obtain STEB ().

6 Prepare ice bath

7 Remove samples from -80°C freezer and place .

8 Add Tris buffer (,) and STEB () to the bead tube. Place tubes .

9 Transfer sample/blank filter to the tube.
Forceps must be cleaned by ethanol.
Blank filter should be the same type as the one with samples collected.

10 Vortex immediately then put back .






20s

11 Disrupt samples on the bead mill at 6.5 m/s.










30s

12 Keep tubes . Check the label on each tube, restore the label if it fades.

1m 30s

- 13 Disrupt samples on the bead mill at 6.5 m/s. 30s
- 14 Keep tubes  **On ice** . Check the label on each tube, restore the label if it fades. 1m 30s
- 15 Disrupt samples on the bead mill at 6.5 m/s 30s
- 16 Keep tubes  **On ice** . Check the label on each tube, restore the label if it fades. 1m 30s
- 17 Disrupt samples on the bead mill at 6.5 m/s. 30s
- 18 Continuously shake homogenate in a multi-head vortex at the highest reasonable speed for  **01:00:00**
 **Room temperature** 1h
- 19 In the biosafety cabinet, aliquot 150 uL homogenates into RNase free 600 uL tubes.
- 20 Store at  **-80 °C** .

Prepare standard and nucleases

- 21 RNA primary standard ( **200 ug/ml**)
- 21.1 Add  **1 mL** Tris buffer ( **5 mM** ,  **8.0**) to the original package and vortex.
- 21.2 Aliquot 30 uL by stepper with sterile tip to 600 uL RNase free microtubes. Keep frozen at  **-80 °C** .
- 22 DNA primary standard ( **500 ug/ml**)
- 22.1 Add  **2 mL** Tris buffer ( **5 mM** ,  **8.0**) to the original package and vortex.

22.2 Aliquot 10 uL by stepper with sterile tip to 600 uL RNase free microtubes. Keep frozen at -80°C .

23 RNase primary stock (10 mg/ml)

23.1 Add 5 mL Tris buffer (5 mM , pH8.0) to the original package and vortex.

23.2 Aliquot 30 uL by stepper with sterile tip to 600 uL RNase free microtubes. Keep frozen at -20°C .

24 DNase primary stock (5 mg/ml = 10,000 U/mL)

24.1 Add 1 mL Tris buffer (5 mM , pH8.0) to the original package and vortex.

24.2 Aliquot 100 uL by stepper with sterile tip to 600 uL RNase free microtubes. Keep frozen at -20°C .

The assay with nucleases

25 Prepare ice bath.

26 Expose biosafety cabinet to UV light for 00:15:00 and clean working surface with decontamination solution.

27 Prepare falcon tubes, microtubes and tube racks in biosafety cabinet

* Maximum number of samples (including blanks) per assay is 20.

Number	Tube type	Usage
6	5 mL falcon tubes	1 M MgCl ₂ 1 M CaCl ₂ WS-A WS-B WS-C SG-II WS
2	15 mL falcon tubes	5 mM Tris, 0.05% STEB
3	600 uL RNase free tubes	RNase 900 mM MgCl ₂ 900 mM CaCl ₂
2	1.7 mL RNase free tubes	RNA standard DNA standard
33	1.7 mL RNase free tubes	standards reaction
N=Samples + blanks	1.7 mL RNase free tubes	sample dilution
3N	1.7 mL RNase free tubes	diluted sample reaction
5	Microtube racks	Microtubes
1	Tube rack	Falcon tubes

28 Organize and label the tubes as shown below

Set 1: Sample dilution

1	2	3	4	5	6	7	8	9	...	Blank
---	---	---	---	---	---	---	---	---	-----	-------

Set 2: Standards

Tubes A	R1A	R2A	R3A	R4A	R5A	R6A	D1A	D2A	D3A	D4A	D5A
Tubes B	R1B	R2B	R3B	R4B	R5B	R6B	D1B	D2B	D3B	D4B	D5B
Tubes C	R1C	R2C	R3C	R4C	R5C	R6C	D1C	D2C	D3C	D4C	D5C

Set 3: Sample reaction

Tubes A	1A	2A	3A	4A	5A	6A	7A	8A	9A	...	BlankA
Tubes B	1B	2B	3B	4B	5B	6B	7B	8B	9B	...	BlankB
Tubes C	1C	2C	3C	4C	5C	6C	7C	8C	9C	...	BlankC

29 Label tubes for reagents as following.

Add Tris buffer ([M] **5 mM** , **pH 8.0**) to the reagent tubes:

Reagent	5 mM Tris (uL)
0.05% STEB	9X1000+500
RNase	190
900 mM MgCl ₂	20
900 mM CaCl ₂	20
WS-A	2X1000+820
WS-B	2X1000+820
WS-C	2X1000+940
RNA	990+495
DNA	998

30 Add  900 µl Tris buffer ([M]5 mM , pH8.0) to each tube in Set 1

31 Add Tris buffer ([M]5 mM , pH8.0) to the tubes in Set 2. The unit of volume is uL.

Tubes A	650	640	625	600	550	500	640	630	610	580	550
Tubes B	650	640	625	600	550	500	640	630	610	580	550
Tubes C	600	590	575	550	500	450	590	580	560	530	500

32 Add Tris buffer ([M]5 mM , pH8.0) to tubes in Set 3. The unit of volume is uL.

Tubes A	650	650	650	650	650	650	650	650	...	650
Tubes B	650	650	650	650	650	650	650	650	...	650
Tubes C	600	600	600	600	600	600	600	600	...	600

33 Prepare STEB ([M]0.05 %)

33.1 Add  500 µl STEB ([M]1 %) to 0.05% STEB tube, and vortex.

34 Add  250 µl STEB ([M]0.05 %) to Set 2.

35 Place RNase, DNase primary stock, RNA primary standard, DNA primary standard and samples  On ice .

36 Turn on refrigerated centrifuge and set the temperature to  4 °C .

- 37 Turn on shaker/incubator and set temperature to **37 °C** .
- 38 Prepare **900 mM** MgCl_2
- 38.1 Pour **1 M** MgCl_2 solution into 5 mL RNase free Falcon tube
- 38.2 Transfer **180 μl** **1 M** MgCl_2 solution into 900 mM MgCl_2 tube
- 39 Add **60 μl** **900 mM** MgCl_2 to WS-A and WS-B
- 40 Prepare **900 mM** CaCl_2
- 40.1 Pour **1 M** CaCl_2 solution into 5 mL RNase free Falcon tube
- 40.2 Transfer **180 μl** **1 M** CaCl_2 solution into 900 mM CaCl_2 tube
- 41 Add **60 μl** **900 mM** CaCl_2 to WS-A and WS-B
- 42 Prepare RNase working solution **0.5 mg/ml**
Add **10 μl** RNase primary stock (**10 mg/ml**) to RNase tube
- 43 Add **60 μl** **0.5 mg/ml** RNase to WS-B and WS-C, keep WS-B and WS-C **On ice** .
- 44 Add **60 μl** DNase primary stock (**5 mg/ml**) to WS-A, keep WS-A **On ice** .
- 45 Centrifuge extracted samples **10000 x g, 4°C 00:04:00**
- 46 Prepare RNA secondary standard **2 $\mu\text{g/ml}$**


Add  **15 µl** RNA primary standard to RNA standard tube and vortex. Keep  **On ice** .

47 Prepare DNA secondary standard  **1 µg/ml**


Add  **2 µl** DNA primary standard to DNA standard tube and vortex. Keep  **On ice** .

48 Load  **50 µl** WS-A to Tubes A in Set 2 and Set 3.

49 Load  **50 µl** WS-A to Tubes C in Set 2 and Set 3.

50 Load  **50 µl** WS-B to Tubes B in Set 2 and Set 3.




51 Load  **50 µl** WS-C to Tubes C in Set 2 and Set 3.

52 Add  **100 µl** centrifuged samples to its corresponding tubes in Set 1. Vortex each tube.

53 From Set 1, transfer  **250 µl** of diluted samples to each corresponding tubes in Set 3.

54 Add RNA and DNA secondary standards to Set 2 tubes. The unit of volume is µL.




Tubes A	0	10	25	50	100	150	10	20	40	70	100
Tubes B	0	10	25	50	100	150	10	20	40	70	100
Tubes C	0	10	25	50	100	150	10	20	40	70	100


55 Vortex each tube for  **00:00:02** and place all tubes into the shaker/incubator at  **37 °C** , continuously shaking at ^{20m} 200 RPM for  **00:20:00** .


Fluorescence reading

56 Prepare SYBR Green II working solution (SG-II WS)

56.1 Each 96-well microplate requires 1 mL of SG-II WS.

56.2 Wrap 5 mL tube with foil, add  **7 µl** SYBR Green II 10,000X concentrate to each one mL Tris buffer ( **5 mM** ,  **8.0**) in biosafety cabinet.

56.3 Load  10 µl SG-II WS to each well in the microplate with 0.5 mL tip of stepper.
Adhere black film on the top of a microplate lid and cover the plate with this lid.

57 After incubation, stop reaction by vortexing each tube for  00:00:02 and then place into the fridge.

58 Allow samples to reach  Room temperature for  00:02:00 before loading the microplate.

59 Load  190 µl working sample to the microplate.

45m

Standards must be included in each microplate.

	1	2	3	4	5	6	7	8	9	10	11	12
A	R1A	R2A	R3A	R4A	R5A	R6A	1A	1A	1B	1B	1C	1C
B	R1B	R2B	R3B	R4B	R5B	R6B	2A	2A	2B	2B	2C	2C
C	R1C	R2C	R3C	R4C	R5C	R6C	3A	3A	3B	3B	3C	3C
D	D1A	D2A	D3A	D4A	D5A		4A	4A	4B	4B	4C	4C
E	D1B	D2B	D3B	D4B	D5B		5A	5A	5B	5B	5C	5C
F	D1C	D2C	D3C	D4C	D5C		6A	6A	6B	6B	6C	6C
G	9A	9A	9B	9B	9C	9C	7A	7A	7B	7B	7C	7C
H	Blank A	Blank A	Blank B	Blank B	Blank C	Blank C	8A	8A	8B	8B	8C	8C

60  Room temperature Shake black film covered microplate for  00:10:00

10m

61 Setup microplate reader:

Plate: Greiner F bottom chimney well PP 96 well;

Endpoint reading: Ex 490 nm/Em 520 nm;

Shake: Continuous 5s at 600 rpm

62 Read fluorescence and export data to excel sheet.

Standard curve and calculation

63 RNA standard curve

63.1 Concentrations of RNA standards in the microplate

Standard	2 ug/mL (uL)	Tris+WS (uL)	0.05% STEB (uL)	SG II WS (uL)	Final (ng/mL)
R1	0	700	250	50	0
R2	10	690	250	50	20
R3	25	675	250	50	50
R4	50	650	250	50	100
R5	100	600	250	50	200
R6	150	550	250	50	300

- 63.2 Slope of fluorescence in Tube A vs concentration of RNA standard gives $m_{RNA+DNase}$
Slope of fluorescence in Tube B vs concentration of RNA standard gives $m_{RNA+RNase}$

- 63.3 Calculate ρ

$$\rho = \frac{m_{RNA+RNase}}{m_{RNA+DNase}}$$

- 64 Total RNA of the samples

$$\mu g \text{ Total RNA (ml assay)}^{-1}$$

$$= 0.001 \times \frac{(RFU_A - RFU_{A_{Blank}}) - (RFU_C - RFU_{C_{Blank}})}{(1 - \rho)} / m_{RNA+DNase}$$

Where,

RFU_A and RFU_C are the fluorescence in Tube A and Tube C of the same sample.

$RFU_{A_{Blank}}$ and $RFU_{C_{Blank}}$ are the fluorescence in Tube A and Tube C of the blank.

- 65 DNA standard curve

- 65.1 Concentrations of DNA standards in the microplate

Standard	1 ug/mL (uL)	Tris+WS (uL)	0.05% STEB (uL)	SG II WS (uL)	Final (ng/mL)
D1	10	690	250	50	10
D2	20	680	250	50	20
D3	40	660	250	50	40
D4	70	630	250	50	70
D5	100	580	250	50	100

- 65.2 Slope of fluorescence in Tube A vs concentration of DNA standard gives $m_{DNA+DNase}$
Slope of fluorescence in Tube B vs concentration of DNA standard gives $m_{DNA+RNase}$

- 65.3 Calculate δ

$$\delta = \frac{m_{DNA+DNase}}{m_{DNA+RNase}}$$

66 Total DNA of the samples

$\mu\text{g Total DNA (ml assay)}^{-1}$

$$= 0.001 \times \frac{(RFU_B - RFU_{BBlank}) - (RFU_C - RFU_{CBBlank})}{(1 - \delta)} / m_{DNA+RNase}$$

Where,

RFU_B and RFU_C are the fluorescence in Tube B and Tube C of the same sample

RFU_{BBlank} and $RFU_{CBBlank}$ are the fluorescence in Tube B and Tube C of the blank.

67 Dilution factor=40

If,

- Sample is extracted by 1 mL extraction reagent
- In Set 1, sample is diluted to 1/10
- In Set 3, diluted by Tris and all working solutions to 250/950
- In microplate, diluted by SG-II WS to 190/200