

Mar 23,
2020

Viral Sequencing, from Gunk to Graph

David Eccles¹

¹Malaghan Institute of Medical Research (NZ)

In Development

dx.doi.org/10.17504/protocols.io.bd3yi8pw

Coronavirus Method Development Community



David A. Eccles

Malaghan Institute of Medical Research (NZ)



ABSTRACT

This is a fast "gunk to graph" protocol for analysing viral RNA from nasopharyngeal swabs. The approach involves swab lysis and inactivation at the point of sampling, uses a cellulose binding / wash protocol to reduce extraction cost, incorporates sample-specific barcodes during first-strand synthesis, nanopore rapid-attachment primers during PCR amplification, and nanopore sequencing with parallel RAMPART analysis for fast assembly and phylogenetics.

MATERIALS

NAME ▾	CATALOG # ▾	VENDOR ▾
Q5 Hot Start High-Fidelity 2X Master Mix - 500 rxns	M0494L	New England Biolabs
MinION sequencer	View	Oxford Nanopore Technologies
ONT MinION Flow Cell R9.4.1	FLO-MIN106D	Oxford Nanopore Technologies

MATERIALS TEXT

Additional materials TBA.

SAFETY WARNINGS

This protocol is UNTESTED, and is in the early stages of development. Do not trust the protocol; question everything.

Assume samples are potentially infectious during extraction, and make sure to use proper sterile technique to avoid cross-contamination.

Swab Lysis

- 1 Prepare a 1.5 ml centrifuge tube with heated lysis buffer and a cellulose disc

1.1 Add **500 µl** lysis / RNase inactivation buffer ([Twitter reference](#)) to 1.5ml centrifuge tube:

- **10 Milimolar (mM) Tris**
- **10 Milimolar (mM) EDTA**
- **0.5 % volume SDS**
- **150 Milimolar (mM) NaCl**

OR **500 µl** extraction buffer #2 (see [paper](#)):

- **800 Milimolar (mM) guanidine hydrochloride**
- **50 Milimolar (mM) Tris [pH 8]**
- **0.5 % volume Triton X100**
- **1 % volume Tween-20**

1.2 Add a **3 mm diameter** punched disc from [Whatman #1 filter paper](#) (see [paper](#))

1.3 Preheat **1.5 ml** tube to **60 °C**

2 Collect sample using a sterile polystyrene swab with a 30mm breakpoint (e.g. [Puritan 25-3606-U; PurFlock Ultra 6" Sterile Elongated Flock Swab w/Polystyrene Handle, 30mm Breakpoint](#)).

RNA Wash

3 Transfer disc to a new **1.5 ml** tube containing **200 µl** wash buffer using a pipette tip to remove contaminants:

- **10 Milimolar (mM) Tris [pH 8.0]**
- **0.1 % volume Tween-20**

4 Incubate tube at **Room temperature** for **00:01:00**

cDNA Synthesis

5 Add the following additional components into the **200 µl** PCR tube (see the [Nanopore protocol for Sequence-specific cDNA-PCR Sequencing \(SQK-PCS109\)](#)) in a **11 µl** reaction:

- **1 µl x 2 Micromolar (µM) reverse primers**
- **1 µl x 10 Milimolar (mM) dNTPs**
- **9 µl RNase-free water**

Reverse primers should be prefixed with sample-specific barcode sequences (if used) and the ONT reverse anchor sequence, i.e. [5' - ACTTGCTGTCGCTCTATCTTC - [barcode] - [sequence-specific] - 3']

- 6 Mix gently *by flicking the tube* and spin down ⌚ 00:00:05
- 7 Denature RNA and anneal primers at **65 °C** for ⌚ 00:05:00 and then snap cool on a pre-chilled freezer block for ⌚ 00:01:00

- 8 In a separate tube, mix together the following in an **8 µl** reaction:
 - **4 µl** 5X RT Buffer
 - **1 µl** RNaseOUT
 - **1 µl** Nuclease-free water
 - **2 µl** x **10 Micromolar (µM)** Strand-switching primer (SSP)

Note: It might be possible to instead carry out only the first-strand synthesis (i.e. excluding SSP), then use a forward primer, tailed with sample-specific barcode sequences (if used) and the ONT forward anchor sequence, i.e. [5' - TTTCTGTTGGTGCTGATATTGC - [barcode] - [sequence-specific] - 3']. One-Step RT-PCR sequencing kits may help with this (e.g. [OneTaq One-Step RT-PCR Kit](#)). For more details about the reverse anchor sequence and four-primer amplicon sequencing, see the Nanopore protocol for [Four-primer PCR \(SQK-PSK004 or SQK-PBK004\)](#).

- 9 Mix gently *by flicking the tube* and spin down ⌚ 00:00:05
- 10 Add the strand-switching buffer to the snap-cooled, annealed RNA, mix by *flicking the tube* and spin down
- 11 Incubate at **42 °C** for ⌚ 00:02:00
- 12 Add **1 µl** of Maxima H Minus Reverse Transcriptase, to a total volume of **20 µl**
- 13 Mix gently by *flicking the tube* and spin down ⌚ 00:00:05
- 14 Incubate using the following protocol:

Cycle step	Temperature	Time	No. of cycles
Reverse transcription and strand-switching	42° C	90 mins	1
Heat inactivation	85° C	5 mins	1
Hold	4° C	∞	

Thermal cycler settings for reverse transcription and strand switching

PCR amplification

15 In four new **200 µl** PCR tubes, prepare the following reaction at **Room temperature** in a **50 µl** reaction:

- **25 µl** 2X Q5 Hot Start High-Fidelity Master Mix
- **1.5 µl** cDNA primer (cPRM)
- **18.5 µl** Nuclease-free water
- **5 µl** Reverse-transcribed cDNA from the previous step (pool, or single sample)

16 Amplify using the following cycling conditions:

Cycle step	Temperature	Time	No. of cycles
Initial denaturation	95 °C	30 secs	1
Denaturation	95 °C	15 secs	10-18*
Annealing	62 °C	15 secs	10-18*
Extension	65 °C	50 secs per kb	10-18*
Final extension	65 °C	6 mins	1
Hold	4 °C	∞	

Thermal cycler settings for PCR amplification

* The recommended starting point is 14 cycles - adjust this depending on experimental needs.

17 Add **1 µl** of NEB Exonuclease 1 (20 units) directly to each PCR tube to remove unextended primers. Mix by **pipetting**.

18 Incubate the reaction at **37 °C** for **00:15:00**, followed by **80 °C** for **00:15:00**

Bead Cleanup

19 Add 160 µl of resuspended AMPure XP beads to the **1.5 ml** tube and mix by **pipetting**

20 Incubate on a gentle agitator (e.g. hula mixer or rotator mixer) for **00:05:00** at **Room temperature**

21 Spin down **00:00:05** the sample and pellet on a magnet. Keep the tube on the magnet, and pipette off the supernatant.

22 Keep the tube on the magnet and wash the beads with **200 µl** of freshly-prepared **70 % volume ethanol** without disturbing the pellet. Remove the ethanol using a pipette and discard.

23 Repeat the previous step: wash with **200 µl 70 % volume ethanol**, and discard the ethanol / wash liquid.

- 24 Spin down ⌚ 00:00:05 and place the tube back on the magnet. Pipette off any residual ethanol. Allow to dry for ⌚ 00:00:30 [at most] but do not dry the pellet to the point of cracking (the magnetic beads should just start to lose their shiny sheen).
- 25 Remove the tube from the magnetic rack and resuspend pellet in 📏 12 µl of Elution Buffer (EB).
- 26 Incubate at 🌡 Room temperature for ⌚ 00:10:00
- 27 Pellet beads on magnet ⌚ 00:05:00 until the eluate is clear and colourless
- 28 While still on the magnet, carefully remove and retain 📏 12 µl of eluate into a clean 📏 1.5 ml Eppendorf DNA LoBind tube
- 29 Quantify 1 µl of the amplified cDNA library using the Quantus Fluorometer using the ONE dsDNA assay (see [ncov 2019 sequencing protocol, step 16](#))

Adapter Addition


- 30 Add 📏 1 µl of Rapid Adapter (RAP) to the amplified cDNA library
- 31 Mix by *pipetting* and spin down ⌚ 00:00:05
- 32 Incubate the reaction for ⌚ 00:05:00 at 🌡 Room temperature
- 33 Store the prepared library 🌡 On ice until ready to load onto a flow cell.

Nanopore Sequencing

- 34 Load 📏 20 ng sequencing library onto a MinION flow cell (see [ncov 2019 sequencing protocol, step 21](#))
- 35 Start the sequencing run using MinKNOW, using SQK-PCS109 as the sample preparation protocol (see [ncov 2019 sequencing protocol, step 22](#))

RAMPART Analysis

- 36 Analyse the run results using RAMPART (see <https://artic.network/ncov-2019/ncov2019-using-rampart.html>)

 This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited