



Nov 20, 2020

Fungal Extraction from Beetle Galleries

You Li¹, Jiri Hulcr¹¹University of Florida

1

Works for me

dx.doi.org/10.17504/protocols.io.bnutmewn

Bark Beetle Mycobiome Research Coordination Network

ABSTRACT

This protocol describes how to extract fungi from beetle galleries.

This protocol is part of the Bark Beetle Mycobiome (BBM) Research Coordination Network. For more information on the BBM international network: Hulcr J, Barnes I, De Beer ZW, Duong TA, Gazis R, Johnson AJ, Jusino MA, Kasson MT, Li Y, Lynch S, Mayers C, Musvuugwa T, Roets F, Seltmann KC, Six D, Vanderpool D, & Villari C. 2020. Bark beetle mycobiome: collaboratively defined research priorities on a widespread insect-fungus symbiosis. *Symbiosis* 81: 101–113 <https://doi.org/10.1007/s13199-020-00686-9>.

DOI

dx.doi.org/10.17504/protocols.io.bnutmewn

DOCUMENT CITATION

You Li, Jiri Hulcr 2020. Fungal Extraction from Beetle Galleries. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bnutmewn>

LICENSE

This is an open access document distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Oct 23, 2020

LAST MODIFIED

Nov 20, 2020

DOCUMENT INTEGER ID

43635

ABSTRACT

This protocol describes how to extract fungi from beetle galleries.

This protocol is part of the Bark Beetle Mycobiome (BBM) Research Coordination Network. For more information on the BBM international network: Hulcr J, Barnes I, De Beer ZW, Duong TA, Gazis R, Johnson AJ, Jusino MA, Kasson MT, Li Y, Lynch S, Mayers C, Musvuugwa T, Roets F, Seltmann KC, Six D, Vanderpool D, & Villari C. 2020. Bark beetle mycobiome: collaboratively defined research priorities on a widespread insect-fungus symbiosis. *Symbiosis* 81: 101–113 <https://doi.org/10.1007/s13199-020-00686-9>.

Materials

(all autoclaved)

- eppendorf vials
- Pestles or some other grinding tool
- scalpel
- Reserve hood time

- strong forceps
- vials
- ethanol in petri dish for dipping the forceps & scalpel into
- vial rack
- (chisel)
- (pruner)
- pipette
- 1000 uL pipette tips
- EtOH bath

For field sampling of beetle associated fungi:

- ethanol
- squirt bottle (or separate vials with ethanol)
- hard forceps
- scalpel
- scalpel blades
- lighter (for forceps & scalpel sterilization)
- ziploc bags
- permanent markers (alcohol-resistant)

Procedure

Scrape piece of 1 gallery in the hood, scrape into

a. 1 ml of 1X PBS AUTOCLAVED phosphate buffer saline solution in tube

3. mortar & pestle
4. plate directly, refrigeration decreases growing success
5. Plate three different dilutions, refer to the dilution protocol.