



Mar 17,
2020

A protocol for rapid detection of the 2019 novel coronavirus SARS-CoV-2 using CRISPR diagnostics: SARS-CoV-2 DETECTR

Mammoth Biosciences¹, James P. Broughton², Wayne Deng^{3,4}, Clare L. Fasching², Jasmeet Singh², Charles Y. Chiu^{3,4,5}, Janice S. Chen²

¹[Mammoth Biosciences, Inc., South San Francisco, California, USA], ²Mammoth Biosciences, Inc., South San Francisco, California, USA, ³Department of Laboratory Medicine, University of California, San Francisco, California, USA, ⁴UCSF-Abbott Viral Diagnostics and Discovery Center, San Francisco, California, USA, ⁵Department of Medicine, Division of Infectious Diseases, University of California, San Francisco, California, USA

1 Works for me dx.doi.org/10.17504/protocols.io.bcmtiu6n

Coronavirus Method Development Community

Mammoth Biosciences


ABSTRACT

DISCLAIMER: This protocol has not been approved by the FDA and should not be used as a clinical diagnostic

Introduction

Given the global health emergency, rapid transmission, and severe respiratory disease associated with the outbreak of the 2019 novel coronavirus (SARS-CoV-2), Mammoth Biosciences has reconfigured our DETECTR platform to rapidly and accurately detect SARS-CoV-2 using a visual lateral flow strip format within 30 minutes from sample to result. To ensure specificity of detection, we selected a high-fidelity CRISPR detection enzyme and designed sets of gRNAs that can either 1) differentiate SARS-CoV-2 or 2) provide multi-coronavirus strain detection. SARS-CoV-2 DETECTR couples CRISPR detection with isothermal pre-amplification using primers based on protocols validated by the US Centers for Disease Control and Prevention (CDC) and World Health Organization (WHO). Currently in the United States, the CDC SARS-CoV-2 real-time RT-PCR diagnostic panel has a laboratory turnaround time of approximately 4-6 hours, with results that can be delayed for >24 hours after sample collection due to shipping requirements. In addition, these tests are only available in CDC-designated public health laboratories certified to perform high-complexity testing.

Mammoth is working to enable point of care testing (POCT) solutions that can be deployed in areas at greatest risk of transmitting SARS-CoV-2 infection, including airports, emergency departments, and local community hospitals, particularly in low-resource countries. Leveraging an "off-the-shelf" strategy to enable practical solutions within a short time frame, we describe here a protocol that is fast (<30 min), practical (available immediately from international suppliers), and validated using contrived samples.



Specifications	
Targets	<ul style="list-style-type: none"> ● N-gene (SARS-CoV-2 specific) ● E-gene (SARS-CoV, bat-SARS-like-CoV, and SARS-CoV-2 coronaviruses) ● RNase P (human sample control)

Limit of detection	70-300 copies/μl input
--------------------	------------------------

Table 1: SARS-CoV-2 DETECTR assay workflow and specifications.




Figure 1: SARS-CoV-2 DETECTR has a limit of detection (n=7) of 156-625 copies per 20 μl reaction (or 70-300 copies per μl input) and generates a clear visible signal on lateral flow strips within 30 minutes sample to result.

Acknowledgements: We thank Vikram Joshi, Nefeli Tsaloglou and Xin Miao for advice and helpful discussions in the preparation of this whitepaper.


Conflicts of Interest: JPB, CLF, JS and JSC are employees of Mammoth Biosciences, CYC is on the Scientific Advisory Board of Mammoth Biosciences, and JSC is a co-founder of Mammoth Biosciences. JPB, CLF, JS, CYC and JSC are co-inventors on CRISPR-related technologies.

EXTERNAL LINK
<https://mammoth.bio/2020/02/15/white-paper-a-protocol-for-rapid-detection-of-sars-cov-2-using-crispr-sars-cov-2-detectr/>

ATTACHMENTS
[SARS-CoV-2.pdf](#)

GUIDELINES
Appendix

While we were preparing this whitepaper, another [protocol for SARS-CoV-2 detection using CRISPR diagnostics \(SHERLOCK, v.20200214\)](#) was published. We compare the assay workflows and specifications between CRISPR diagnostics and established CDC/WHO protocols below. (Note: as of this publication, CRISPR diagnostics workflows have not yet been approved by the FDA)



Appendix Figure 1: Comparison of SARS-CoV-2 assay workflows for DETECTR, SHERLOCK, and CDC/WHO

	SARS-CoV-2 DETECTR	SARS-CoV-2 SHERLOCK	CDC SARS-CoV2 qRT-PCR
Target	N gene & E gene (N gene gRNA compatible with CDC N2 amplicon, E gene compatible with WHO protocol)	S gene & Orf1ab gene	N-gene (3 amplicons)
Sample control	RNase P	None	RNase P
Limit of Detection	70-300 copies/μl input	10-100 copies/μl input	1 copy/μL input
Assay reaction time	~30 min	~60 min	~45-60 minutes
Assay components	RT-LAMP (62 °C, 20 min), Cas12 (37 °C, 10 min), Lateral flow (RT, 2 min)	RT-RPA (42 °C, 25 min), IVT + Cas13 (37 °C, 30 min), Lateral flow (RT, 2 min)	UDG digestion (25 °C, 2 min), reverse transcription (50 °C, 15 min), denature (95 °C, 2 min), amplification (95 °C, 3 min; 55 °C 30 sec; 45 cycles)
Heavy instrumentation required	No	No	Yes
FDA EUA approval	No	No	Yes

Appendix Table 1: Comparison of SARS-CoV-2 specifications for CRISPR diagnostic protocols to the current CDC assay.

MATERIALS TEXT

SARS-CoV-2 DETECTR Reagents

Step 1: Isothermal amplification (62°C, 20 min)

- RT-LAMP Master Mix (Supplier: NEB)


WarmStart LAMP Kit (DNA and RNA) - 100 rxns
by New England Biolabs
Catalog #: E1700S

- [DNA oligos \(Supplier: IDT\)](#)

- Primer sequences:

Name	Sequence (5' → 3')
N-gene F3	AACACAAGCTTTTCGGCAG
N-gene B3	GAAATTTGGATCTTTGTTCATCC
N-gene FIP	TGCGGCCAATGTTTGTAAATCAGCCAAGGAAATTTGGGGAC
N-gene BIP	CGCATTGGCATGGAAGTCACTTTGATGGCACCTGTGTAG
N-gene LF	TTCCTTGTCTGATTAGTTC
N-gene LB	ACCTTCGGGAACGTGGTT
E-gene F3	CCGACGACGACTACTAGC
E-gene B3	AGAGTAAACGTAAAAAGAAGGTT
E-gene FIP	ACCTGTCTCTTCCGAAACGAATTTGTAAGCACAAAGCTGATG
E-gene BIP	CTAGCCATCCTTACTGCGCTACTCACGTTAACAATATTGCA
E-gene LF	TCGATTGTGTGCGTACTGC
E-gene LB	TGAGTACATAAGTTCGTAC
RNaseP POP7 F3*	TTGATGAGCTGGAGCCA
RNaseP POP7 B3*	CACCCTCAATGCAGAGTC
RNaseP POP7 FIP*	GTGTGACCCTGAAGACTCGGTTTTAGCCACTGACTCGGATC
RNaseP POP7 BIP*	CCTCCGTGATATGGCTCTTCGTTTTTTCTTACATGGCTCTGGTC
RNaseP POP7 LF*	ATGTGGATGGCTGAGTTGTT
RNaseP POP7 LB*	CATGCTGAGTACTGGACCTC


* RNaseP POP7 primers published in Curtis et al. , (2018) .



Curtis KA, Morrison D, Rudolph DL, Shankar A, Bloomfield LSP, Switzer WM, Owen SM (2018). A multiplexed RT-LAMP assay for detection of group M HIV-1 in plasma or whole blood.. Journal of virological methods.
<https://doi.org/10.1016/j.jviromet.2018.02.012>

Step 2: Cas12 detection (37°C, 10 min)

- LbCas12a (Supplier: NEB)



EnGen Lba Cas12a (Cpf1) - 70 pmol
 by New England Biolabs
 Catalog #: M0653S

- [crRNA \(Supplier: Synthego\)](#)
-
- [Reporter \(Supplier: IDT\)](#)

Name	Sequence (5' → 3')
N gene gRNA (SARS-CoV-2 specific)	UAAUUUCUACUAAGUGUAGAUAUCCCCAGCGCUUCAGCGUUC
E gene gRNA (pan-coronavirus)	UAAUUUCUACUAAGUGUAGAUGUGGUAUUCUUGCUAGUUAC
RNase P gRNA (Sample control)	UAAUUUCUACUAAGUGUAGAUGACCGCGAGCGGGUUCUGA

Step 3: Lateral flow (RT, 2 min)

[Milenia HybriDetect 1 lateral flow strips \(Supplier: TwistDx\)](#)

Sample Equipment

- Pipette tips
- 37 °C heat block
- 62°C heat block
- Microcentrifuge
- Eppendorf tubes
- Pipettes
- Lateral flow strips
- Sample collection device (nasopharyngeal swab)
- Timer

SAFETY WARNINGS

Please see SDS (Safety Data Sheet) for hazards and safety warnings.

Prepare nucleic acid sample and CRISPR reagents

1 

Extract patient RNA following [CDC recommendations](#).

2 Prepare *LbCas12a* RNP complexes for the samples to be tested. **One** complex for N-gene, E-gene, and RNase P gRNAs is needed **for each sample**.

Reagent	Volume	Final Concentration
Nuclease-free water	15.75 µl	
10X NEBuffer 2.1	2 µl	1X
1 µM LbCas12a	1 µl	50 nM
1 µM gRNA	1.25 µl	62.5 nM
TOTAL VOLUME	20 µl	

3 

Incubate LbCas12a with gRNA to generate RNP complexes for  **00:30:00** at  **37 °C**.

4 

Add *reporter substrate* to final concentration of  **500 Nanomolar (nM)**.

- 5 Place reactions  **On ice** until ready to proceed.



Complexes are stable at 4°C for at least 24 hours.

Run DETECTR reaction

- 6 

 **On ice**, prepare **three** *RT-LAMP reactions*, **one each** for N-gene, E-gene, and

RNase P primer sets:

Reagent	Volume	Final Concentration
10X Isothermal Amplification Buffer (NEB)	2.5 µl	
100 mM MgSO ₄ (NEB)	1.13 µl	6.5 mM (4.5 mM added, 2 mM in 1X IsoAmp Buffer)
10 mM dNTPs (NEB)	3.5 µl	1.4 mM
10X Primer Mix	2.5 µl	0.2 µM F3 / 0.2 µM B3 / 1.6 µM FIP / 1.6 µM BIP / 0.8 µM LF / 0.8 µM LB
Bst 2.0 polymerase (NEB)	1 µl	8 units / rxn
Warmstart RTx (NEB)	0.5 µl	7.5 units / rxn
Nuclease-free water	3.87 µl	
Nucleic acid sample	5 µl	
TOTAL VOLUME	25 µl	



- 7 

Incubate at  **62 °C** for  **00:20:00**.



Note: Use precaution when opening amplification tubes to prevent amplicon contamination.

- 8 

Combine  **2 µl** of the *RT-LAMP reaction* with  **20 µl** of the *LbCas12a RNP complex* with the appropriate gRNA.

- 9 

Add  **80 µl** **1X NEBuffer 2.1**.

- 10 

Incubate at  **37 °C** for  **00:10:00**.

11 Insert *Milenia HybriDetect 1 (TwistDx) lateral flow strip* directly into reaction.

12 

Allow the lateral flow strip to run for **00:02:00** at **Room temperature** and observe the result.

Test interpretation

13 



Note: The *line closest to the sample pad* is the *control line* and the line that appears farthest from the sample pad is the *test line* (see Figure 1). A sample with complete cleavage of the reporter molecule may appear to have no signal at the control line.




Figure 1 | SARS-CoV-2 DETECTR has a limit of detection (n=7) of 156-625 copies per 20 µl reaction (or 70-300 copies per µl input) and generates a clear visible signal on lateral flow strips within 30 minutes sample to result.

N-gene	E-gene	RNase P	Interpretation
+	+	+/-	SARS-CoV-2 positive
+	-	+/-	Indeterminate
-	+	+/-	Indeterminate
-	-	+	SARS-CoV-2 negative
-	-	-	QC failure



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited