



DEC 13, 2022

WORKS FOR ME

1

Western blotting to detect ATP13A2 and ATP13A3

COMMENTS 0

DOI

dx.doi.org/10.17504/protocols.io.81wgbyzqovpk/v1[Marine Houdou](#)¹, [Peter Vangheluwe](#)¹¹KU Leuven

Marine Houdou

ABSTRACT

Protocol to detect ATP13A2 and ATP13A3 via Western Blotting.

DOI

dx.doi.org/10.17504/protocols.io.81wgbyzqovpk/v1

PROTOCOL CITATION

Marine Houdou, Peter Vangheluwe 2022. Western blotting to detect ATP13A2 and ATP13A3.

protocols.io<https://dx.doi.org/10.17504/protocols.io.81wgbyzqovpk/v1>

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Dec 02, 2022

LAST MODIFIED

Dec 13, 2022

PROTOCOL INTEGER ID

73481

Citation: Marine Houdou, Peter Vangheluwe Western blotting to detect ATP13A2 and ATP13A3

<https://dx.doi.org/10.17504/protocols.io.81wgbyzqovpk/v1>

This is an open access protocol distributed under the terms of the **Creative Commons Attribution License** (<https://creativecommons.org/licenses/by/4.0/>), which permits unrestricted use, distribution, and reproduction in any medium,

■ **Antibodies:**

- Goat anti-mouse IgG (H+L) secondary antibody HRP conjugated: Thermo Scientific, 31430
- Goat anti-rabbit IgG (H+L) secondary antibody HRP conjugated: Thermo Scientific, 31460
- Mouse monoclonal anti-GAPDH (lot #067M4785V, dilution 1:5,000): Sigma, G8795.
- Rabbit anti-ATP13A2 antibodies (lot #0000102992, dilution 1:1,000): Sigma, A3361.
- Rabbit anti-ATP13A3 antibody (lot # 000035781, dilution 1:2,000): Atlas Antibodies, HPA029471.

- 0.25% Trypsin-EDTA: Gibco, 25200056

- Dulbecco's Phosphate Buffered saline modified without calcium chloride and magnesium chloride (DPBS): Gibco, D8537

- Micro-BCA Protein Assay Kit: Pierce BCA Protein Assay Kit, Thermo Scientific, 23225

- NuPAGE LDS sample buffer: Invitrogen, NP0007

- Ponceau staining: Sigma, P7170

- Pre-cast 4-12% Bis-Tris gels: Invitrogen, NP0321BOX

- PVDF membranes: Thermo Scientific, 88518

- RIPA Lysis and Extraction Buffer: Invitrogen: 89900

- SIGMAFAST Protease Inhibitor Cocktail Tablets, EDTA-Free: Sigma: S8830

- SuperSignal West Pico PLUS chemiluminescent Substrate: Thermo Scientific, 34095

Harvesting cells

- 1 Depending on cell type, collect the cells by scrapping them with a scrapper in Dulbecco's Phosphate Buffered saline modified without calcium chloride and magnesium chloride (DPBS) (SH-SY5Y) or, using 0.25% Trypsin-EDTA (HMEC-1) for which stop enzymatic reaction by adding culture medium.

2 Centrifuge cell suspensions at 450 g (SH-SY5Y) or 2500 rpm (HMEC-1), 4°C for 00:05:00 .

5m

3 Resuspend cell pellets with DPBS and centrifuge following the same indications as in 2. Repeat once.

4 Discard supernatants and keep cell pellets on ice.

Cell lysis and protein concentration determination

5 Resuspend cell pellets in RIPA buffer (RIPA Lysis and Extraction Buffer supplemented with protease cocktail inhibitors).

6 Vortex 00:00:30 and keep on ice for 00:30:00 .

30m 30s

7 Centrifuge at 20,000g, 4°C for 00:30:00 .

30m

8 Keep supernatants on ice to proceed with protein concentration determination using the micro-BCA Protein Assay Kit.

7h 10m

SDS-PAGE



9 Loading

9.1 Mix 20 µg of protein with NuPAGE LDS sample buffer and 5% β-mercaptoethanol final.

9.2 For this specific protocol to detect ATP13A2 and ATP13A3, do not boil samples.


9.3 Load protein on pre-cast 4-12% Bis-Tris gels. Include at least one lane with a protein ladder.

10 Running

10.1 Run for  00:10:00 at 100V and  01:30:00 at 110-130V.

1 h 40m

11 Transfer

11.1 Transfer onto PVDF membranes using a liquid transfer and following settings: 100V,  01:15:00, 4°C.

1 h 15m

12 Ponceau staining



12.1 Rinse membrane with distilled water.

12.2 Incubate membrane with Ponceau staining for  00:05:00, 19 rpm.

5m


12.3 Scan membrane if necessary.

13 Blocking


13.1  01:00:00 Block membranes with blocking buffer (5% milk powder in 1X TBS and 0.1% Tween20 (REF)) for  01:00:00 at room temperature, 19 rpm.

2h

14 Primary antibodies


14.1 Incubate membrane with primary antibodies in solution (1% bovine serum albumin in 1X TBS-Tween20 (TBS-T) buffer),  Overnight at 4°C, 19 rpm.

1h


14.2 Wash membrane three times for  00:05:00 in TBS-T, 19 rpm.

5m

15 Secondary antibodies

15.1 Incubate membrane with peroxidase-conjugated secondary antibodies in solution (1% milk powder in 1X TBS-T) for  01:00:00 at room temperature and, 19 rpm.

1h

15.2 Wash membrane five times for  00:05:00 in TBS-T, 19 rpm.

5m

16 Detection

- 17 Use a chemiluminescence reagent to detect signal and acquire with a Biorad Camera (Vilber Lourmat) and its software (ImageLab).