

♦ Thawing, Expanding, and Freezing Human Fibroblasts

Jacob Marsh¹, Rj Martinez², Celeste Karch³

¹Washington University in Saint Louis - WUSTL (MO), ²Washington University, Saint Louis, ³Washington University in St Louis

May 26, 2020

1 Works for me dx.doi.org/10.17504/protocols.io.bgunjwve

 $Neurode generation \ Method \ Development \ Community \ Tech. \ supporternal: \ ndcn-help@chanzuckerberg.com$



Celeste Karch Washington University in St Louis

Thawing Human Fibroblasts

- 1 Prior to thawing cells, prepare culture media
 - 1.1 hFibroblast Media:

DMEM/F12, 10-20%

FBS

1% L-Glut (or Glutamax)

1% Antibiotic

2 Remove cells from liquid nitrogen storage and thaw quickly by placing vial in a $837\,^{\circ}$ C water bath for approximately

© 00:01:00



Be sure to monitor cells during thaw and make sure they are removed from water bath just prior to completely thawing

- Quickly transfer thawed cells into conical tube containing $\Box 9$ ml of room temperature DMEM/F12
- 4 **3750 rpm, Room temperature 00:03:00**
- 5 Aspirate supernatant, be careful not to distrub cell pellet
- 6 Re-suspend cells in **5 ml** of hFibroblast Media by pipetting cells up and down 2-3 times to achieve a single cell suspension
- 7 Transfer cell suspension to T25 cell culture flask

protocols.io
1
05/26/2020

- 8 Incubate at § 37 °C and 5% CO2 overnight
- 9 Replace hFibroblast Media every 2-3 days

Expanding Human Fibroblasts

- 10 Upon reaching 80-90% confluency, fibroblasts will create a swirling pattern. This is the appropriate time to split the cells
 - Confluent fibroblasts can be split from one T25 cell culture flask into one T75 cell culture flask OR two T25 cell culture flasks.
- 11 Aspirate hFibroblast Media from T25 cell culture flask
- 12 Wash cells with **5 ml** 1X PBS
- Add 1 ml of 0.25% Trypsin to cells. Incubate at 37 °C for 00:05:00 to allow cells to dissociate from cell culture flask
 - After incubation period, check cells under a microscope to determine the extent to which cells have detached from the flask. If many cells remain attached, incubate at § 37 °C for an additional © 00:05:00.
- Re-suspend cells in **5 ml** of hFibroblast Media. Pipet cells up and wash flask 4-6 times to collect all cells and to achieve a single cell suspension.
- 15 Add cell suspension to **77 ml** of fresh hFibroblast Media and transfer to T75 cell culture flask
- 16 Incubate at § 37 °C at 5% C02 overnight
- 17 Change media every 2-3 days until fibroblasts become confluent

At this point fibroblasts can be passaged again, used for experiments, or frozen down and banked

Freezing Human Fibroblasts

- 18 Upon reaching 80-90% confluency, fibroblasts will create a swirling pattern. This is the appropriate time to split the cells
 - It is important to keep fibroblast cultures at a low passage number, thus when thawing cells it should be a priority to expand cells for freezing over expanding cells for experiments.
 - One T75 cell culture flask will produce enough cells to freeze 2-3 **1 ml** aliquots.
- 19 Aspirate hFibroblast Media from T25 cell culture flask
- 20 Wash cells with 5 ml 1X PBS
- 21 Add 2 ml of 0.25% Trypsin to cells. Incubate at 37 °C for 00:05:00 to allow cells to dissociate from cell culture flask
- Re-suspend cells in **7 ml** of hFibroblast Media. Pipet cells up and wash flask 4-6 times to collect all cells and to achieve a single cell suspension.
- 23 Transfer cell suspension to conical tube for centrifugation
- 24 **3500 rpm, Room temperature 00:05:00**
- 25 Aspirate supernatant, be careful not to distrub cell pellet
- 26 Re-suspend cell pellet in **1.5 ml** hFibroblast Media.
- 27 Add 1.5 ml 2X Freezing Media (FBS, 10% DMSO) to cell suspenion from Step 26.



 This should result in a final volume of 3 ml which can be equally distributed into 11 ml aliquots across three cryovials.

28 Transfer 1 ml to each pre-labeled cryovial

29 Immediately place cells in & -80 °C freezer for © 72:00:00

30 Transfer cells to long-term liquid nitrogen storage