



May 20, 2021

# Gel electrophoresis

Victoria Jackson<sup>1</sup><sup>1</sup>Living Systems Institute, University of Exeter

1 Works for me

Share

This protocol is published without a DOI.

Victoria Jackson  
Living Systems Institute, University of Exeter

## ABSTRACT

Gel electrophoresis

## PROTOCOL CITATION

Victoria Jackson 2021. Gel electrophoresis. **protocols.io**  
<https://protocols.io/view/gel-electrophoresis-buy6nxze>

## LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

May 13, 2021

## LAST MODIFIED

May 20, 2021

## PROTOCOL INTEGER ID

49918

## PARENT PROTOCOLS

In steps of

[Standard PCR protocol](#)

### Making the gel

8m 30s

1 For **200 mL** of a 1.2% agarose gel, add **2.4 g** of agarose to **200 mL** TAE or TBE in a conical flask and mix by swirling the flask.

2 Microwave on a high power for about **00:03:00**, or until boiling, swirling the flask approximately every **00:00:30** during heating to ensure all the agarose dissolves.

3m 30s




Use heat proof gloves.

3 Allow the agarose to cool until you are able to comfortably touch the outside of the flask (approximately **00:05:00**<sup>5m</sup>).

4 Add **10 µl** SafeView stain to the hot agarose and swirl to mix.

5 Pour the agarose into the casting tray. Pop any bubbles or move them to the edges of the mould using a pipette tip. Place the comb as desired.

6 Allow the gel to set for approximately  **00:40:00** .



40m

#### Loading the gel

7 Place the gel (still in the casting tray) gently into the tank, ensuring there is enough TAE or TBE running buffer to cover it.

8 If the buffer used for PCR already contains loading dye, skip to step 11.

9 On a piece of parafilm, dot  **1 µl** of 6x loading dye using a pipette.

10 Pipette  **5 µl** ladder or  **6 µl** PCR product directly onto the loading dye and mix by pipetting up and down several times.

11 Load  **6 µl** into each well.

#### Running the gel

1h

12 Place the lid on the tank, making sure that the wells are at the negative electrode end of the tank. The samples should run from black to red (negative to positive).

13 Set the power supply to between **80-120V**, depending on the expected amplicon size and desired duration of electrophoresis.

14 Run the electrophoresis until the visible dye bands are about 75% of the way down the gel, approximately **1-1.5 hours**.



Once complete, ensure you disconnect the power supply from the gel tank, before removing the lid.