



VERSION 1

MAR 18, 2024

OPEN ACCESS



DOI:

[dx.doi.org/10.17504/protocols.io.x54v92e8ml3e/v1](https://dx.doi.org/10.17504/protocols.io.x54v92e8ml3e/v1)

**Protocol Citation:** Wenxin Zhao, Sheng Zhong 2024. Protocol for Preparing Brain Samples for MUSIC. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.x54v92e8ml3e/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol and it's working

## Protocol for Preparing Brain Samples for MUSIC V.1

Wenxin Zhao<sup>1</sup>, Sheng Zhong<sup>1</sup>

<sup>1</sup>Shu Chien-Gene Lay Department of Bioengineering, University of California San Diego

Human BioMolecular Atlas Program (HuBMAP) Method Development Community

MUSIC



Wenxin Zhao

Shu Chien-Gene Lay Department of Bioengineering, University ...

### DISCLAIMER

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The [protocols.io](#) team notes that research involving animals and humans must be conducted according to internationally-accepted standards and should always have prior approval from an Institutional Ethics Committee or Board.

### ABSTRACT

Here states the detailed procedure to prepare brain samples for MUSIC study.

Created: Mar 15, 2024

Last Modified: Mar 19, 2024

PROTOCOL integer ID: 96794

## Tissue pulverization and crosslinking

- 1 Cut a portion of post-mortem human brain frontal cortex sample on dry ice with heavy razor blades, and collect 50 mg of the sample in a 1.5 mL LoBind tube.
- 2 Thaw the 50 mg of brain sample on ice, and chop the tissue into smaller pieces by pestle. Store the rest of the sample at -80°C.
- 3 Incubate the sample with 10 mL of 2 mM disuccinimidyl glutarate (DSG) in 1X PBS in a 15 mL LoBind tube at room temperature for 45 min with gentle rotation.
- 4 Wash once with 10 mL of 1X PBS by centrifugation at 1,000 x g for 4 min.
- 5 Resuspend the sample in 15 mL of 1X PBS containing 3% formaldehyde, and incubate for 10 min with a gentle rotation.
- 6 Quench the crosslinking reaction by the addition of 5 mL of 1.25 M glycine followed by an incubation of 5 min with a rotation.

- 7 Centrifuge the sample at 1,000 x g for 4 min, and wash the sample twice with ice-cold 1X PBS containing 0.3% BSA (wt/vol).

## Nuclei isolation

- 8 Use Chromium Nuclei Isolation kit (10X genomics, 1000494) to isolate nuclei from crosslinked cortex samples.

- 9 Transfer 50 mg frozen tissue into pre-chilled sample dissociation tube.

- 10 Add 400 µL of Lysis Buffer to Sample Dissociation Tube. Dissociate tissue with plastic pestle until homogeneous.

A	B
Component	Volume (µL)
Lysis Reagent	1000
Reducing Agent B	1
Surfactant A	10
Total Volume	1011

### Lysis Buffer

- 11 Add 600 µL of lysis buffer into the tube, and mix 10 times by pipetting. Incubate on ice for 10 min.
- 12 Equally load the solution into two nuclei isolation column, and centrifuge the tubes at 16,000 x g for 20 sec at 4°C.

- 13 Vortex the flowthrough in the collection tube that contains nuclei for 10 sec at 3,200 rpm or max speed to resuspend nuclei.
- 14 Centrifuge the collection tubes for 3 min at 500 x g at 4°C to pellet nuclei. Carefully discard the supernatant.
- 15 Resuspend the nuclei in 500 µL of Debris Removal Buffer provided by the kit by pipetting 15 times.

A	B
Component	Volume (µL)
Debris Removal Reagent	500
Reducing Agent B	0.5
Total Volume	500.5

#### Debris Removal Buffer

- 16 Centrifuge the nuclei at 700 x g for 10 min at 4°C. Carefully discard the supernatant.
- 17 Resuspend the nuclei in 1 mL of Wash and Resuspension Buffer.

A	B
Component	Volume (µL)
1X PBS	1750
10% BSA	200
RNase Inhibitor (40X)	50
Total Volume	2000

#### Wash and Resuspension Buffer

- 18 Centrifuge the nuclei at 500 x g for 5 min at 4°C. Carefully discard the supernatant.

- 19** Resuspend the nuclei again in 1 mL of Wash and Resuspension Buffer.
- 20** Centrifuge the nuclei at 500 x g for 5 min at 4°C. Carefully discard the supernatant as much as possible.
- 21** The nuclei are subjected to nuclei counting and the following procedures.