



Version 2

Dec 11, 2020

Leaf Protein Extraction for Immunoblot (Soybean, Cowpea, Tobacco) V.2

Steven J Burgess¹¹University of Illinois at Urbana-Champaign

1

Works for me

dx.doi.org/10.17504/protocols.io.bqnmvmc6

Ag SynBio Lab UIUC



Steven Burgess

University of Illinois at Urbana-Champaign

ABSTRACT

This protein extraction protocol was developed for analysis of protein abundance in leaf tissue by immunoblot. It was optimized for dicot species including *Glycine max*.

Quantities - 450 μ L protein extraction buffer (PEB) per three 13.4 mm diameter leaf disks (size #7 Humboldt Cork Borer)

Note from Agrisera: Keeping sample volumes in a range of 0.2-0.5 mL has been found to contribute to better extraction results, an upscale in volume is not recommended, if no cork borer available the suggestion is 500 μ L PEB per 100 mg of plant tissue.

- Final volume ~400 μ L

- Expected yield: 1.5 - 6 μ g / μ L

- Total yield: 6 - 24 mg protein

- Amount of total soluble protein to load per lane 30 μ g

Note: when using the TissueLyser II it is recommended to use 2 mL centrifuge tubes in conjunction with 4mm beads. 1.5mL tubes are narrow at the end and samples will not grind properly if you use the 4mm beads.

DOI

dx.doi.org/10.17504/protocols.io.bqnmvmc6

PROTOCOL CITATION

Steven J Burgess 2020. Leaf Protein Extraction for Immunoblot (Soybean, Cowpea, Tobacco).

protocols.io

<https://dx.doi.org/10.17504/protocols.io.bqnmvmc6>

Version created by Steven Burgess

KEYWORDS

Glycine max, soybean, protein extraction

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

IMAGE ATTRIBUTION

Image reproduced from the QIAGEN website <https://www.qiagen.com/us/products/human-id-and-forensics/automation/tissuelyser-ii/#orderinginformation>

CREATED

Dec 11, 2020

LAST MODIFIED

Dec 11, 2020

PROTOCOL INTEGER ID

45485

MATERIALS TEXT

- TissueLyser II ([QIAGEN; 85300](#))
- 2 mL centrifuge tubes
- 4mm SPEX™ stainless steel grinding beads ([SPEX; 2150](#))
- Humboldt brass cork borer set (07-865-10B; Fisher Scientific)
- 13.4 mm diameter, flash-frozen leaf disks
- 4x Protein Extraction Buffer (PEB) (8 % SDS (w/v); 40 % glycerol (v/v); 0.25 M Tris HCl (pH 6.8))

(from Abcam) SDS grade is important for high-quality protein separation: a protein stained background along individual gel tracks with indistinct or slightly distinct protein bands are indicative of old or poor quality SDS. This buffer is essentially the same as Laemmli buffer so proteins can be directly loaded onto PAGE gels.

SAFETY WARNINGS

Perform steps with protein extraction buffer in a fumehood

ABSTRACT

This protein extraction protocol was developed for analysis of protein abundance in leaf tissue by immunoblot. It was optimized for dicot species including Glycine max.

Quantities - 450 µL protein extraction buffer (PEB) per three 13.4 mm diameter leaf disks (size #7 Humboldt Cork Borer)





Note from Agrisera: Keeping sample volumes in a range of 0.2-0.5 mL has been found to contribute to better extraction results, an upscale in volume is not recommended, if no cork borer available the suggestion is 500 µL PEB per 100 mg of plant tissue.

- Final volume ~400 µL
- Expected yield: 1.5 - 6 µg / µL
- Total yield: 6 - 24 mg protein
- Amount of total soluble protein to load per lane 30 µg

Note: when using the TissueLyser II it is recommended to use 2 mL centrifuge tubes in conjunction with 4mm beads. 1.5mL tubes are narrow at the end and samples will not grind properly if you use the 4mm beads.

BEFORE STARTING

Prepare a working solution of 1x Protein Extraction Buffer

- 1x PEB (10 mL) ( **2.5 mL** 4x PEB;  **250 µl** 2-mercaptoethanol;  **100 µl** protease inhibitor cocktail;  **7.3 mL** dH₂O)

Buffer Preparation

20m

20m

1 Incubate 4x PEB (if stored at 4 °C) at  **50 °C** for  **00:20:00** to re-suspended precipitated SDS in buffer.

2 Make  **500 µl** 1x PEB per sample by diluting 4x stock with dH₂O.

It is advised to make more 1x PEB than necessary to avoid running out.

Tissue Lysis

1h 3m

3

Pre-cool components for tissue lysis

3.1 Pre-cool 4mm SPEX stainless steel grinding beads at -80°C for 00:30:00

30m

3.2 Pre-cool the TissueLyser inserts for LN_2 for 00:30:00

30m

4

Prepare sample tubes for tissue lysis

4.1 Using forceps cooled in LN_2 , add one SPEX bead per sample tube.

4.2 Insert tubes into pre-cooled TissueLyser II cassettes, ensuring a balanced number of samples between cassettes.

5

Grind tissue for 00:01:30 at 20 Hz

1m 30s

Warning: Do not exceed this frequency, higher frequencies increase the number of cases where the steel beads will break the lid of centrifuge tubes resulting in sample loss. If using the SPEX 2150 beads it is necessary to use 2 mL centrifuge tubes to ensure proper grinding, in 1.5 mL tubes the bead will not reach the tapered bottom leaving samples unground

TissueLyser II
Bead Mill

QIAGEN

85300




6

Remove cassettes from TissueLyser II and submerge them in LN_2 to prevent thawing

1m 30s


- 7 Grind tissue for  00:01:30 at 20 Hz

This repeat is to ensure all tissue is correctly grounded. There will be odd instances where leaf tissue has not properly ground for reasons such as samples sticking to the edge of tubes, check for this on removal, it may be necessary to repeat again.

- 8 Remove tubes with forceps cooled in LN₂. Samples can be stored at  -80 °C for several months before processing further.

It is advised not to proceed directly from samples in LN₂ to the next step as they take longer to thaw and resuspend in extraction buffer. A period of 30 mins on dry ice can be beneficial.

Dissolve ground powder 1h 3m

- 9 Add  450 µl 1x PEB to ground powder (3x 1cm leaf disk).

- 10 

1m


Vortex immediately (maximum speed) to mix. Approximately  00:01:00

No lumps should be visible at this stage if the sample is ground well.

- 11 Heat samples at  95 °C for  00:05:00 to denature proteins and inactivate proteases.

5m

Prolonged heating can cause cleavage of peptide bonds leading to artifacts. Do not delay heating after sample buffer addition as not all proteases are denatured by SDS and partially denatured peptides are sensitive to protease degradation.

- 12 Vortex heated samples to shear nucleic acids ( 00:00:05 ; max speed)

5s

- 13 Spin samples for  00:03:00 at  10000 x g, Room temperature to pellet insoluble material

3m

the pellet should be white/light-grey an intense green color of the pellet can indicate that disruption was not

optimal and extraction conditions need to be adjusted (e.g. improved grinding, or adjusting buffer volume)

- 14 Transfer supernatant to a new 1.5 mL centrifuge tube and transfer to ice or store at $-20\text{ }^{\circ}\text{C}$ for up to a month

Be careful not carry over debris