



DEC 28, 2022

## Processing Stack-of-Stars DCE Data

Rong Zhou<sup>1,2</sup>

<sup>1</sup>Department of Radiology University of Pennsylvania;

<sup>2</sup>mAbramson Cancer Center University of Pennsylvania



Miguel Romanello Joaquim

### ABSTRACT

Step-wise protocol for reconstruction of Stack-of-stars acquired DCE series, T1 and B1 maps and PK (pharmacokinetic) modeling of DCE data using a reference region model is provided.

OPEN ACCESS

#### DOI:

[dx.doi.org/10.17504/protocols.io.5qpvobb2bl4o/v1](https://dx.doi.org/10.17504/protocols.io.5qpvobb2bl4o/v1)

**Protocol Citation:** Rong Zhou 2022. Processing Stack-of-Stars DCE Data.

**protocols.io**

<https://dx.doi.org/10.17504/protocols.io.5qpvobb2bl4o/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working  
We use this protocol and it's working

**Created:** Feb 08, 2022

**Last Modified:** Dec 28, 2022

**PROTOCOL integer ID:**  
57940

## Reconstruction of Stack-of Stars (SoS) radial k-space sampled DCE data

### 1 Image reconstruction

#### 1.1 VFA and AFI images

The same reconstruction procedure is used for images acquired by the variable flip angle (VFA) and actual flip angle (AFI) stack-of-stars (SoS) pulse sequences.

1. Apply a Fourier transform in the slice direction to separate slices
2. Shift image in the slice direction to match geometry of reference T2-weighted image
3. Apply the following corrections to each view:
  - 3.1. Center each view in k-space by moving its peak to the center of k-space
  - 3.2. Phase normalize
  - 3.2. Correct for off-resonance frequency using the average phase difference between views in opposite directions
4. Re-grid radial k-space data to a 128x128 Cartesian grid as described by O'Sullivan et al. To summarize, for each slice:
  - 4.1. Multiply signal of each point by its respective area on a Voronoi diagram of the points (including zerofill points) in k-space
  - 4.2. Re-grid each radially defined point to its nearest Cartesian coordinates using its Kaiser-Bessel index
5. Apply Fourier transform to now Cartesian-defined k-space

#### CITATION

O'Sullivan, JD (1985). A Fast Sinc Function Gridding Algorithm for Fourier Inversion in Computer Tomography. IEEE Transactions on Medical Imaging.

LINK

[10.1109/TMI.1985.4307723](https://doi.org/10.1109/TMI.1985.4307723)

## 1.2 DCE images

The k-space weighted image contrast (KWIC) method described by Song et al (2000 and 2004). was used to reconstruct the DCE images. To summarize, for DCE-MRI, the KWIC method uses radial acquisition's inherent oversampling of the k-space center by only using a subset of the acquired views in that region. By using a sliding window to select the views that are included in that central region, multiple images can be created from a single acquisition of k-space, thus increasing the temporal resolution.

After applying the KWIC, the resulting k-space data is reconstructed using the same method as **1.1- VFA and AFI Images**.

#### CITATION

Song HK, Dougherty L (2000). k-Space weighted image contrast (KWIC) for contrast manipulation in projection reconstruction MRI. Magnetic Resonance in Medicine.

LINK

[https://doi.org/10.1002/1522-2594\(200012\)44:6<825::AID-MRM2>3.0.CO;2-D](https://doi.org/10.1002/1522-2594(200012)44:6<825::AID-MRM2>3.0.CO;2-D)

## CITATION

Song HK, Dougherty L (2004). Dynamic MRI with projection reconstruction and KWIC processing for simultaneous high spatial and temporal resolution. Magnetic Resonance in Medicine.

LINK

<https://doi.org/10.1002/mrm.20237>

## 2 B1 field map generation from AFI images

1. Compute pixel-wise actual flip angles using the equation

$$\alpha = \arccos \left( \frac{n-1}{n-r} \right), \text{ where}$$

$$n = TR_2/TR_1, \text{ and}$$

$$r = S_1/S_2, \text{ where}$$

$S_i$  is the signal intensity from the image acquired at  $TR_i$ .

2. Divide resulting actual flip angle maps by nominal flip angle to yield normalized pixel-wise B1 field maps.
3. Fit normalized B1 field maps to 3D 3rd degree polynomial using a least-squares fit to yield final B1 field maps used for T1 correction (step 3.1)

## 3 T1 map generation from VFA images

1. Compute T1 values using a non-linear least squares fit of the VFA image signal intensity to the Ernst equation

$$S(\alpha_{\text{act}}) = M_0 \left( \frac{\sin \alpha_{\text{act}}}{(1 - E_1)} \{1 - E_1 \cos \alpha_{\text{act}}\} \right), \text{ where}$$

$$E_1 = e^{-TR/T_1}, \text{ and}$$

$$\alpha_{\text{act}} = B1 * \alpha_{\text{nom}}.$$

## 4 Generate tissue masks for reference region (muscle), kidney and tumor.

1. Open DCE images in an imaging processing software (eg. ImageJ)

Note: Images immediately following contrast agent injection tend to outline structures well. Anticipating slight movement of the slices either due to respiratory motion or due to other reasons, a T2-weighted scan is usually acquired before the DCE sequence and used for ROI masking.

2. Create a mask accessible by your analysis software which defines skeletal muscle and any other ROIs of interest.

For example, on ImageJ:

1. Create a new empty image with same dimensions as DCE image in File -> New -> Image
2. Draw ROIs by hand and add to ROI Manager by pressing shortcut "t"
3. Set ROIs in empty image to a values for each tissue by Process -> Map -> Set
4. Save as raw image by File -> Save As -> Raw Data

## 5 DCE metric map generation from DCE images

### 5.1 Compute contrast agent concentration time-course from signal time-course for each voxel and for spinal muscle (reference region)

For each voxel:

1. Compute actual (B1-corrected) flip angle at voxel using  $\alpha = B1 * \alpha_{\text{nom}}$ .
2. Normalize signal time course by mean signal prior to bolus injection
3. Compute T1 time-course using equation

$$T_1(t) = -TR \cdot \ln \left( \frac{A-S_R(t)}{A-S_R(t) \cos \alpha} \right), \text{ where}$$

$$A = \frac{1 - e^{-TR/T_{1,0}} \cos \alpha}{1 - e^{-TR/T_{1,0}}} \text{ and}$$

$T_{1,0}$  is the baseline T1, obtained from the previously generated T1 map (step 3)

4. Estimate contrast agent concentration time-course using equation

$$C(t) = \frac{1}{R_{Gd}} \left( \frac{1}{T_1(t)} - \frac{1}{T_{1,0}} \right), \text{ where}$$

$$R_{Gd} = 4.6 \text{ s}^{-1} \cdot \text{mM}^{-1} \text{ (check for your specific contrast agent)}$$

### 5.2 Compute spinal muscle (reference region) concentration time-course

1. Using a manually defined tissue mask for the spinal muscle, compute the mean concentration time-course for the entire muscle ROI

### 5.3 Compute quantitative DCE parametric maps

For each voxel:

1. Using the reference region model (Jones et al.) and the muscle as a reference tissue, fit for  $K^{trans}$  and  $v_e$  using the following equation:

$$\frac{dC_t}{dt} = \frac{K^{trans}_t}{K^{trans}_m} \frac{dC_m}{dt} + \frac{K^{trans}_t}{v_{e,t}} C_m - \frac{K^{trans}_t}{v_{e,t}} C_t, \text{ where}$$

the reference  $K^{trans}$  of muscle,  $K^{trans}_m = 0.1/\text{min}$ ,  
the reference  $v_e$  of muscle,  $v_{e,m} = 0.1$ , and

C is the concentration of contrast agent, and the subscripts m and t refer to muscle (the reference region) and the tissue in the voxel being analyzed.

Note:  $K^{trans}_m$  and  $v_{e,m}$  are reference values, therefore voxel  $K^{trans}$  and  $v_e$  values are relative to those assigned when fitting for the equation above. The values of  $0.1/\text{min}$  and 0.1 are from Cardenas-Rodriguez et al.

#### CITATION

Jones KM, Pagel MD, Cárdenas-Rodríguez J (2018). Linearization improves the repeatability of quantitative dynamic contrast-enhanced MRI. Magnetic resonance imaging.

LINK

<https://doi.org/10.1016/j.mri.2017.11.002>

#### CITATION

Cardenas-Rodriguez J, Howison CM, Pagel MD (2013). A linear algorithm of the reference region model for DCE-MRI is robust and relaxes requirements for temporal resolution. Magnetic Resonance Imaging.

LINK

<https://doi.org/10.1016/j.mri.2012.10.008>

## 6 Obtaining ROI metrics from maps

While pixel-wise parametric maps allow us to assess heterogeneity of a specific metric within the tumor, we also compute all parameters ( $K^{trans}, v_e$ ) from the ROI of interest (tumor, kidney, and phantom):

1. Extract  $K^{trans}$  and  $v_e$  values for each voxel in the ROI.
2. Compute ROI metrics of choice (eg. mean, median, percentile values, standard deviation)

Note: To mitigate the impact of pixels whose  $K^{trans}$  and  $v_e$  values are outliers, we prefer median instead of mean of all pixels in the ROI.