

Jul 29, 2024 Version 1

Opentrons Pipeline: DNA Extraction with the Mag-Bind Blood & Tissue DNA HDQ 96 Kit Tissue Protocol V.1

DOI

dx.doi.org/10.17504/protocols.io.bp2l62rqrgqe/v1

Thalía Silvestre¹

¹UNMSM



Thalía Silvestre

UNMSM

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.bp2l62rqrgqe/v1

Protocol Citation: Thalía Silvestre 2024. Opentrons Pipeline: DNA Extraction with the Mag-Bind Blood & Tissue DNA HDQ 96 Kit Tissue Protocol. [protocols.io <https://dx.doi.org/10.17504/protocols.io.bp2l62rqrgqe/v1>](https://dx.doi.org/10.17504/protocols.io.bp2l62rqrgqe/v1)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development

We are still developing and optimizing this protocol

Created: July 28, 2024

Last Modified: July 29, 2024

Protocol Integer ID: 104196

Keywords: Opentrons, OT-2, Mag-Bind Blood & Tissue DNA HDQ 96 Kit Tissue Protocol



Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment.

Abstract

This protocol is an automated pipeline to extract a full 96-well plate of DNA from tissue lysates.

This protocol was developed and optimized for the following:

Platform: Opentrons OT-2 automated pipetting robot

Kit: Mag-Bind® Blood & Tissue DNA HDQ Prefilled 96 Kit

Recommended number of samples: 96

Two ellutions: Firts ellution(60μL) , Second ellution (30μL)



Samples Lysis

- 1 Prepare 96 samples lysis following the Mag-Bind® Blood & Tissue DNA HDQ 96 Kit Tissue Protocol:
 - 1.1 Mince up to 10 mg tissue and transfer to a 1.5 mL tubes (one for each sample). (Cut each 10 mg sample with sterilized cutters as small as possible)
 - 1.2 Add 250 µL TL Buffer to each sample.
 - 1.3 Add 20 µL Proteinase K Solution to each sample. Vortex to mix.
 - 1.4 Incubate at 55°C overnight in a shaking water bath.
 - 1.5 Centrifuge at maximum speed ($\geq 4,000 \times g$) for 5 minutes to pellet undigested tissue debris and transfer 200 µL of the supernatant to the VWR 96-Well Deep Well Plate (1000 µL) without disturbing the undigested pellet

List of materials to start to use OT-2

- 2 It is necessary to have this equipment to use this pipeline.
Opentrons Equipment List:


A	B
OT-2	
OT-2 8-Chann el Pipette P300	
OT-2 Magnetic Module GEN2	Slot 7

- 3 Tips & Labware:

A	B
Four Opentrons 200µL Filter Tips	Slot: 5,6,8,9
Two Nest 1-Well Reservoirs (195 mL)	Slot: 10 and 11
VWR 96-Well Deep Well Plate (1000 µL)	Slot 7 on the Magnetic module
Two Nest Well plate (100 µL)	Slot 2 and 3
Two Nest 12 Well Reservoirs (15 mL)	Slot 1 and 4

Star OT-2 run

- 4 The pipeline start on the step number 7 of the Mag-Bind [®] Blood & Tissue DNA HDQ 96 Kit Tissue Protocol

Load the pipeline on the Opentrons app:  Extracción.py 11KB

The VWR 96-Well Deep Well Plate (1000 µL) labware can be download from:
dx.doi.org/10.17504/protocols.io.dm6gp39njvzp/v1

4.1 Order of the reagents on the Two Nest 12 Well Reservoirs (15 ml)

Slot 1:

A	B	C	D	E	F	G	H	I	J	K	L
Well 1	Well 2	Well 3	Well 4	Well 5	Well 6	Well 7	Well 8	Well 9	Well 10	Well 11	Well 12
11.5 mL of AL buffer	11.5 mL of AL buffer	Mag-Bind mixed with HBQ buffer (20:340) in a total of 12 mL	Mag-Bind mixed with HBQ buffer (20:340) in a total of 12 mL	Mag-Bind mixed with HBQ buffer (20:340) in a total of 12 mL	15 ml of VH B buffer	15 ml of VH B buffer	15 ml of VH B buffer	15 ml of VH B buffer			

Slot 4:

	A	B	C	D	E	F	G	H	I	J	K	L
	Well 1	Well 2	Well 3	Well 4	Well 5	Well 6	Well 7	Well 8	Well 9	Well 10	Well 11	Well 12
	15 ml of VH B buffer	15 ml of VH B buffer	15 ml of VH B buffer	15 ml of VH B buffer	15 ml of SP M buffer	15 ml of SP M buffer	15 ml of SP M buffer	15 ml of SP M buffer				Ellution Buffer