



Oct 01, 2024

## snASE - MAS-Seq protocol

DOI

**[dx.doi.org/10.17504/protocols.io.ewov19w7plr2/v1](https://dx.doi.org/10.17504/protocols.io.ewov19w7plr2/v1)**

Asa Shin<sup>1</sup>, Aziz Al'Khafaji<sup>2</sup>

<sup>1</sup>None; <sup>2</sup>Broad Institute



Asa Shin

None

OPEN  ACCESS



DOI: **[dx.doi.org/10.17504/protocols.io.ewov19w7plr2/v1](https://dx.doi.org/10.17504/protocols.io.ewov19w7plr2/v1)**

**Protocol Citation:** Asa Shin, Aziz Al'Khafaji 2024. snASE - MAS-Seq protocol. **protocols.io**

**<https://dx.doi.org/10.17504/protocols.io.ewov19w7plr2/v1>**

**License:** This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** September 19, 2024

**Last Modified:** October 01, 2024

**Protocol Integer ID:** 108045

## Abstract

Modified version of PacBio's MAS-seq protocol for snASE samples.



## Quality Control

- 1 1 Bring the Qubit 1X dsDNA HS working solution and standards to room temperature.
- 2 2 Pulse vortex or pipette mix each sample to homogenize the DNA in solution.
- 3 3 Quick spin each sample to collect liquid.
- 4 4 Take a 1  $\mu$ L aliquot from each sample.
- 5 5 Measure DNA concentration with a Qubit fluorometer using the 1X dsDNA HS kit.
- 6 6 Dilute each sample to 1.0-1.5 ng/ $\mu$ L in elution buffer or water, based on the Qubit reading.
- 7 7 Measure DNA size distribution with a Bioanalyzer system using the High Sensitivity DNA Kit.
- 8 8 Proceed to the next step of the protocol if sample quality is acceptable.

## TSO PCR

- 9 Prepare PCR mix per sample:

A	B
Component	Volume
Nuclease-free water	Make up volume
MAS PCR mix	25 $\mu$ L
MAS 5' capture primer mix	10 $\mu$ L
10x 3' cDNA library	Up to 15 $\mu$ L
Total Volume	50 $\mu$ L



- 10 Pipette-mix RM1.
- 11 Quick spin RM1 in a microcentrifuge to collect liquid.
- 12 Select the TSO PCR program based on cDNA input.
- 13 Cleanup with 1.5X SMRTbell cleanup beads
- 14 Add 1.5X v/v (75  $\mu$ L) of resuspended, room-temperature SMRTbell cleanup beads to each tube of amplified cDNA.
- 15 Pipette mix the beads until evenly distributed.
- 16 Quick spin the tube strip in a microcentrifuge to collect liquid.
- 17 Leave at room temperature for 10 minutes to allow DNA to bind beads.
- 18 Place tube strip in a magnetic separation rack until beads separate fully from the solution.
- 19 Slowly pipette off the cleared supernatant without disturbing the beads. Discard the supernatant.
- 20 Slowly dispense 200  $\mu$ L, or enough to cover the beads, of freshly prepared 80% ethanol into each tube. After 30 seconds, pipette off the 80% ethanol and discard.
- 21 Repeat the previous step.
- 22 Remove residual 80% ethanol: • Remove tube strip from the magnetic separation rack. • Quick spin tube strip in a microcentrifuge. • Place tube strip back in a magnetic separation rack until beads separate fully from the solution. • Pipette off residual 80% ethanol and discard.



- 23 Remove tube strip from the magnetic rack. Immediately add 42  $\mu$ L of elution buffer to each tube and resuspend the beads by pipetting 10 times or until evenly distributed.
- 24 Quick spin the tube strip in a microcentrifuge to collect liquid.
- 25 Leave at room temperature for 5 minutes to elute DNA.
- 26 Place tube strip in a magnetic separation rack until beads separate fully from the solution.
- 27 Slowly pipette off the cleared supernatant without disturbing the beads. Transfer supernatant to a new tube strip. Discard old tube strip with beads.
- 28 Recommended: Evaluate sample concentration. Take a 1  $\mu$ L aliquot from each tube. Measure DNA concentration with a Qubit fluorometer using the 1x dsDNA HS kit.
- 29 Proceed to the next step of the protocol if sample quantity is acceptable (at least 150 ng) and not exceeding 1  $\mu$ g. If the sample quantity is higher than 1  $\mu$ g, only carry forward to Step 3 with a maximum input of 1  $\mu$ g

## TSO artifact removal

- 30 Bring MAS capture beads kit to room temperature. Resuspend the beads by vortexing.
- 31 Transfer 10  $\mu$ L resuspended MAS capture beads per sample to a PCR tube. Scale up the amount of beads if processing more than 4 samples (with 10% overage). If preparing more than 40  $\mu$ L of beads, use a 1.5 mL LoBind tube instead of PCR tube.
- 32 Place the tube on the magnet until the beads separate fully from the solution.
- 33 Carefully remove and discard the supernatant while the tube remains on the magnet. Avoid touching the bead pellet with the pipette tip.
- 34 Remove the tube from the magnet.  
Add 40  $\mu$ L MAS bead binding buffer along the insidewall of the tube where the beads are collected and gently resuspend by pipetting using wide bore tips.  
DO NOT VORTEX.

Note: the solution may be viscous. Highly recommend using wide bore tips to avoid foaming. When excess bubbles are present, lower cDNA recovery is expected.

Quick-spin the tube in a microcentrifuge if needed.

Note: Scale up the volume of MAS capture binding buffer accordingly if preparing more than 40  $\mu\text{L}$  of beads.

- 35 Place the tube on the magnet until the beads separate fully from the solution and remove the supernatant.

- 36 Resuspend the beads in 40  $\mu\text{L}$  MAS bead binding buffer by pipetting slowly using wide bore tips.

DO NOT VORTEX.

Note: the solution may be viscous. Highly recommend using wide bore tips to avoid foaming. When excess bubbles are present, lower cDNA recovery is expected.

Note: Scale up the volume of MAS capture binding buffer accordingly, if preparing more than 40  $\mu\text{L}$  of beads.

Distribute 40  $\mu\text{L}$  of resuspended MAS capture beads into the appropriate number of PCR tubes before proceeding to Step 3.8.

- 37 Add 40  $\mu\text{L}$  of a solution containing the biotinylated DNA-fragments (from Step 2.18) to the resuspended beads. Mix carefully using wide bore tips to avoid foaming of the solution.

- 38 Incubate the tube at room temperature for 15 minutes on a rotator to keep the beads in suspension. Quick-spin the tube in a microcentrifuge to collect liquid.

- 39 Place the tube on the magnet until the beads separate fully from the solution and remove the supernatant.

- 40 Resuspend the MAS capture beads/DNA-complex in 80  $\mu\text{L}$  MAS bead washing buffer by pipette mixing until evenly distributed.

- 41 Place the tube on the magnet until the beads separate fully from the solution and remove the supernatant.

- 42 Remove the tube from the magnet. Resuspend the MAS capture beads/DNA-complex in 80  $\mu\text{L}$  MAS bead washing buffer by pipette mixing until evenly distributed.

- 43 Place the tube on the magnet until the beads separate fully from the solution and remove the supernatant.

- 44 Remove the tube from the magnet. Resuspend the MAS capture beads/DNA complex in 80  $\mu\text{L}$  nuclease free water by pipette mixing until evenly distributed.



- 45 Place the tube on the magnet until the beads separate fully from the solution and remove the supernatant.
- 46 Resuspend the capture beads/DNA-complex in 40  $\mu\text{L}$  of elution buffer by pipette mixing until evenly distributed.
- 47 Add 2  $\mu\text{L}$  MAS enzyme to the sample with capture beads to cleave the captured DNA products from MAS capture beads.
- 48 Pipette-mix each sample and a very quick spin in a microcentrifuge to collect liquid.

- 49 Run the TSOartifact removal program. Heated lid set at 47°C

A	B	C
Step	Time	Temperature
1	30 minutes	37°C
2	Hold	4°C

- 50 Place the tube on the magnet for 1 minute and move the supernatant containing the library to a fresh tube.
- 51 Cleanup with 1.5X SMRTbell cleanup beads.
- 52 Add 1.5X v/v (63  $\mu\text{L}$ ) of resuspended, room-temperature SMRTbell cleanup beads to each sample.
- 53 Pipette-mix the beads until evenly distributed.
- 54 Quick-spin the tube strip in a microcentrifuge to collect liquid.
- 55 Leave at room temperature for 10 minutes to allow DNA to bind beads.
- 56 Place tube strip in a magnetic separation rack until the beads separate fully from the solution.



- 57 Slowly pipette off the cleared supernatant without disturbing the beads. Discard the supernatant.
- 58 Slowly dispense 200  $\mu$ L, or enough to cover the beads, of freshly prepared 80% ethanol into each tube. After 30 seconds, pipette off the 80% ethanol and discard.
- 59 Repeat the previous step.
- 60 Remove residual 80% ethanol: • Remove the tube strip from the magnetic separation rack. • Quick-spin the tube strip in a microcentrifuge. • Place the tube strip back in a magnetic separation rack until the beads separate fully from the solution. Pipette off residual 80% ethanol and discard.
- 61 Remove the tube strip from the magnetic rack. Immediately add 46  $\mu$ L of elution buffer to each tube and resuspend the beads by pipetting 10 times or until evenly distributed.
- 62 Quick-spin the tube strip in a microcentrifuge to collect liquid.
- 63 Leave at room temperature for 5 minutes to elute DNA.
- 64 Place the tube strip in a magnetic separation rack until beads separate fully from the solution.
- 65 Slowly pipette off the cleared supernatant without disturbing the beads. Transfer supernatant to a new 0.5 mL LoBind tube. Discard the old tube strip with beads.
- 66 Recommended: Evaluate sample concentration. • Take a 1  $\mu$ L aliquot from each tube. Measure DNA concentration with a Qubit fluorometer using the 1x dsDNA HS kit.
- 67 Proceed to the next step of the protocol if sample quantity is acceptable (maximum 50 ng). If cDNA amount is >50 ng, dilute the cDNA to 50 ng using elution buffer in a total volume of 45 $\mu$ L

## MAS PCR

- 68 Set up the following PCR reaction mix per sample on ice.  
Reaction Mix 2 (RM2):



A	B
<b>Master mix components</b>	<b>Volume for 16X concatenation*</b>
PCR Grade Water	176–X $\mu$ L
MAS PCR mix	220 $\mu$ L
55 ng of amplified cDNA from Step 67	X $\mu$ L
Total volume	396 $\mu$ L

- 69 Quick-spin RM2 in a microcentrifuge to collect liquid.
- 70 Add 22.5  $\mu$ L of RM2 to a new PCR tube on ice. Repeat this step to prepare a total of 16 tubes per sample (each containing 22.5  $\mu$ L of RM2)
- 71 Add 2.5  $\mu$ L of MAS primers premix A-PQ into each of 16 PCR tubes on ice, with each tube getting only one primer. Tube 1= Primer mix A, Tube 2=Primer mix B,... Tube 16= Primer mix PQ. DO NOT ADD MORE THAN ONE PRIMER PREMIX TO A TUBE
- 72 Pipette-mix each sample. The total volume of each tube should be 25.0  $\mu$ L.
- 73 Quick-spin the strip tubes in a microcentrifuge to collect liquid.
- 74 Run the MAS PCR program. Reactions can be held overnight in the cycler.  
Note: if the total sample quantity is less than 50 ng, follow the table below for cycle number.

A	B
cDNA Input amount	Cycle Number
30–50 ng	9
12.5–29.9 ng	10

- 75 Clean up with 1.5X SMRTbell cleanup beads.
- 76 Pool entire volume of all 16 reactions into a single 1.5 mL LoBind tube.





- 77 Add 1.5X v/v (600  $\mu$ L) of resuspended, room-temperature SMRTbell cleanup beads to the PCR pool.
- 78 Pipette-mix the beads until evenly distributed.
- 79 Quick-spin the tube strip in a microcentrifuge to collect all liquid from the sides of the tubes.
- 80 Leave at room temperature for 10 minutes to allow DNA to bind beads.
- 81 Place 1.5mL LoBind tube in a magnetic separation rack until beads separate fully from the solution.
- 82 Slowly pipette off the cleared supernatant without disturbing the beads. Discard the supernatant.
- 83 Slowly dispense 1 mL, or enough to cover the beads, of freshly prepared 80% ethanol into each tube. After 30 seconds, pipette off the 80% ethanol and discard.
- 84 Repeat the previous step.
- 85 Remove residual 80% ethanol: • Remove the LoBind tube from the magnetic separation rack. • Quick-spin the LoBind tube in a microcentrifuge. • Place the LoBind tube back in a magnetic separation rack until beads separate fully from the solution. Pipette off residual 80% ethanol and discard.
- 86 Remove the LoBind tube from the magnetic rack. Immediately add 50  $\mu$ L of elution buffer to each tube and resuspend the beads.
- 87 Quick-spin the LoBind tube in a microcentrifuge.
- 88 Incubate at room temperature for 5 minutes to elute DNA.



- 89 Place the LoBind tube in a magnetic separation rack until beads separate fully from the solution.
- 90 Slowly pipette off the cleared supernatant without disturbing the beads. Transfer supernatant to a new LoBind tube. Discard old tube with beads.
- 91 Recommended: Evaluate sample concentration. • Take a 1  $\mu\text{L}$  aliquot from each tube, dilute with 9  $\mu\text{L}$  of elution buffer. Using 1  $\mu\text{L}$  of the dilution, measure DNA concentration with a Qubit fluorometer using the 1x dsDNA HS kit.
- 92 Proceed to the next step of the protocol if sample quantity is acceptable (required input: 10  $\mu\text{g}$ ). Do not proceed if less than 6  $\mu\text{g}$  is available.

## MAS Array Formation

- 93 In a 0.2 mL PCR tube, add 10  $\mu\text{g}$  of sample from Step 4.22, in 47  $\mu\text{L}$  of volume. Dilute with elution buffer going into this step if sample is too concentrated.
- 94 Add 10  $\mu\text{L}$  of MAS enzyme to create single-stranded extensions on PCR-amplified cDNA fragments to enable subsequent directional assembly of 16 PCR products.
- 95 Pipette-mix each sample.
- 96 Run the MAS primer digestion program.
- 97 Add 3  $\mu\text{L}$  of MAS adapter bc01–04 (use a single barcode per sample) and 20  $\mu\text{L}$  of MAS ligation additive to each sample for a total volume of 80  $\mu\text{L}$ .
- 98 Pipette-mix each sample.
- 99 Add the following components in the order and volume listed below to a new microcentrifuge tube.  
Adjust component volumes for the number of samples being prepared, plus 10% overage. For individual preps, add components directly to each sample in the order and volume listed below.  
Reaction Mix 3 (RM3):

A	B
Component	Volume



A	B
MAS Ligase Buffer	10 µL
MAS Ligase	10 µL
<b>TOTAL volume</b>	<b>20 µL</b>

100 Pipette-mix RM3 with wide bore tips.

101 Quick-spin RM3 in a microcentrifuge to collect liquid.

102 Add 20 µL of RM3 to each sample.

103 Pipette-mix each sample with wide bore tips

104 Run the MAS array ligation program. Heated lid set at 47°C

A	B	C
Step	Time	Temperature
1	30 Minutes	37°C
2	Hold	4°C

105 Cleanup with 1.2X SMRTbell cleanup beads

106 Add 1.2X v/v (120 µL) of resuspended, room-temperature SMRTbell cleanup beads to each sample.

107 Pipette-mix the beads with wide bore tips until evenly distributed.

108 Quick-spin the tube strip in a microcentrifuge to collect liquid.

109 Leave at room temperature for 10 minutes to allow DNA to bind beads.



- 110 Place the tube strip in a magnetic separation rack until beads separate fully from the solution.
- 111 Slowly pipette off the cleared supernatant without disturbing the beads. Discard the supernatant.
- 112 Slowly dispense 200  $\mu$ L, or enough to cover the beads, of freshly prepared 80% ethanol into each tube. After 30 seconds, pipette off the 80% ethanol and discard.
- 113 Repeat the previous step.
- 114 Remove residual 80% ethanol: • Remove the tube strip from the magnetic separation rack.
- 115 Add 1.2X v/v (120  $\mu$ L) of resuspended, room-temperature SMRTbell cleanup beads to each sample.
- 116 Pipette-mix the beads with wide bore tips until evenly distributed.
- 117 Quick-spin the tube strip in a microcentrifuge to collect liquid.
- 118 Leave at room temperature for 10 minutes to allow DNA to bind beads.
- 119 Place the tube strip in a magnetic separation rack until beads separate fully from the solution.
- 120 Slowly pipette off the cleared supernatant without disturbing the beads. Discard the supernatant.
- 121 Slowly dispense 200  $\mu$ L, or enough to cover the beads, of freshly prepared 80% ethanol into each tube. After 30 seconds, pipette off the 80% ethanol and discard.
- 122 Repeat the previous step.



- 123 Remove residual 80% ethanol:  
Remove the tube strip from the magnetic separation rack.  
Quick-spin the tube strip in a microcentrifuge.  
Place the tube strip back in a magnetic separation rack until beads separate fully from the solution.  
Pipette off residual 80% ethanol and discard.
- 124 Remove the tube strip from the magnetic rack. Using a wide bore pipette tip, immediately add 43  $\mu\text{L}$  of elution buffer to each tube and resuspend the beads by pipetting 10 times or until evenly distributed.
- 125 Quick-spin the tube strip in a microcentrifuge to collect liquid.
- 126 Leave at room temperature for 5 minutes to elute DNA.
- 127 Place the tube strip in a magnetic separation rack until beads separate fully from the solution.
- 128 Slowly pipette off the cleared supernatant without disturbing the beads. Transfer supernatant to a new 0.5 mL LoBind tube. Discard old tube strip with beads.
- 129 Recommended: Evaluate sample concentration.  
Take a 1  $\mu\text{L}$  aliquot from each tube, dilute with 4  $\mu\text{L}$  of elution buffer.  
Measure DNA concentration with a Qubit fluorometer using the 1x dsDNA HS kit.  
The required amount of purified MAS array products to proceed with the DNA damage repair step is 5 $\mu\text{g}$

## DNA Damage repair

- 130 In a 0.2 mL PCR tube, add 5  $\mu\text{g}$  of sample from Step 5.26, in 42  $\mu\text{L}$  of volume. Dilute with elution buffer  
going into this step if sample is too concentrated
- 131 Add the following components in the order and volume listed below to a new microcentrifuge tube.  
Adjust component volumes for the number of samples being prepared, plus 10% overage. For individual preps, add components directly to each sample in the order and volume listed below.  
Reaction Mix 4 (RM4):

A	B
Component	Volume
Repair buffer	6 $\mu\text{L}$



A	B
DNA Damage repair mix	2 $\mu$ L

132 Pipette-mix RM4.

133 Quick-spin RM4 in a microcentrifuge to collect liquid.

134 Add 8  $\mu$ L of RM4 to each sample. Total volume should equal 50  $\mu$ L.

135 Pipette-mix each sample with wide bore tips.

136 Quick-spin the strip tube in a microcentrifuge to collect liquid.

137 Run the DNA damage repair program. Heated Lid set to 47°C

A	B	C
Step	Time	Temperature
1	30 Minutes	37°C
2	Hold	4°C

138 Clean up with 1.2X SMRTbell cleanup beads

139 Add 1.2X v/v (60  $\mu$ L) of resuspended, room-temperature SMRTbell cleanup beads to each sample.

140 Pipette-mix the beads with wide bore tips until evenly distributed.

141 Quick-spin the tube strip in a microcentrifuge to collect liquid.

- 142 Leave at room temperature for 10 minutes to allow DNA to bind beads
- 143 Place the tube strip in a magnetic separation rack until beads separate fully from the solution.
- 144 Slowly pipette off the cleared supernatant without disturbing the beads. Discard the supernatant.
- 145 Slowly dispense 200  $\mu$ L, or enough to cover the beads, of freshly prepared 80% ethanol into each tube. After 30 seconds, pipette off the 80% ethanol and discard.
- 146 Repeat the previous step.
- 147 Remove residual 80% ethanol: • Remove the tube strip from the magnetic separation rack. • Quick-spin the tube strip in a microcentrifuge. • Place the tube strip back in a magnetic separation rack until beads separate fully from the solution. • Pipette off residual 80% ethanol and discard.
- 148 Remove the tube strip from the magnetic rack. Using a wide bore pipette tip, immediately add 40  $\mu$ L of elution buffer to each tube and resuspend the beads by pipetting 10 times or until evenly distributed.
- 149 Quick-spin the tube strip in a microcentrifuge to collect liquid.
- 150 Leave at room temperature for 5 minutes to elute DNA.
- 151 Place the tube strip in a magnetic separation rack until beads separate fully from the solution.
- 152 Slowly pipette off the cleared supernatant without disturbing the beads. Transfer supernatant to a new tube strip. Discard old tube strip with bead

## Nuclease Treatment

- 153 Add the following components in the order and volume listed below to a new microcentrifuge tube.  
Adjust component volumes for the number of samples being prepared, plus 10% overage. For



individual preps, add components directly to each sample from the previous step in the order and volume listed below.

Reaction Mix 5 (RM5):

A	B
Component	Volume
Nuclease Buffer	5 $\mu$ L
Nuclease Mix	5 $\mu$ L
Total Volume	10 $\mu$ L

- 154 Pipette-mix RM5.
- 155 Quick-spin RM5 in a microcentrifuge to collect liquid.
- 156 Add 10  $\mu$ L of RM5 to each sample. Total volume should equal 50  $\mu$ L.
- 157 Pipette-mix each sample with wide bore tips.
- 158 Quick-spin the strip tube in a microcentrifuge to collect liquid.
- 159 Run the nuclease treatment program.

A	B	C
Step	Time	Temperature
1	60 Minutes	37°C
2	Hold	4°C

## Final cleanup with SMRTbell cleanup beads

- 160 Add 60  $\mu$ L SMRTbell cleanup beads to each sample from the previous step. Using wide bore tips, pipette-mix the beads until evenly distributed.
- 161 Quick-spin the tube strip in a microcentrifuge to collect all liquid.



- 162 Leave at room temperature for 10 minutes to allow DNA to bind beads.
- 163 Place the tube strip in a magnetic separation rack until beads separate fully from the solution.
- 164 Slowly pipette off the cleared supernatant without disturbing the beads. It is recommended to save the supernatant in another tube strip in case of poor DNA recovery.
- 165 Slowly dispense 200  $\mu$ L, or enough to cover the beads, of freshly prepared 80% ethanol into each tube. After 30 seconds, pipette off the 80% ethanol and discard.
- 166 Repeat the previous step.
- 167 Remove residual 80% ethanol: • Remove tube strip from the magnetic separation rack. • Quick spin tube strip in a microcentrifuge. • Place tube strip back in a magnetic separation rack until beads separate fully from the solution. Pipette off residual 80% ethanol and discard.
- 168 Remove the tube strip from the magnetic rack. Immediately add 20  $\mu$ L of elution buffer to each tube and resuspend the beads by pipetting 10 times or until evenly distributed with wide bore tips.
- 169 Quick-spin the tube strip in a microcentrifuge to collect liquid.
- 170 Leave at room temperature for 5 minutes to elute DNA.
- 171 Place the tube strip in a magnetic separation rack until beads separate fully from the solution.
- 172 Slowly pipette off the cleared supernatant without disturbing the beads. Transfer supernatant to a new 0.5 mL LoBind tube. Discard old tube strip with beads.
- 173 Take a 1  $\mu$ L aliquot from each tube. Measure DNA concentration with a Qubit fluorometer using the 1x dsDNA HS kit. Calculate the total mass. Recommended: Further dilute each aliquot to 250 pg/ $\mu$ L with Femto Pulse dilution buffer. Measure final SMRTbell library size distribution with a Femto Pulse system.
- 174 Proceed to SMRT Link Sample Setup to prepare the SMRTbell library for sequencing.



- 175 Store SMRTbell libraries at 4°C if sequencing within the week. Long-term storage should be at -20°C. Minimize freeze-thaw cycles when handling SMRTbell libraries