



Apr 21, 2021

# Human Islet Cell Culture

Bridget Wagner<sup>1</sup><sup>1</sup>Broad Institute

1

Works for me

This protocol is published without a DOI.

HIRN



bwagner

## ABSTRACT

Human islet cell culture

## PROTOCOL CITATION

Bridget Wagner 2021. Human Islet Cell Culture. **protocols.io**  
<https://protocols.io/view/human-islet-cell-culture-budyns7w>



## LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Apr 21, 2021

## LAST MODIFIED

Apr 21, 2021

## OWNERSHIP HISTORY

Apr 21, 2021  avuApr 21, 2021  bwagner

## PROTOCOL INTEGER ID

49304

### Cell-culture formulations

- 1 The HTB-9 bladder carcinoma cell line (ATCC) is cultured in medium containing RPMI-1640 supplemented with 10% fetal bovine serum (FBS) and 1 x penicillin-streptomycin.
- 2 Islet culture medium: Islets are cultured in CMRL-1066 supplemented with 10% FBS, 1 x penicillin-streptomycin, and 2mM L-glutamine.

### Conditioned medium: Preparation of matrix-coated 384-well plates

- 3 Culture HTB-9 cells in T-175 flasks until 100% confluent, then incubate for an additional 3 days after changing media.
- 4 Collect the conditioned HTB-9 supernatant, and filter with 0.22 mm filter and store at -20C until use.

- 5 The day before you plan to seed islet cells, coat each well of a 384-well plate (PerkinElmer, CellCarrier-384 Ultra) with 10 mL/well HTB-9 supernatant and incubate overnight at 37C, 5% CO<sub>2</sub>. If you are short of time, 1-2 h is acceptable.
- 6 Wash each well twice with 50 mL/well PBS.
- 7 Add 10 mL/well islet culture medium (see above) and return to the incubator in advance of receiving islets.
- 8 As an alternative, coated plates filled with 50 mL/well PBS (after the same two washes) can be stored at 4C for up to 1 month.

#### Islet culture in 384-well plates

- 9 Upon receipt, pellet islets by centrifugation in a 50-mL Falcon tube and wash once with PBS.
- 10 Dissociate islets in Accutase (5000 IEQ/mL) in a Falcon tube for 20 min at 37C in 5% CO<sub>2</sub> incubator, keeping cell clumps down to ~10 cells. After ~10 min, agitate the Falcon tube and return to the incubator.
- 11 Pipette 5-7 times gently.
- 12 Add CMRL medium to dilute Accutase (9 mL for every 1 mL Accutase).
- 13 Resuspend islets in islet culture medium (see above) at a density of 250,000 cells/mL.
- 14 Dispense 40 mL/well into coated 384-well plates, which will give you 10,000 cells/well. We use the MultiDrop Combi (Thermo Fisher) for semi-automated dispensing. NOTE: this section can be adapted to 96-well format by altering the volumes and cell density. NOTE: this section can be adapted to 96-well format by altering the volumes and cell density.

#### Cell proliferation

- 15 After incubating for 24hr at 37C, 5% CO<sub>2</sub>, aspirate medium, leaving 10 mL/well. We use the BioTek ELx-405 automated plate washer for this step. Unless otherwise indicated, we aspirate medium so that 10 mL remains, and then add 40 mL/well to bring the total to 50 mL.
- 16 Add 40 mL/well medium containing 12.5 mM EdU (final concentration = 10 mM). Total volume is 50 mL/well.
- 17 Treat and incubate according to your experiment. We transfer compounds from compound plate (Abgene 384 PP, V-bottom) to cell plates using the Cybio Well Vario (Analytik Jena). We have found that the final DMSO concentration should remain  $\leq 0.3\%$ .

Change the culture medium containing 10 mM EdU (and add compounds as appropriate) every two days.

18

19 Aspirate the medium, then add 40 mL/well 4% PFA/PBS (final concentration will be 3.2%) at room temperature for 20 min.

20 Wash cells twice with PBS.

21 Permeabilize islet cells with 0.2% TritonX-100/PBS (final concentration) at room temperature for 20min.

22 Wash twice in PBS.

23 Stain cells using the Click-iT EdU HCS assays (catalog number C10351, Invitrogen) according to the manufacturer's protocol.

24 Block with 2% BSA/PBS at room temperature for 2 hr while gently shaking. No need to wash after this step.

25 Incubate primary antibody (C-peptide) at 1:100 dilution in 2% BSA/PBS overnight at 4C with gentle shaking.

26 Wash cells three times with PBS, then once with 2% BSA/PBS.

27 Incubate with Alexa Fluor 594-conjugated secondary antibodies (1:1000 dilution) and HCS NuclearMask Blue stain for 1 hr at room temperature.

28 Wash cells five times in 80 mL/well PBS. On final wash, leave wells filled with 50 mL/well PBS.

29 Seal plates with foil seal and store at 4C before imaging.