



Jan 07, 2022

Detection of Cryptosporidium in stool samples by Taq-Man qPCR

botchiesenyo¹¹University of Ghana

2

dx.doi.org/10.17504/protocols.io.bwghpbt6

Ayi and Pawlowic Lab Collaborations



mcpawlowic

Detection of *Cryptosporidium parvum* and/or *Cryptosporidium hominis* from fecal samples by qPCR.

DOI

dx.doi.org/10.17504/protocols.io.bwghpbt6

botchiesenyo 2022. Detection of *Cryptosporidium* in stool samples by Taq-Man qPCR. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bwghpbt6>



protocol ,

Jul 09, 2021

Jan 07, 2022

51433

Generating samples for a standard curve 25m

- 1 To 100 mg of fecal material (known to not be infected with *Cryptosporidium*), spike in a known number of *Cryptosporidium parvum*.^{25m}

DNA Extraction 2h 30m

- 2 Follow protocol from Quick-DNA Fecal/Soil Microbe Miniprep Kit from Zymo Research to extract total genomic DNA from stool sample (both standard curve samples and experimental samples).^{2h 30m}

Stool samples (0.1 g) should be thoroughly vortexed in lysis buffer and then subjected to five freeze-thaw cycles prior to DNA extraction.

qPCR

- 3 Combine the following for a single qPCR sample (9µl per sample). Scale accordingly for number of samples and technical replicates: 45m

5 µl	10 x Luna Universal Probe qPCR Master Mix
0.25 µl	10 µM Forward Primer (ATGACGGGTAACGGGGAAT)
0.25 µl	10 µM Reverse Primer (CCAATTACAAAACCAAAAGTCC)
1 µl	1 µM Probe ([FAM]CGCGCCTGCTGCCTTCCTTAGATG[BHQ1])
2.5 µl	Ultrapure water

- 4 Load mastermix in 96-well PCR plate using a multichannel pipette. 1h

Add 1µl extracted DNA to corresponding wells.

Keep a detailed plate map of which wells contain standard curve DNA, sample DNA, positive control DNA, no DNA control.

- 5 Spin plate in a centrifuge at low speed to make sure the samples are at the bottom of the well. 5m

- 6 Run qPCR cycling as described for Luna Universal qPCR Mixture. Make sure qPCR machine is measuring FAM at the correct cycle: 1h

6.1 95°C - 3 minutes 3m

6.2 57m
Repeat 40x:
95°C - 10 seconds
60°C - 30 seconds (record FAM in this step)

Calculating number of oocysts/g 2h

- 7 Calculate delta Ct values for samples that contain the standard curve. 1h

Using Microsoft excel (or a similar graphing program) plot the delta Ct values on the Y-axis and the number of oocysts/g on the x-axis.

Use the analysis tools to determine the linear equation of your standard curve samples.

- 8 Use the equation determined from your standard control samples, you can estimate the number of oocysts/g in your experimental samples. 1h

Y = delta Ct value for experimental sample

Solve for X = oocysts/g for your experimental sample