



Feb 28, 2022

# 🌐 Untargeted Top-down Proteomics by LC-MS/MS on Lumos

 [Untargeted Top-down Proteomics by LC-MS/MS on Eclipse](#)

Bryon Drown<sup>1</sup>, Jeannie Camarillo<sup>1</sup>, Neil Kelleher<sup>1</sup>

<sup>1</sup>Northwestern University

1



[dx.doi.org/10.17504/protocols.io.bttnnnme](https://dx.doi.org/10.17504/protocols.io.bttnnnme)

Human BioMolecular Atlas Program (HuBMAP) Method Development Community

Kelleher Research Group



Kelleher KRG Research Group

Northwestern University, National Resource for Translational...

Describes the LC-MS/MS data acquisition procedure for top-down proteomics samples using the Thermo Scientific Orbitrap Fusion Lumos Tribrid mass spectrometer

DOI

[dx.doi.org/10.17504/protocols.io.bttnnnme](https://dx.doi.org/10.17504/protocols.io.bttnnnme)

Bryon Drown, Jeannie Camarillo, Neil Kelleher 2022. Untargeted Top-down Proteomics by LC-MS/MS on Lumos. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.bttnnnme>



Office of the Director of the National Institutes of Health

Grant ID: UH3 CA24663

National Institute of General Medical Sciences of the National Institutes of Health

Grant ID: P41 GM108569

National Institute of Cancer of the National Institute of Cancer of the National Institutes of Health

Grant ID: F32 CA246894

[Untargeted Top-down Proteomics by LC-MS/MS on Eclipse, Kelleher KRG Research Group](#)

\_\_\_\_\_ protocol ,

Mar 30, 2021

Feb 28, 2022

48718

In steps of

[Overall protocol for top-down LC-MS/MS of human heart tissue](#)

[Overall protocol for top-down LC-MS/MS of human small intestine tissue](#)

[Overall protocol for top-down LC-MS/MS of human spleen tissue](#)

PLRP-S 5- $\mu$ m particles 1,000-Å pore size (Agilent Technologies)

Water Optima LC/MS Grade (Fisher Scientific #W64)



Acetonitrile Optima LC/MS Grade (Fisher Scientific #A955-4)

Formic Acid LC/MS Grade (Thermo Scientific #28905)

15  $\mu$ m SilicaTip PicoTip Emitter (New Object #FS360-50-15-N-20-C12)

Buffer A: 94.8 % water, 5 % acetonitrile, 0.2 % formic acid

Buffer B: 4.8 % water, 95 % acetonitrile, 0.2 % formic acid

- 1 Samples were analyzed on a Thermo Scientific Orbitrap Fusion Lumos Tribrid mass spectrometer in line with a Dionex Ultimate 3000 RSLCnano system
- 2 Samples (  6  $\mu$ L ) were injected via the autosampler and loaded onto a self-packed trap<sup>10m</sup> column (150  $\mu$ m i.d. x 2 cm length packed with PLRP-S 5- $\mu$ m particles 1,000-Å pore size) for  00:10:00 with 100% loading buffer (94.8% water:5% acetonitrile:0.2% formic acid)
- 3 Following a valve switch and initiation of the nanopump at 300 nL/min (buffer A: 94.8 % water, 5 % acetonitrile, 0.2 % formic acid; buffer B: 4.8 % water, 95 % acetonitrile, 0.2 % formic acid), proteins were separated on a self-packed analytical column (75  $\mu$ m i.d. x 25 cm length packed with PLRP-S 5- $\mu$ m particles 1,000-Å pore size) according to the following gradient for fractions 1-4:

A	B	C
Time (min)	%B	Valve Position
0	5	10_1
10	5	1_2
13	15	
70	45	
72	95	
76	95	
80	5	
90	5	

For fraction 5 and later, nanopump used the following gradient:

A	B	C	D	E	F	G	H	I	J	K	L	M	N	O	P	Q	R	S
Time (min)	%B	Valve Position																
0	5	10_1																
10	5	1_2																
13	15																	
70	50																	
72	95																	
76	95																	
80	5																	
90	5																	

- 4 Eluted proteins were ionized in positive ion mode nanoelectrospray ionization (nESI) using a pulled tip nanospray emitter (15- $\mu$ m i.d.  $\times$  125 mm) packed with 1mm of PLRP-S 5- $\mu$ m particles 1,000-Å pore size with a custom nano-source ([https://proteomicsresource.washington.edu/docs/protocols05/UWPR\\_NSI\\_Source.pdf](https://proteomicsresource.washington.edu/docs/protocols05/UWPR_NSI_Source.pdf)).

A	B
	<b>High-High</b>
<b>Spray voltage</b>	1800
<b>Sweep gas</b>	0
<b>Ion transfer tube temp</b>	320
<b>Application mode</b>	Intact Protein
<b>Pressure mode</b>	Low Pressure
<b>Advanced Peak Determination</b>	True
<b>Default charge state</b>	15
<b>S-lens RF</b>	30
<b>Source fragmentation</b>	15 eV

Global MS parameters

## 5 Precursor (intact protein) spectra were acquired at 120k FTRP

A	B
	<b>High-High</b>
<b>Detector type</b>	Orbitrap
<b>Resolving power</b>	120000
<b>m/z RP measured</b>	200 m/z
<b>Scan range</b>	600-2000
<b>Mass range</b>	Normal
<b>AGC target</b>	1000000
<b>Normalized AGC target</b>	250%
<b>Max Injection Time</b>	100 ms
<b>Microscans</b>	4
<b>Data type</b>	Profile
<b>Polarity</b>	Positive
<b>Use wide quad isolation</b>	True

Parameters for MS1 acquisition

- The mass spectrometer was operated using a Top2 data-dependent acquisition mode. Precursor ions were filtered by intensity, charge state, and dynamic exclusion.

<b>A</b>	<b>B</b>
<b>Intensity minimum</b>	20000
<b>Intensity maximum</b>	1E20
<b>Included charge states</b>	6-60
<b>Include undetermined charge states</b>	False
<b>Dynamic exclusion after n times</b>	1
<b>Dynamic exclusion duration</b>	60 s
<b>Mass tolerance</b>	1.5 m/z
<b>Exclude isotopes</b>	True

Precursor selection filters for DDA

## 7 Ions for fragmentations were isolated and fragmented via higher energy dissociation (HCD)

A	B
	<b>High-High</b>
<b>Detector type</b>	Orbitrap
<b>Isolation mode</b>	Quadrupole
<b>Resolving power</b>	60000
<b>m/z RP measured</b>	200 m/z
<b>Scan range</b>	350-2000
<b>AGC target</b>	1000000
<b>Normalized AGC target</b>	2000%
<b>Max injection time</b>	400 ms
<b>Microscans</b>	4
<b>Isolation window</b>	3 m/z
<b>Activation type</b>	HCD
<b>Collision energy</b>	27
<b>Collision energy mode</b>	Fixed
<b>Polarity</b>	Positive

Parameters for MS2 acquisition