



May 25, 2022

Expression and purification protocol of GST-mCh-FYVE

Expression and purification protocol of GST-NDP52

Liv Jensen¹, Chunmei Chang²

¹James Hurley Lab, University of California Berkeley; ²James Hurley Lab, UC, Berkeley

1



dx.doi.org/10.17504/protocols.io.ewov1n39pgr2/v1



Liv Jensen

This protocol details the expression and purification of GST-mCherry-FYVE.

[246-484.docx](#)

DOI

dx.doi.org/10.17504/protocols.io.ewov1n39pgr2/v1

Liv Jensen, Chunmei Chang 2022. Expression and purification protocol of GST-mCh-FYVE. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.ewov1n39pgr2/v1>



Expression and purification protocol of GST-NDP52, Chunmei Chang

ASAPCRN

protocol ,

May 02, 2022

May 25, 2022

61821

General information:

A	B
Expression system	E.Coli BL21DE3
Medium	Luria Bertani
Plasmid origin	
Backbone	pGST2
Resistance	Amp
Insert	
Tags & cleavage sites	N-term GST
Ext coeff	82320 M-1cm-1, MW 77.9 kDa

Lysis Buffer:

A	B
Hepes pH=7.5	50 mM
NaCl	300 mM
TCEP	1 mM
Protease Inhibitors (Roche)	

Wash Buffer:

A	B
Hepes pH=7.5	50 mM
NaCl	150 mM
TCEP	1 mM

Elution Buffer:

A	B
Hepes pH=7.5	50 mM
NaCl	150 mM
TCEP	1 mM
Glutathione	50 mM

SEC Buffer:

A	B
Hepes pH=8	20 mM
NaCl	150 mM
TCEP	1 mM





Columns/Resin:

- Glutathione Sepharose 4B

Protein expression 20m

- 1 Transform the E.Coli BL21DE3 cells with plasmid encoding for GST-mCh-FYVE and plate them on Amp plate.





2 

Carry out protein expression in  **1 L** medium, induce with  **200 micromolar (μM)** IPTG (isopropyl- β-d-thiogalactopyranoside) to an OD₆₀₀ of 0.8 and grow at  **18 °C**  **Overnight** .

3






20m

Harvest the cells by spinning at  **4500 x g** for  **00:20:00** at  **4 °C** and stock at  **-80 °C** until purification.

Protein purification 5h



- 4 Follow the GST batch purification by Size Exclusion Chromatography.


- 5 Resuspend the pellets in Lysis Buffer, sonicate for cell lysis and clear at  **16000 rpm** at ^{1h}
 **4 °C** for  **01:00:00**

6



30m

Incubate the supernatant with Glutathione Sepharose 4B (GE Healthcare) at  **4 °C** with gentle shaking for  **00:30:00** , apply to a gravity column, and wash extensively with Wash Buffer.

- 7 Elute the protein of interest with Elution Buffer and then apply onto a Superdex 6 column (10/300 Increase) pre-equilibrated in SEC Buffer at  **4 °C** .

- 8 Pool the peak fractions containing pure protein, snap-frozen in liquid nitrogen, and store at -80°C .