



Oct 05, 2022

Obtain LB broth medium

An.Huang¹¹XJTLU

1 Works for me

 Sharedx.doi.org/10.17504/protocols.io.5qpveryxdv4o/v1 An.Huang

DISCLAIMER

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

ABSTRACT

LB broth medium is a basic culture media for microorganisms. This protocol helps make an ordinary LB broth medium.

DOI

dx.doi.org/10.17504/protocols.io.5qpveryxdv4o/v1

PROTOCOL CITATION

An.Huang 2022. Obtain LB broth medium. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.5qpveryxdv4o/v1>



LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Oct 05, 2022

LAST MODIFIED

Oct 05, 2022

PROTOCOL INTEGER ID

70864

DISCLAIMER:

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

1 Dissolve  **25 g** LB broth powder in  **1000 mL** distilled water.

2 Autoclave the medium.



Use autoclave safely.

3 Store the medium at  **4 °C**.