



6

Jun 01, 2022

2-step PCR mixture and conditions (Barcoded-head primers for seqs pooling) V.6

Yin-Tse Huang¹¹Kaohsiung Medical University

protocol .



Yin-Tse Huang

PCR mixture and condition (2X SUPERGREEN PCR MASTER MIX)

Yin-Tse Huang 2022. 2-step PCR mixture and conditions (Barcoded-head primers for seqs pooling).

protocols.io

<https://protocols.io/view/2-step-pcr-mixture-and-conditions-barcoded-head-pr-cadfsa3n>

Yin-Tse Huang

protocol

Herbold CW, Pelikan C, Kuzyk O, Hausmann B, Angel R, Berry D, Loy A. 2015. A flexible and economical barcoding approach for highly multiplexed amplicon sequencing of diverse target genes. Front. Microbiol. [Internet] 6:731. Available from: <http://dx.doi.org/10.3389/fmicb.2015.00731>

protocol ,

Jun 01, 2022

Jun 01, 2022

63623

- 1 Wear glove, clean up the working bench w. 1% bleach

For 1' PCR head-primers

- 2 Prepare 1' PCR master mixutre for **head-primers** (prepare 1.2X of solutions for pipetting error if needed)

PCR mixture for head-primers for each reaction

A	B	C	D
Component	Volume	Volume (1.2X)	Final conc.
Forward Primer (10 µM)	1.6 µl	1.9 µl	1 µM
Reverse Primer (10 µM)	1.6 µl	1.9 µl	1 µM
2X Supergreen PCR Master Mix	7.8 µl	9.4 µl	-
ddH2O	4.1 µl	4.9 µl	-
Total volume	15 µl	18 µl	-

Negative control ALWAYS NEEDED! For example, if you have 5 PCR reactions to run, prepare master mixture for 6 reactions (5 DNA template + 1 negative control).

- 3 Mix the 1' PCR master mixture gently by pippeting. Quick spin the tube.

- 4 Transfer **15 µL** 1' PCR master mixutre in 8-strip PCR tubes.

- 5

Add **0.6 µL** DNA template in 8-strip PCR tubes, resulting in a **15.6 µL** reaction mixture for 1' PCR.

Negative control contains only **15 µL** master mixture but not DNA template

6 Mix the reaction mixture gently by tapping the tubes. Quick spin the tubes.

7 Carry out PCR using the following condition:

1' PCR condition for **head-primers**

A	B	C	D
Step	Temp	Sec	Cycle
Initial denaturation	95 °C	30-180 (a)	
Denaturation	98 °C	15	20-25 cycles
Annealing	64-68 °C varied (b)	15	
Extension	72 °C	60-180 (c)	
Final extension	72 °C	210	
Preservation	Preservation	4 °C	∞

a. Varied depend on template complexity


b. Annealing varied, **62-65°C** is working based on test on 220530; Refer to 1' PCR primers for annealing temperature

c. 1kb ~ 1min extension; enough time allow full extension of sequence

7.1 1' hear-primers used in Huang lab


A	B	C
Name	Tm°C	CG%
NS1B1ngs_H1	71.7	60
ITS4ngs_H1	69	55.6
LR5_H1	70.2	60.6
EF1-526F_H1	71	59.3
EF1-1567R_H1	70.6	56
EF1-2218R_H1	72.2	60.3
Ben2f_H1	70.5	56.8
Bt2b_H1	74.5	62.5
T22_H1	70.3	56.8
RPB2-3bF_H1	65.6	48.7
RPB2-7cR_H1	72.3	59.7
fRPB2-11aR_H1	71.7	60.8

8 Carry out **electrophoresis** for inspection of DNA products



Gel before markdown

9 Markdown wells and upload the pictures to the Lab Google drive



Marked gel picture go to the Lab Google drive

For 2' PCR barcoded-head primers

- 10 Prepare 2' PCR master mixutre for **barcoded-primers** (prepare 1.2X of solutions for pipetting error if needed)

PCR mixture for barcoded-primers for each reaction (**NO PRIMERS at this point!!**)

A	B	C	D
Component	Volume	Volume (1.2X)	Final conc.
2X Supergreen PCR Master Mix	10.75 µL	12.9 µL	-
ddH2O	10.75 µL	12.9 µL	-
Total volume	21.5 µL	25.8 µL	-

Negative control ALWAYS NEEDED! For example, if you have 5 PCR reactions to run, prepare master mixture for 6 reactions (5 DNA template + 1 negative control).

- 11 Mix the 2' PCR master mixture gently by pipetting. Quick spin the tube.
- 12 Transfer **21.5 µL** of the 2' PCR master mixture to PCR tubes.
- 13 Add **2.5 µL pre-mixed barcoded-head primers** (Forward + Reverse) to each PCR tube

14

Add **1 µL** of **1' PCR product as template**, resulting in **25 µL** reaction mixture for 2' PCR.

Negative control contains only **24 µL** master mixture and premixed barcoded-head primers but not DNA template

- 15 Mix gently by tapping the tubes. Quick spin the tubes.

- 16 Carry out 2' PCR using the following condition:

2' PCR condition for **barcoded-head primers**


A	B	C	D
Step	Temp	Sec	Cycle
Initial denaturation	98 °C	30	
Denaturation	98 °C	15	10-15 cycles
Annealing	64-68 °C varied (a)	15	
Extension	72 °C	60 (b)	
Final extension	72 °C	210	
Preservation	Preservation	4 °C	∞

a. Annealing varied, **65 °C** is working based on test on 220531; Refer 2' PCR primers for annealing temperature
b. 1kb ~ 1min extension; enough time allow full extension of sequence

16.1 2' barcoded-head primers used in Huang lab


A	B	C	D	E	F	G	H	I	J
Name (1st)	Sequence	Tm °C	CG%	Name (2nd)	Sequence	Tm °C	CG%	Name (3rd)	Sequence
F1-1	aaatctcattttctgtagtgcgctatgctgcgagctgc	71.7	53.8	F2-1	ggccattctgtagtaagtagttcagctatgctgcgagctgc	70	53.8	F3-1	gggctagagtttaaccttaattag
F1-2	aaccgcgacgtctattttaagaagctatgctgcgagctgc	70.1	51.3	F2-2	cagatactaatctcgtcacatggctatgctgcgagctgc	68.5	51.3	F3-2	taatttatggctacagacgttag
F1-3	aactgtttattttcaaatatagctatgctgcgagctgc	61.3	35.9	F2-3	agagattgagcttattcgtttttgctatgctgcgagctgc	68.7	48.7	F3-3	gtaagcgtctacacagcaaaact
F1-4	aatatgaagctccgacatatctggctatgctgcgagctgc	69.5	51.3	F2-4	ccattgattccagattatcatgtgctatgctgcgagctgc	67.7	48.7	F3-4	atccgcatactcaactcgaaagc
F1-5	acagtacaaacacggtcttaagctatgctgcgagctgc	70	51.3	F2-5	cgttaggactctgtagtttaaacgctatgctgcgagctgc	70.1	53.8	F3-5	gtgctgtcaaatcgttttctatgc
F1-6	acagttgacctcgcgacattatgctatgctgcgagctgc	71.2	53.8	F2-6	gtcatgaccgctacatatctcagctatgctgcgagctgc	70	53.8	F3-6	gtgggataagcttgacattttag
F1-7	caccataaatgagattgctgaggtatgctgcgagctgc	70.7	53.8	F2-7	cgataggacaagcaatgtactcagctatgctgcgagctgc	70.1	53.8	F3-7	gttgataggaggtcaaacaga
F1-8	accgtgtatgtattcgtgttaccgctatgctgcgagctgc	71	53.8	F2-8	cgatataatcgatccgcataccgctatgctgcgagctgc	71.4	56.4	F3-8	gttaaccataagtgcgacctat
R1-1	acgtcaatgctattcccagtcagctatgctgcgagctgc	71.4	53.8	R2-1	cagcgtttccaaagacattattgctatgctgcgagctgc	69.7	51.3	R3-1	taaaaaattcgtggaactccacag
R1-2	acgtcgaggtatcataaatacttgcctatgctgcgagctgc	67.9	48.7	R2-2	cgctttatgtcttaagtagccgctatgctgcgagctgc	69.1	51.3	R3-2	taatacattgtgtctattcttagagc
R1-3	cgcaagtatgctcttctcatagctgctatgctgcgagctgc	70.1	53.8	R2-3	ctgctgagttataccacagtgacgctatgctgcgagctgc	71.5	56.4	R3-3	gtaaatctaggtgtaaaatgagt
R1-4	agtcattgctttgggtacataaagctatgctgcgagctgc	70.1	51.3	R2-4	tacgttatattaactctagccgagctatgctgcgagctgc	67.4	48.7	R3-4	ctttctataataatccgggctaagc
R1-5	atatcttacacaaaagtagctgctatgctgcgagctgc	65.4	43.6	R2-5	tagaatgtcaacacaaagtaggacgctatgctgcgagctgc	69.5	51.3	R3-5	ctttgttaattggtgttccgttgt
R1-6	cgactatttaacttcgcgaacagctatgctgcgagctgc	69.3	51.3	R2-6	gaaatcgacaaaattctgctctgctatgctgcgagctgc	69.7	51.3	R3-6	taggtatcattctcatcctatcggc
R1-7	atcgttggtctatgttcaggtatgctgcgagctgc	69.5	51.3	R2-7	gaatcatcaagaagggaacacagctatgctgcgagctgc	69.3	51.3	R3-7	tatagagacgggtttcggtaaaa
R1-8	atgagatctatctagtagccgttgctatgctgcgagctgc	68.8	51.3	R2-8	gactatagtgaaaaatcacatagctatgctgcgagctgc	65.8	46.2	R3-8	tattagtagatagacactcggg
R1-9	gaggagccatagatgaatgagctatgctgcgagctgc	69.7	53.8	R2-9	atgatcctacggagatttacctgctatgctgcgagctgc	70.7	53.8	R3-9	tcagggttagaagactagttgta
R1-10	attgtagcattgaataggagcagctatgctgcgagctgc	68.7	48.7	R2-10	gagttctggatctatctgggtgctatgctgcgagctgc	70	53.8	R3-10	tcgggttaaatgctaagcgtaat
R1-11	caaatctcatgctagctgacgctatgctgcgagctgc	70.2	53.8	R2-11	gcaacgaaacattcgttaagtagtctatgctgcgagctgc	68.1	48.7	R3-11	tgctgaaaaacaggaagtctcac
R1-12	caagaactagacgctgctcttaagctatgctgcgagctgc	70.2	53.8	R2-12	tgcccaatatgttagcaccctaagctatgctgcgagctgc	71.4	53.8	R3-12	accattcctaataagccaataggg

17 Carry out electrophoresis for inspection of DNA products



Gel before markdown

18 Markdown wells and upload the pictures to the Lab Google drive



Marked gel picture go to the Lab Google drive

